

Denture 4U System

Denture 4U KIT

# Denture 4U Kit




| HIOSSEN'S FIXED DENTURE KIT

| USER MANUAL

**HIOSSEN**  
IMPLANT

# Denture 4U Kit

- Accurate Surgery with Precision Guide Positioning
- Safe Surgery without Nerve Damage



## Perform Denture 4U Treatment

When there is vertical bone loss due to alveolar bone resorption

## Denture 4U KIT 0° Posterior Guide can be used to treat patients

who do not lack of vertical bone mass

## Up to 6 implants can be placed by using the Denture 4U KIT

to acquire high fixation power on the maxillary bone with soft osseous tissue

## Place 4~6 Implants in tilted manner

- Avoid inferior alveolar nerve and gain stability by placing 4~6 implants in a tilted manner in cases where the patient lacks of bone volume.
- In cases where the implant are tilted, the cantilever length can be reduced, which disperses the load efficiently on just 4 implants and thus making denture treatment possible.

※ It is advised to place 6 implants in the maxilla to secure stability.

## Semi-Permanent use owned to Fixed Full-Denture


- No need for re-lining which is usually needed due to gum recession.
- Unlike removable dentures, there is no need to replace abutment components.

## Excellent Aesthetics Compared to Conventional Dentures

- Denture 4U enables placement of 4 implants in a way that they can properly disperse the pressure from masticatory movements, and therefore prevents alveolar bone resorption and involution.
- Maintains shape and volume of the jawbone, which results into better esthetics than conventional dentures.

## Denture 4U Kit Line-up

### Implant




10mm 11.5mm 13mm 15mm  
ET

### Prosthetic




Multi ABT Multi Angled ABT Esthetic- low Cylinder Multi Ti- Base

### Denture 4U KIT




# Denture 4U KIT

KIT for Denture 4U Surgery : Enables accurate and safe Drilling



## Why is Denture 4U Kit Essential

- Most edentulous patients lack of vertical bone volume due to alveolar bone resorption.



### Denture 4U Treatment

- Placing long implants is possible in order to gain stability.
- Place Implants in tilted manner in order to reduce cantilever length.

Beware of the alveolaris inferior nerve since the long implants are inserted with an inclined angle.

### Denture 4U KIT

- Guides the placement site of the 4~6 implants in edentulous cases.
- Adjusts the angle and distance between implants.

Implant placement is safe while avoiding the inferior alveolar nerve.

# Denture 4U KIT SURGICAL SEQUENCE


## STEP 1 | Preparation

### STEP 2 | 1-point Fixation (refer to p.5)

### STEP 3 | 2-point Fixation (refer to p.8)

### STEP 4 | Drilling (refer to p.8)

### STEP 5 | Reaming (refer to p.10)




## STEP 1 | Preparation

※ Before the procedure, check the location and shape of inferior alveolar nerve, and involution of alveolar bone.

### 1 Bone Flattening


- Flatten the bone with the crest remover in order to set the conditions for Guide Positioning



### 2 Check the median line

- Find and set the median line by checking the labial frenum or the mid line of the nose and chin.

#### Guide 1 | Check labial frenum



Set the median line by checking the labial frenum.

#### Guide 2 | Check the midline of the nose and chine




Connect the midline of the nose and chin. The line will run over the alveolar bone, which will be the median line.

## STEP 2 | 1-point Fixation

### 1 1-point fixation in Anterior region

- Use an Anchor Screw to firmly fix the guide on the bone.



#### User Guide TIP

- First, try to fix the guide with an Anchor Screw. If the Anchor Screw can not be inserted because the bone quality is hard, use the Anchor Drill before placing the Anchor Screw.
- Soft Bone : Possible to fix guide with Anchor Screw.
- Normal/Hard Bone : Fix the guide with Anchor Screw after drilling a hole with the Anchor Drill.
- ※ Stop the engine when the mount driver reaches the guide in order to prevent tickover of the Anchor Screw.



#### Select specifications for Anchor Screw

- Perform drilling with the 3mm Anchor Drill first, before drilling with 11mm Anchor Drill.

※ There is no contact between the drill and the guide, if the surgeon performs the initial drill with the 11mm Anchor Drill.

### 2 Guide Positioning in Anterior region (adjust Anterior Guide)


- Position the guide according to the patient's dental arch.



- Position the guide and manually adjust ① according to the patient's dental arch
- In order to fix position ①, tighten ② with hand driver.

#### User Guide TIP


- Based on first premolar, the most narrow space is 31mm, and from that point on, the angle can be widened up to 30°
- Guide can be correctly positioned even on asymmetric dental arches, because each left & right side, anterior & posterior region can be adjusted separately.



# Denture 4U KIT SURGICAL SEQUENCE

## 3 Guide Positioning in Posterior region (adjust Posterior Guide)

- Adjust and fix the Posterior Guide according to the patient's dental arch.



Posterior Guide  
Option : 0°, 17°, 30°


※ The Posterior Guide enables adjustment of the Implant's placement angle, distance between Implants and the buccal/lingual angle. The surgeon can therefore place the Implants in the desired way without damaging the alveolar nerves.

### Usage Guide **TIP**

- Adjust placement angle, distance between implants and buccal/lingual angle, and then tighten screw with hand driver to fix the adjustments.


#### 1 | Adjust Placement Angle

Posterior guide can be changed during the surgery, but it is advised to select appropriate specification with CT image before the surgery.



#### 2 | Adjust Distance between Implants

Adjust the distance with the help of the laser markings that come in 2mm units.



#### 3 | Adjust the Buccal/Lingual Angle

Buccal/lingual angle can be adjusted up to  $\pm 35^\circ$




## 4 Check Surgery Safety (check alveolar nerve)

※ Needs to be checked before drilling stage.


- Before the drilling stage, the location of the alveolar nerve needs to be checked after positioning the guide.
- Denture 4U Guide knobs can be removed. (Better Panorama images can be acquired with CT checker, when knobs are removed.)

### Guide 1 | Check with Indicator



Perform a full flap surgery in order to spot the mental foramen with naked eye. Safety can be checked with the indicator.

### Guide 2 | Check with Path Checker




In case the mental foramen is not visible with the naked eye, place Path Checker and check location of the nerve with a CT image.


### Usage Guide **TIP**

#### Guide 1 | Locate with indicator


Case for visual confirmation of the mental foramen (Hold the indicator with instruments such as hemostat or needle holder.)



The alveolar nerves goes through the mental foramen, and therefore the drilling path should be more in the mesial direction.




Check drilling path by placing the indicator in correct direction depending on the guide that will be used.




Check the drilling path with naked eye and adjust distance between implants.

#### Guide 2 | Locate with path checker


Case for ordinary flap surgery



Place path checker inside the drilling hole and check drilling path on panorama or CT image.



Readjust the guide in case the drilling path passes through the alveolar nerve. (Laser markings come in 2mm units)



Posterior Guide adjustment is completed so that drilling path does not pass over alveolar nerves.


# Denture 4U KIT SURGICAL SEQUENCE

## STEP 3 | 2-point Fixation

### 1 2-point Fixation in Posterior region

- After placing the Guide positioning taking into consideration the alveolar nerve, the guide needs to be fixed on 2 points fixation in order to prevent movements of the guide.

※ Use Anchor Screws to fix Posterior Guide on 2 points. (The guide is then firmly fixed on 3 points, including the fixation in the Anterior region, and therefore drilling can be performed in a stable manner.)



#### User Guide TIP

- First, try to fix the guide with an Anchor Screw. If the Anchor Screw can not be inserted because the bone quality is hard, use the Anchor Drill before placing the Anchor Screw.

- Soft Bone : Possible to fix guide with Anchor Screw.
- Normal/Hard Bone : Fix the guide with Anchor Screw after drilling a hole with the Anchor Drill.

※ Stop the engine when the mount driver reaches the guide in order to prevent tickover of the Anchor Screw.

#### Choosing Anchor Screw Specification



- When 2-point fixing guide in Posterior region choose the appropriate Anchor Screw, depending on the severity of the bone resorption. (11mm Anchor Screws are available, in order to provide stable fixation in regions with severe bone resorption).
- Perform drilling with the 3mm Anchor Drill first, before drilling with 11mm Anchor Drill.

※ If the surgeon performs the initial drill with the 11mm Anchor Drill, there is no contact between the Drill and the Guide.

## STEP 4 | Drilling



### 1 Drilling in Posterior region (Ø3.0)

- Perform Drilling in Posterior region with Ø3.0 Twist Drill.



#### Usage Guide TIP

- Place the drill carefully into the guide hole by referring to the marking line which is marked in the lateral side of the guide.
- Control the drilling depth by referring to the drill's marking line in the mesial direction. Use the 5mm drill first and then the 18mm drill in case the surgeon uses a 0° Guide or experiences interference from occluding teeth.




Check mesial direction when referring to the markings.

#### Tips for Preventing jumps of the drills

Set angle of the drill by taking the guide angle into consideration, and press the pedal as you advance carefully with the drill. (If your hand is relaxed and the drill angle matches with the guide hole angle, the drill will glide into the hole and drill as planned.)

### 2 Drilling in Anterior region (Ø2.2)

- Perform drilling in Anterior region with Ø2.2 Twist Drill.




#### Checklist Before Anterior Drilling TIP

- ① Check whether the Dental arch's curve is the same in the Anterior and Posterior region.
  - In case the Guide does not fit due to the curve difference, re-position the guide before performing drilling in the Anterior region.
- ② Check whether the Posterior Guide is blocking the guide hole for the Anterior region.
  - In case it's blocked, remove the Posterior Guide first, and then perform drilling in Anterior region.



Posterior Guide is blocking the guide hole for Anterior region.



Perform drilling after removal of Posterior Guides.


→ Since the Posterior Guides need to be removed in the 2 cases above, firmly hold the Guide, which has then only 1-point fixation, and perform drilling.

# Denture 4U KIT SURGICAL SEQUENCE

## STEP 5 | Reaming

### 1 Countersink Drilling (to prevent interference from Stoppers and Prostheses)

- Remove Denture 4U Guide and perform Countersink Drilling in order to prevent interference from Taper Drill Stoppers and Prostheses.



**Interference from Stoppers and Prostheses**

When implants need to be placed in a tilted way, there will be interference from the drill's stopper and prostheses.


Initiate drilling performed with Denture 4U Guide.

Countersink drilling In order to prevent interferences.

**Placement Guide for Marking Line TIP**


Drill Type	Distal	Mesial
When using 17° Guide	Lower line (Bone level Implantation)	Upper line (1mm subcrestal Implantation)
When using 30° Guide	Lower line (Bone level Implantation)	Upper line (1mm subcrestal Implantation)

※ Refer to the lower line in the distal direction when Implant needs to be placed at bone level. Refer to the upper line in the distal direction when Implant needs to be placed 1mm subcrestal.




### 2 Drill hole expansion with Taper Drills

- Perform additional drilling with Taper KIT or 122 Taper KIT in order to have appropriate drill holes for the implant.



#### Taper KIT

The tapered drills form optimal drill holes for tapered implants that gain good initial stability in the alveolar bone.



#### 122 Taper KIT

A Kit with simple drill protocol: 1 drilling in soft bone, and 2 drillings in normal and hard bone.

# Denture 4U Prosthetic PROCESS


## Temporary Denture

- Abutment Placement**
- Impression Taking & Try in**
- Create through holes**
- Check Abutment Location**
- Place Temporary Cylinder**
- Seat Temporary Denture**
- Inject Resin to attach**
- Cut out excessive part of the Cylinder**
- Cut excessive parts of the Temporary Denture**
- Final seating and finishing**

※ Please refer to ET Prosthetics Manual for detailed fabrication protocol.

# Denture 4U Prosthetic PROCESS

## Final Denture



※ Please refer to ET Prosthetics Manual for detailed fabrication protocol.

# Denture 4U KIT ORDER CODE

## Denture 4U KIT | OD4UK

### Denture 4U Guide



### Posterior Guide


#### Degree

0°	D4UPG0
17°	D4UPG17
30°	D4UPG30



### Crest Remover

L	D Ø5.0
29	CERM50A
45	CERM50S



### Anchor Screw

L	D Ø1.65
5	D4UAS5
8	D4UAS8
11	D4UAS11



### Anchor Drill

L	D Ø1.65
3	D4UAD3
11	D4UAD11



### Guided Initial Drill

L	D Ø2.2
5	GD2208NC



### Twist Drill

L	D Ø3.0
5	D4U2D3005
18	D4U2D3018




### Countersink

D4UCS



### Indicator

D4UI



### Path Checker

D4UPC



### Simple Mount Driver

L  
Short ASMDS




### Multi Abutment Machine Driver

MAMD



### Multi Abutment Outer Driver

MAOD



# Denture 4U KIT ORDER CODE

## Prosthetic | ET



Multi Abutment




Multi Angled Abutment



Multi Abutment NP-Cast Cylinder



Esthetic-low Plastic Cylinder



Esthetic-low Temporary Cylinder




Multi Ti-Base


### Multi Abutment

G/H	1.0	2.0	3.0	4.0	5.0	
<b>M</b>	ETMTA501M	ETMTA502M	ETMTA503M	ETMTA504M	ETMTA505M	
<b>R</b>	ETMTA501R	ETMTA502R	ETMTA503R	ETMTA504R	ETMTA505R	

### Multi Angled Abutment

17°	G/H 2.5	3.0	4.0	5.0	
<b>M</b>	ETMA217MHW	ETMA317MHW	ETMA417MHW	-	
<b>R</b>	ETMA217SHW	ETMA317SHW	ETMA417SHW	-	
30°	G/H 3.5	4.0	5.0		
<b>M</b>	ETMA330MHW	ETMA430MHW	ETMA530MHW		
<b>R</b>	ETMA330SHW	ETMA430SHW	ETMA530SHW		


### Multi Ti-Base

H	5°	10°	
4	TSMTB0405GTH	TSMTB0410GTH	
6	TSMTB0605GTH	MGW100	

### Multi Ti-Base


TSMSBC	
--------	---

### Esthetic-low Temporary Cylinder

Hex	Non-hex	
MTR200	MTR100	

※ Regular Non-hex


### Esthetic-low Temporary Cylinder (Narrow type)

Hex	Non-hex	
NMTR200	NMTR200	


※ Regular Non-hex

※ Please refer to the Product Catalog for information on ET & EK system.

### Multi Abutment NP-Cast Cylinder

Hex	Non-hex	
TSMN500	TSMN500N	

### Esthetic-low Plastic Cylinder

Hex	Non-hex	
MGR200	MGR100	





# **HIOSSEN** IMPLANT

Smiles that last a lifetime

Please contact your local sales representative or  
visit our website today to learn more about Hiossen and its products.

270 Sylvan Ave. Ste 1130, Englewood Cliffs, NJ 07632  
Marketing@hiossen.com | 888.678.0001 | www.hiossen.com

 @hiossen\_implants

 @HiossenImplants

 @Hiossen

Copyright © 2022 Hiossen Inc. All rights reserved.

All information included herein is subject to change without notice. No part of this brochure may be reproduced in any manner without the prior written permission of Hiossen Inc. Hiossen and the Hiossen logo are trademarks of Hiossen Inc. All other brand or product names are trademarks of their respective companies or organizations.

PB22D4KLTR1.0