

Full Arch Treatment

# Multi Abutment

Multi Abutment


Multi Angled Abutment

**HÍOSSEN**  
IMPLANT




# Multi Abutment

Maximize your practice’s potential with the Multi Abutment’s exceptional clinical adaptability. This groundbreaking solution ensures a superior outcome, by overcoming the limitations of full arch treatments. Elevate your practice by incorporating the the Multi Abutment, the ultimate solution to enhance your capabilities.



## Case



Bridge

Fixed full denture

Bar with removable full denture

# The concept of Full-arch treatment

The full-arch treatment option offers the opportunity to craft extensive dental restorations with minimal implant placement, resulting in reduced costs and shorter treatment times.

The selection of an implant system can be customized to align with the surgical plan. Moreover, the use of components such as the Multi Abutment and Multi Angled Abutment provides flexibility, allowing for a range of customized restoration options tailored to each unique case.

Furthermore, advanced digital libraries and CAD/CAM technology streamline the fabrication of diverse restorations.






# Multi Abutment


## Gingival Heights

- A range of gingival heights, from 1mm to 5mm, are available.



## Combination Options

- For each case, a wide range of prosthetic components and cylinders can be utilized, offering numerous options for various types of prosthetic restoration.




- The Multi Abutment is attached to an internal submerged implant, which operates at its optimal performance.



# Multi Angled Abutment


## Two Angulation Options

- The Multi Abutment provides two angle options, 17 degrees and 30 degrees, serving as compensation angles for implants implanted during the prosthetic engagement stage. This versatility enables it to address a wide spectrum of clinical cases.



## Gingival Heights



- A wide range of gingival heights from 2.5mm to 5mm are available.



\* Options available in both 17 and 30 degrees

## Divergence Compensation


- Additional compensation can be made with the angle at the top of the abutment to compensate for up to 108 degrees.




# Components & Tools

## Conventional


### Cylinder




**Esthetic-low Temporary Cylinder**  
Used for fabrication of temporary prosthesis



**Multi Combination Cylinder**  
Used for producing a combination of prostheses



**Esthetic-low Gold Cylinder**  
Used for fabrication of customized prosthesis by casting with gold alloys



**Esthetic-low Gold Cylinder**  
Used for fabrication of customized prosthesis by casting with nonprecious metal alloys

### Impression



**Esthetic-low Transfer Coping**




**Esthetic-low Pick-up Coping**


### Lab Analog




**Esthetic-low Lab Analog**



**Esthetic-low Polishing Protector**  
Protects the joint in the polishing procedure after producing a prosthesis using esthetic GoldCast/plastic cylinder.



**Esthetic-low Healing Cap**  
Protection cap for esthetic-low Abutment.




**Multi Abutment Outer driver**  
A dedicated torque Driver for Straight Multi Abutment

## Digital

\*Digital libraries are available for the following manufacturers

### Cylinder



**Multi Ti Base**


- EK/ET Multi Abutment is utilized in the fabrication of combination-retained prostheses
- Used by connecting it to the EK/ET Multi Scan Body

3shape

exocad

dental wings  
A Straumann Group Brand

### Impression




**Multi Scan Body**  
Used for oral scanning by connecting to EK/ET Multi Abutment

### Lab Analog




**Multi Digital Lab Analog**



**Reamer Drill**  
Used for correction of insertion path of 3D-printed working model



**Positioning Jig**  
Used for inserting Lab Analog to the working model in the accurate hex direction








# Optional Components

## Temporary Restoration

Temporary prosthetics can be set immediately after surgery by the use of a Temporary Cylinder. It results in a high treatment satisfaction by minimizing dental patient’s functional, aesthetic, and psychosocial discomfort.



### Esthetic-low Temporary Cylinder





The temporary cylinder used during the temporary restoration stage has a total length of 12mm. It is conveniently adjustable in length due to its segmentation into 1mm sections from top to bottom.

## Digital Restoration

\*Digital libraries are available for the following manufacturers **3shape** **exocad** **dental wings**  
A Straumann Group Brand

The Multi Ti Base offers the advantages of digital manufacturing, which include enhanced precision and customization in prosthetics. It is compatible with a wide range of Multi/Multi Angled Abutments, and various software options provide access to digital libraries for streamlined workflows















# Multi Abutment

## Product Overview














|                         |    |
|-------------------------|----|
| Abutment                | 11 |
| Prosthetic Components   | 12 |
| Impression Components   | 12 |
| Placement Tool          | 13 |
| For Digital Prosthetics | 14 |

## Abutment

### Multi Abutment Packing unit: Abutment + Carrier\*)

| <br>D Ø 4.8<br>M | G/H | 1.0  | 2.0  | 3.0  | 4.0  | 5.0  |
|---|-----|--|--|--|--|--|
|   |     | <br>ETMTA501M | <br>ETMTA502M | <br>ETMTA503M | <br>ETMTA504M | <br>ETMTA505M |
| <br>D Ø 4.8<br>R | G/H | 1.0  | 2.0  | 3.0  | 4.0  | 5.0  |
|   |     | <br>ETMTA501R | <br>ETMTA502R | <br>ETMTA503R | <br>ETMTA504R | <br>ETMTA505R |


### Multi Angled Abutment Packing unit: Abutment + Carrier\*) + EbonyGold Screw(\*)

| <br>D Ø 4.8<br>M | Angle | G/H | 2.5   | 3.0   | 4.0   |
|---|-------|-----|---|---|---|
|   | 17°   |     | <br>ETMA217MHW | <br>ETMA317MHW | <br>ETMA417MHW |
|   | Angle | G/H | 3.5   | 4.0   | 5.0   |
|   | 30°   |     | <br>ETMA330MHW | <br>ETMA430MHW | <br>ETMA530MHW |
|   | Angle | G/H | 2.5   | 3.0   | 4.0   |
|   | 17°   |     | <br>ETMA217SHW | <br>ETMA317SHW | <br>ETMA417SHW |
|   | Angle | G/H | 3.5   | 4.0   | 5.0   |
|   | 30°   |     | <br>ETMA330SHW | <br>ETMA430SHW | <br>ETMA530SHW |

Prosthetic Components

Packing unit: Cylinder + Ti Cylinder Screw(\*)

Esthetic-low Temporary Cylinder

| D ∅ 4.8 \ Type | Non-Hex  |
|----------------|--|
|                | <div><div>Narrow Type</div><div></div><div>MTR100TH</div></div> |


Multi Combination Cylinder

| D ∅ 4.8 \ Type | Hex   | Non-Hex  |
|----------------|---|--|
|                | <div></div> <div>TSMC500TH</div> | <div></div> <div>TSMC500NTH</div> |

Multi Abutment NP-Cast Cylinder

| D ∅ 4.8 \ Type | Hex   | Non-Hex  |
|----------------|---|--|
|                | <div></div> <div>TSMN500TH</div> | <div></div> <div>TSMN500NTH</div> |

Esthetic-low Gold Cylinder

| D ∅ 4.8 \ Type | Non-Hex  |
|----------------|--|
|                | <div></div> <div>MGR100TH</div> |

Esthetic-low Plastic Cylinder

| D ∅ 4.8 \ Type | Non-Hex  |
|----------------|--|
|                | <div></div> <div>MEPRI00TH</div> |









\*Screw (MTS200)


Impression Components

Esthetic-low Pick-up Impression Coping


Packing unit: Impression coping body + Guide pin(\*)

| D ∅ 4.8 \ Type | Hex  | Non-Hex  | 5  | 10   | 12  | 15  |
|----------------|--|--|--|--|---|---|
|                | <div><div></div><div>MSR200</div></div> | <div></div> <div>MSR100</div> | <div></div> <div>GP100</div> | <div></div> <div>*GP150</div> | <div></div> <div>GP170</div> | <div></div> <div>GP200</div> |


Esthetic-low Transfer Impression Coping

| D ∅ 4.8   |
|---|
| <div></div> <div>MTTR100</div> |

Esthetic-low Lab Analog




| D ∅ 4.8   |
|---|
| <div></div> <div>MERR300</div> |

Esthetic-low Polishing Protector

| D ∅ 4.8   |
|---|
| <div></div> <div>MRCR100</div> |



Placement Tool

Torque Driver (1.2 Hex)

|   |  |  |
|---|--|--|
| <div></div> <div>Short</div> | <div></div> <div>Long</div> | <div><div></div><div>Extra long</div></div> |
|---|--|--|



| L \ Type   | 1.2 Hex  |
|------------|----------|
| 8 Ex.Short | TRHD12MS |
| 13 Short   | TRHD12S  |
| 25 Middle  | TRHD12M  |
| 20 Long    | TRHD12L  |
| 25 Ex.Long | TRHD12E  |

Hand Driver (1.2 Hex)



|   |  |
|---|--|
| <div></div> <div>Extra Short</div> | <div></div> <div>Long</div> |
|---|--|

| L \ Type   | 1.2 Hex  |
|------------|----------|
| 8 Ex.Short | AHD12MSH |
| 13 Short   | AHD12SH  |
| 25 Middle  | AHD12MH  |
| 20 Long    | AHD12LH  |
| 25 Ex.Long | AHD12EH  |



Multi Abutment Outer Driver

|  |  |
|--|--|
| <div></div> <div>MAOD</div> | <div></div> |
|--|--|

Multi Abutment Machine Driver

|  |  |
|--|--|
| <div></div> <div>MAMD</div> | <div></div> |
|--|--|

Angled Torque Driver

|   |  |
|---|--|
| <div></div> <div>Short</div> | <div></div> <div>Long</div> |
|---|--|


| L \ Type | 1.2 Hex | 1.2 Hex (Set) |
|----------|---------|---------------|
| 13 Short | ATD12S  | ATD12S3S      |
| 20 Long  | ATD12L  | ATD12L3S      |

# For Digital Prosthetics

\*Digital libraries are available for the following manufacturers



**TS Multi Ti Base** Packing unit: Ti Base + Ti Cylinder Screw(\*)



## Multi Digital Lab Analog

Packing unit :  
Digital Lab Analog  
+ Screw



TSMDLA

## Positioning Jig



TSMDLAPJ

## Multi Scan Body



TSMSBC


## Screw Bulk (10ea)



TSDLASB



**Reamer Drill** \*Use A Type for Multi Abutment




DLARDR



# **Hiossen** IMPLANT

Smiles that last a lifetime






Please contact your local sales representative or visit our website today to learn more about Hiossen and its products.



All Hiossen Implants are processed and  
Manufactured in the USA

270 Sylvan Ave. Ste 1130, Englewood Cliffs, NJ 07632  
Marketing@hiossen.com | 888.678.0001 | www.hiossen.com

 @Hiossen     @HiossenImplants     @Hiossen

Copyright © 2024 Hiossen Inc. All rights reserved.

All information included herein is subject to change without notice. No part of this brochure may be reproduced in any manner without the prior written permission of Hiossen Inc. Hiossen and the Hiossen logo are trademarks of Hiossen Inc. All other brand or product names are trademarks of their respective companies or organizations.

PB24MUALTR1.1