



# HIOSSEN SURGICAL KITS

PRODUCT CATALOG

**Hiossen**  
IMPLANT



## **Hiossen Surgical Kit**

Version: PC26HSKLTR1.0



105 Challenger Road, 7th Floor, Ridgefield Park, NJ 07660  
[www.hiossen.com](http://www.hiossen.com) Email: [marketing@hiossen.com](mailto:marketing@hiossen.com)

📷 @Hiossen    📺 @HiossenImplants    📺 @Hiossen

# Hiossen Surgical Kits

## GUIDED SURGERY KITS

ET/EK OneGuide Kit	004
OneCAS Kit	016
OneEM Kit	022
One485 Kit	028

## IMPLANT GURGERY KITS

122 Taper Kit	032
Taper Kit & Taper Ultra Kit	040
485 Kit	052
Assist Kit	056
Ultra Kit	060

## PROSTHETIC & MAINTENANCE KITS

Prosthetic Simple Kit & Prosthetic Kit	068
ESR Kit & ESR Full Kit	076
EIR Kit & EIR Full Kit	082
IM-Cure Kit	086

## SINUS SURGERY KITS

CAS Kit	090
LAS Kit	094

## NARROW RIDGE SURGERY KITS

ESSET Kit	100
EM Kit	

## SURGICAL INSTRUMENTS

Surgical Instruments Kit	104
--------------------------	-----

# ET OneGuide Kit (HOGK)

For

ETIII/IV

SSII/III

Top panel components



**Torque Wrench**

TQWCB



**Depth Gauge**

ODG





# EK OneGuide Kit (HKOGK)

For

EKIII

SSII/III

Top panel components


**Torque Wrench**

TQWCB





**Depth Gauge**


ODG





# OneGuide Kit Surgical Kit Instruments

OneGuide Template	
Description	Image
<ul style="list-style-type: none"> <li>Two sizes of guide holes depending on the diameter of the implant               <ul style="list-style-type: none"> <li>D5.0 for implant diameters: F3.5/4.0/4.5</li> <li>D5.8 for implant diameter: F5.0</li> </ul> </li> <li>Dual contact feature ensures excellent accuracy in positioning and stability</li> <li>Simple drilling sequence by adopting 122 Taper Kit concepts</li> <li>Packing unit: surgical guide (optional: SmartFit abutment, Temporary Crown)</li> </ul>	



Tissue Punch				
Description	D/Ø	Regular Hole (Ø5.1)	Wide Hole (Ø5.8)	Image/Guide
<ul style="list-style-type: none"> <li>Used to remove gingiva flaplessly</li> <li>7 types according to the OneGuide guide holes</li> <li>Drills except two types are not included in the Kit (HOGTP35R, HOGTP40WC) and are sold separately</li> <li>Recommended speed: 800 ~ 1,200 rpm</li> </ul>	Ø3.0	HOGTP30R	-	
	Ø3.5	HOGTP35R	-	
	Ø4.0	HOGTP40R	HOGTP40WC	
	Ø4.5	HOGTP45R	HOGTP45WC	
	Ø5.0	-	HOGTP50WC	


Flattening Drill				
Description	D/Ø	Regular Hole (Ø5.1)	Wide Hole (Ø5.8)	Image/Guide
<ul style="list-style-type: none"> <li>Used for flattening narrow or irregular bone ridges, before initial drill</li> <li>Multiple cutting edge designed to prevent drill bouncing</li> <li>2 Types (Below F5.0/for F5.0)</li> <li>Recommended speed: 800 ~ 1,200 rpm</li> </ul>	Ø4.5	HOGFD45	-	
	Ø5.2	-	HOGFD52WC	

Initial Drill				
Description	D/Ø	Regular Hole (Ø5.1)	Wide Hole (Ø5.8)	Image/Guide
<ul style="list-style-type: none"> <li>Used after Tissue Punch for initial drilling</li> <li>Secures depth for subsequent drills for more stability</li> <li>Available in 3 types: (F3.5 soft bone/below F5.0/for F5.0)</li> <li>Recommended speed: 800 ~ 1,200 rpm</li> </ul>	Ø3.5 (Soft)	HOGID2206II	-	
	Ø4.0/Ø4.5	HOGID2206I	-	
	Ø5.0 (W)	-	HOGID2206IWC	


Initial Drill (Short Type)				
Description	D/Ø	Regular Hole (Ø5.1)	Wide Hole (Ø5.8)	Image/Guide
<ul style="list-style-type: none"> <li>Short type Initial Drill (5.3mm shorter)</li> <li>Used for limited intermaxillary space</li> <li>Available in 3 types: (F3.5 soft bone/below F5.0/for F5.0)</li> <li>Recommended speed: 800 ~ 1,200 rpm</li> </ul>	Ø3.5 (Soft)	HOGD2206IIS	-	
	Ø4.0/Ø4.5	HOGD2206IS	-	
	Ø5.0 (W)	-	HOGD2206ISWC	


# OneGuide Kit Surgical Kit Instruments


OneGuide Drill							
Description/Item code							
<ul style="list-style-type: none"><li>• Taper Drill optimized for III/IV type implant</li><li>• Used for placing F3.5 ~ F5.0 &amp; 6 ~ 13mm implants</li><li>• Multi-step drill design allows for stable drilling</li><li>• Drills for 6mm and F5.5(W) types are sold separately</li><li>• Recommended speed: Soft Bone (800 ~ 1,200 rpm) / Normal, Hard Bone (1,200 ~ 1,500 rpm)</li></ul>							
Below F5.0 Regular Hole (Ø5.1)		D/Ø	Ø3.5	Ø4.0	Ø4.5	Image/Guide 	
		Y-Dim.	0.7	0.9	1.0		
L	TL	GD	5.0	5.0	5.0		
6	36.1		HOGTPD3506	HOGTPD4006	HOGTPD4506		
7	36.1		HOGTPD3507	HOGTPD4007	HOGTPD4507		
8.5	36.1		HOGTPD3508	HOGTPD4008	HOGTPD4508		
10	36.1		HOGTPD3510	HOGTPD4010	HOGTPD4510		
11.5	37.6		HOGTPD3511	HOGTPD4011	HOGTPD4511		
13	39.1		HOGTPD3513	HOGTPD4013	HOGTPD4513		
Below F5.0 Wide Hole (Ø5.8)		D/Ø	Ø3.5 (w)	Ø4.5 (w)	Ø5.0 (w)	Ø5.5 (w)	Image/Guide 
		Y-Dim.	0.7	0.9	1.0	1.0	
L	TL	GD	5.7	5.7	5.7	5.7	
6	36.1		HOGTPD3506WC	HOGTPD4506WC	HOGTPD5006WC	HOGTPD5506WC	
7	36.1		HOGTPD3507WC	HOGTPD4507WC	HOGTPD5007WC	HOGTPD5507WC	
8.5	36.1		HOGTPD3508WC	HOGTPD4508WC	HOGTPD5008WC	HOGTPD5508WC	
10	36.1		HOGTPD3510WC	HOGTPD4510WC	HOGTPD5010WC	HOGTPD5510WC	
11.5	37.6		HOGTPD3511WC	HOGTPD4511WC	HOGTPD5011WC	HOGTPD5511WC	
13	39.1		HOGTPD3513WC	HOGTPD4513WC	HOGTPD5013WC	HOGTPD5513WC	




Twist Drill					
<div>Description</div> <div><ul style="list-style-type: none"><li>Used to place F3.5 implants in soft bone</li><li>Available in 5 types</li><li>Recommended speed: 800 ~ 1,200 rpm</li></ul></div>	F3.5 Soft Bone		D/Ø	Ø2.2	<div>Image/Guide</div> <div></div>
	Regular Hole (Ø5.1)		Y-Dim.	0.7	
	L	TL	GD	5.0	
	7	36.1		HOGTD2207	
	8.5	36.1		HOGTD2208	
	10	36.1		HOGTD2210	
	11.5	37.6		HOGTD2211	
13	39.1		HOGTD2213		

# OneGuide Kit Surgical Kit Instruments




OneGuide Vertical Twist Drill				
Description	D/Ø	Regular Hole (Ø5.1)	Wide Hole (Ø5.8)	Image/Guide
<ul style="list-style-type: none"> <li>Used for drilling before OneGuide Anchor</li> <li>Sold as an individual item</li> <li>Recommended speed: 800 ~1,200 rpm</li> </ul>	Ø1.5	HOGTD1506	HOGTD1506W	




OneGuide Vertical Bone Anchor				
Description	D/Ø	Regular Hole (Ø5.1)	Wide Hole (Ø5.8)	Image/Guide
<ul style="list-style-type: none"> <li>Used for fixing the OneGuide in place (e. g. edentulous case)</li> <li>Mounted on the alveolar bone vertically to fix OneGuide in place</li> <li>Soft bone: placed directly</li> <li>Normal/hard bone: placed after using the OneGuide Vertical Drill</li> <li>Tighten at 20 rpm with Anchor Driver</li> <li>Sold as an individual item</li> </ul>	Ø2.0	HOGBAR	HOGBAW	



OneGuide Vertical Implant Anchor				
Description	D/Ø	Regular Hole (Ø5.1)	Wide Hole (Ø5.8)	Image/Guide
<ul style="list-style-type: none"> <li>Used for fixing OneGuide in place (e. g. edentulous case)</li> <li>Mounted on the implant vertically to fix OneGuide in place</li> <li>Tighten with 1.2 hex Hand Driver</li> <li>Only used for Regular connection of F4.0 or greater</li> <li>Sold as an individual item</li> </ul>	Ø2.0	HOGFAR	HOGFAW	

OneGuide NoMount Driver for ET				
Description	*C	D/Ø	Item code	Image/Guide
<ul style="list-style-type: none"> <li>Used to place a NoMount implant</li> <li>It is recommended to place the implant ~80% of the planned implant depth with this driver</li> <li>*C = Connection</li> </ul>	F3.5	Mini Regular Hole (Ø5.1)	HOGNMDM50	
	F4.0/4.5	Regular Regular Hole (Ø5.1)	HOGNMDR50	
	F5.0	Regular Wide Hole (Ø5.8)	HOGNMDR57	




# OneGuide Kit Surgical Kit Instruments




OneGuide Implant Driver for ET				
Description	*C	D/Ø	Item code	Image/Guide
<ul style="list-style-type: none"> <li>Used with a wrench for finishing the final implant placement</li> <li>Yellow groove aligns the abutment hex direction</li> <li>Match the grooves on the OneGuide template with the grooves on the driver</li> <li>*C = Connection</li> </ul>	F3.5	<b>Mini</b> Regular Hole (Ø5.1)	HOGFDM50	
	F4.0/4.5	<b>Regular</b> Regular Hole (Ø5.1)	HOGFDR50	
	F5.0 (W)	<b>Regular</b> Wide Hole (Ø5.8)	HOGFDR57	




Implant Driver (Stopper Type)				
Description	*C	D/Ø	Item code	Image/Guide
<ul style="list-style-type: none"> <li>Featuring stopper design to prevent entry below the upper surface of OneGuide hole</li> <li><b>Sold as an individual item</b></li> <li>*C = Connection</li> </ul>	F3.5	<b>Mini</b> Regular Hole (Ø5.1)	HOGFDSM50	
	F4.0/4.5	<b>Regular</b> Regular Hole (Ø5.1)	HOGFDSR50	
	F5.0 (W)	<b>Regular</b> Wide Hole (Ø5.8)	HOGFDSR57	




OneGuide Taper Cortical Drill			
Description	L	Regular Hole (Ø5.1)	Image/Guide
<ul style="list-style-type: none"> <li>Used for placing F4.5 and F5.0 implants in hard bone</li> <li>Optimize placement by cutting cortical bone</li> <li>Drills for 13 mm diameter implant is sold separately</li> <li>Drilling up to the first black line for 6mm placement</li> <li>Recommended speed: 800 ~ 1,200 rpm</li> </ul>	6 / 7 / 8.5mm	HOGTCD4507	
	10 / 11.5mm	HOGTCD4510	
	13mm	HOGTCD4513	
	L	Wide Hole (Ø5.8)	
	6 / 7 / 8.5mm	HOGTCD5007WC	
	10 / 11.5mm	HOGTCD5010WC	
	13mm	HOGTCD5013WC	

# OneGuide Kit Surgical Kit Instruments



OneGuide No Mount Driver for EK				
Description	*C	D/Ø	Item code	Image/Guide
<ul style="list-style-type: none"> <li>Used for placing EK NoMount Implants</li> <li>Recommended to place ~80% of the planned implant depth with this driver</li> <li><b>Sold as an individual item</b></li> <li>*C = Connection</li> </ul>	F3.5	Regular Regular Hole (Ø5.1)	HOGNMDM50K	
	F4.0/4.5		HOGNMDR50K	
	F5.0	Regular Wide Hole (Ø5.8)	HOGNMDR57K	


OneGuide Implant Driver for EK				
Description	*C	D/Ø	Item code	Image/Guide
<ul style="list-style-type: none"> <li>Used with a wrench for finishing the final implant placement</li> <li>Yellow groove aligns the abutment hex direction</li> <li>Match the grooves on the OneGuide template with the grooves on the driver</li> <li>Sold as an individual item</li> <li>*C = Connection</li> </ul>	F3.5	Regular Regular Hole (Ø5.1)	HOGFDM50K	
	F4.0/4.5		HOGFDR50K	
	F5.0	Regular Wide Hole (Ø5.8)	HOGFDR57K	


OneGuide NoMount Driver for SS				
Description	*P	D/Ø	Item code	Image
<ul style="list-style-type: none"> <li>Used for placing SS NoMount Implants</li> <li>Recommended to place ~80% of the planned implant depth with this driver</li> <li>*P = Platform</li> </ul>	F3.5/F4.0/4.5	Mini Regular Hole (Ø5.1)	HOGNMDR50S	
	F5.0	Regular Wide Hole (Ø5.8)	HOGNMDR57S	
	F5.0	Wide Extra Wide Hole (Ø6.8)	HOGNMDW67S	


OneGuide Implant Driver for SS				
Description	*P	D/Ø	Item code	Image
<ul style="list-style-type: none"> <li>Used with a wrench for finishing the final implant placement</li> <li>Yellow groove aligns the abutment hex direction</li> <li>Match the grooves on the OneGuide template with the grooves on the driver</li> <li><b>Sold as an individual item</b></li> <li>*P = Platform</li> </ul>	F3.5/F4.0/4.5	Mini Regular Hole (Ø5.1)	HOGFDR50S	
	F5.0	Regular Wide Hole (Ø5.8)	HOGFDR57S	
	F5.0	Wide Extra Wide Hole (Ø6.8)	HOGFDW67S	

# OneGuide Kit Surgical Kit Instruments

OneGuide Path Drill					
Description	L		Regular Hole (Ø5.1)	Wide Hole (Ø5.8)	Image/Guide
<ul style="list-style-type: none"> <li>Drill for correcting path deviation during OneGuide surgery</li> <li>Used for creating implant placement path for extraction cases</li> <li>Flat blade design optimized for cutting inclined bone</li> <li>4 types for each OneGuide hole diameter, 8 types in total: Regular hole (Ø5.1) / Wide hole (Ø5.8)</li> <li>For 13mm drills, depth is adjusted according to the black lines</li> <li>Recommended speed: 1,200 ~ 1,500 rpm</li> </ul>	7	Ø2.5	HOGSD2507	HOGSD2507WC	
		Ø3.0	HOGSD3007	HOGSD3007WC	
	13	Ø2.5	HOGSD2513	HOGSD2513WC	
		Ø3.0	HOGSD3013	HOGSD3013WC	

Anchor Screw		
Description	Item Code	Image
<ul style="list-style-type: none"> <li>Used to affix the OneGuide firmly</li> <li>Selectable at the preoperative planning stage</li> </ul>	HQGAS18	

Anchor Drill		
Description	Item Code	Image
<ul style="list-style-type: none"> <li>Used for drilling before using anchor screw</li> <li>Recommended speed: 800 ~ 1,200 rpm</li> </ul>	HQGATD13	

OneGuide Anchor Driver (Mount Driver)		
Description	Item Code	Image
<ul style="list-style-type: none"> <li>Used to place anchor screw</li> </ul>	HASMDS	




# Drilling Sequence OneGuide Drill

**EKIII | ETIII**


(Length: 10mm)

**Ø3.5**




Bone Quality	Tissue Punch	Flattening Drill	Initial Drill	Drill (Ø2.2)	Drill (F3.5)	Drill (F4.0)	Nomount Driver	Implant Driver
Soft	►	(►)	(F3.5 Soft)►	►				
Normal	►	(►)	►		►		Implant Placement (Up to 80%)	Implant Placement
Hard	►	(►)	►		►	►		

**Ø4.0**



Bone Quality	Tissue Punch	Flattening Drill	Initial Drill	Drill (F3.5)	Drill (F4.0)	Drill (F4.5)	Nomount Driver	Implant Driver
Soft	►	(►)	►	►				
Normal	►	(►)	►	►	►		Implant Placement (Up to 80%)	Implant Placement
Hard	►	(►)	►	►		►		

**Ø4.5**



Bone Quality	Tissue Punch	Flattening Drill	Initial Drill	Drill (F3.5)	Drill (F4.0)	Drill (F4.5)	Cortical (F4.5)	Nomount Driver	Implant Driver
Soft	►	(►)	►	►	►				
Normal	►	(►)	►	►		►		Implant Placement (Up to 80%)	Implant Placement
Hard	►	(►)	►	►		►	►		




# Drilling Sequence OneGuide Drill

**EKIII | ETIII**

(Length: 10mm)


**Ø5.0**



Bone Quality	Tissue Punch (W)	Flattening Drill (W)	Initial Drill (W)	Drill (W) (F3.5)	Drill (W) (F4.5)	Drill (W) (F5.0)	Cortical (W) (F5.0)	Nomount Driver	Implant Driver
Soft	►	(►)	►	►	►				
Normal	►	(►)	►	►		►		Implant Placement (Up to 80%)	Implant Placement
Hard	►	(►)	►	►		►	►		


※ For extra short implants (Ø4.0, Ø4.5 / 6mm) only

**Ø4.0**



Bone Quality	Tissue Punch	Flattening Drill	Initial Drill	Drill (F3.5)	Drill (F4.0)	Drill (F4.5)	Nomount Driver	Implant Driver
Soft	►	(►)	►	F3.5x6				
Normal	►	(►)	►	F3.5x6	F4.0x6		Implant Placement (Up to 80%)	Implant Placement
Hard	►	(►)	►	F3.5x6		F4.5x6		

**Ø4.5**




Bone Quality	Tissue Punch	Flattening Drill	Initial Drill	Drill (F3.5)	Drill (F4.0)	Drill (F4.5)	Cortical (F4.5)	Nomount Driver	Implant Driver
Soft	►	(►)	►	F3.5x6	F4.0x6				
Normal	►	(►)	►	F3.5x6		F4.5x6		Implant Placement (Up to 80%)	Implant Placement
Hard	►	(►)	►	F3.5x6		F4.5x6	F4.5x6~8.5		

# ETIV


(Length: 10mm)

Ø4.0




Bone Quality	Tissue Punch	Flattening Drill	Initial Drill	Drill (F3.5)	Drill (F4.0)	Nomount Driver	Implant Driver
Soft	►	(►)	►	►		Implant Placement (Up to 80%)	Implant Placement
Normal	►	(►)	►	►	►		

Ø4.5



Bone Quality	Tissue Punch	Flattening Drill	Initial Drill	Drill (F3.5)	Drill (F4.0)	Drill (F4.5)	Nomount Driver	Implant Driver
Soft	►	(►)	►	►	►		Implant Placement (Up to 80%)	Implant Placement
Normal	►	(►)	►	►		►		

Ø5.0



Bone Quality	Tissue Punch (W)	Flattening Drill (W)	Initial Drill (W)	Drill (W) (F3.5)	Drill (W) (F4.5)	Drill (W) (F5.0)	Nomount Driver	Implant Driver
Soft	►	(►)	►	►	►		Implant Placement (Up to 80%)	Implant Placement
Normal	►	(►)	►	►		►		



# OneCAS Kit (HOCK)

For

EKIII

ETIII/IV

**Hydraulic Membrane Lifter Tube**  
SNMT



**Bone Condenser**  
SNBC1114



Top panel components

**Depth Gauge**  
HCDG



**Depth Gauge(W)**  
HCDGW




Lower panel components

**Bone Carrier Head**  
SNBCH30



**Bone Carrier**  
SNBCS35






# OneCAS Kit Surgical Kit Instruments



## OneCAS Stopper




- Stopper number indicates the length when engaged
- On the Kit, protruding length marked in blue for 7mm drills and in red for 10mm drills
- Color coded by length
- Recommended number of use: 50 times


Diameter	1	2	3	4	5	6	7	8	9
									
<b>F4.0/4.5</b>	HOCDS01	HOCDS02	HOCDS03	HOCDS04	HOCDS05	HOCDS06	HOCDS07	HOCDS08	HOCDS09
<b>F5.0 (W)</b>	HOCDS01W	HOCDS02W	HOCDS03W	HOCDS04W	HOCDS05W	HOCDS06W	HOCDS07W	HOCDS08W	HOCDS09W
<b>Color</b>	Yellow	Purple	Blue	Yellow	Purple	Blue	Yellow	Purple	Blue

## Depth Gauge

Description		F4.0/4.5	F5.0 (W)	Image/Guide
	GD	5.0	5.7	
<ul style="list-style-type: none"> <li>• Checks for internal sinus lift</li> <li>• Measures residual bone depth</li> <li>• Use with a stopper for safer approach</li> <li>• Same depth marking lines as a 10mm drill</li> </ul>		HCDG	HCDGW	

Hydraulic Membrane Lifter			
Description	D/Ø	Item code	Image/Guide
<ul style="list-style-type: none"> <li>Dedicated surgical instrument for the OneCAS KIT</li> <li>Used by inserting into the OneGuide hole after marking drilling</li> <li>Thin-wall design optimizes sealing during flapless procedures</li> </ul>	<b>Regular</b> Hole (Ø5.1)	OCHML	
	<b>Wide</b> Hole (Ø5.8)	OCHMLW	


Bone Carrier Head			
Description	D/Ø	Item code	Image/Guide
<ul style="list-style-type: none"> <li>Dedicated bone delivery instrument for the OneCAS KIT</li> <li>Inserted and fixed up to the end of the OneGuide hole during use</li> <li>OCBCH30: Used after OneCAS drill Ø3.1</li> <li>OCBCH32: Used after OneCAS drill Ø3.3 / Ø3.6</li> <li>CBCH37W: Used after OneCAS drill Ø3.8 / Ø4.1</li> <li>After filling bone material up to the head marking line, use with a bone condenser to completely deliver bone into the sinus cavity, then repeat as needed</li> </ul>	<b>Regular</b> Hole (Ø5.1)	OCBCH30	
		OCBCH32	
	<b>Wide</b> Hole (Ø5.8)	OCBCH37W	

Bone Carrier		
Description	Item Code	Image
<ul style="list-style-type: none"> <li>Dedicated bone delivery instrument for the OneCAS KIT</li> <li>The head is secured by tightening the rear handle of the body</li> <li>The head can be replaced and used interchangeably</li> </ul>	OCBCS30	




## Drilling Sequence OneCAS Kit

Residual Bone **3mm**  
**Ø4.0**




Bone Quality	Tissue Punch	Flattening Drill	Twist Drill (Ø2.2x7)	OneCAS Drill	OneCAS Drill	Depth Gauge	OneCAS Drill	Depth Gauge	Hydraulic Lifter	Bone Carrier	Implant
Soft	►	►	►	Ø2.8x7	Ø2.8x7	►	Ø2.8x7	►	►	►	►
Normal	►	►	►	Ø3.1x7	Ø3.1x7	►	Ø3.1x7	►	►	►	►
Stopper			5	5	4	7	3	6			

Residual Bone **3mm**  
**Ø4.5**



Bone Quality	Tissue Punch	Flattening Drill	Twist Drill (Ø2.2x7)	OneCAS Drill	OneCAS Drill	OneCAS Drill	Depth Gauge	OneCAS Drill	Depth Gauge	Hydraulic Lifter	Bone Carrier	Implant
Soft	►	►	►	Ø2.8x7	Ø3.3x7	Ø3.3x7	►	Ø3.3x7	►	►	►	►
Normal	►	►	►	Ø3.1x7	Ø3.6x7	Ø3.6x7	►	Ø3.6x7	►	►	►	►
Stopper			5	5	5	4	7	3	6			

Residual Bone **3mm**  
**Ø5.0**




Bone Quality	Tissue Punch (W)	Flattening Drill (W)	Twist Drill (Ø2.2x7)(W)	OneCAS Drill (W)	OneCAS Drill (W)	OneCAS Drill (W)	Depth Gauge (W)	OneCAS Drill (W)	Depth Gauge (W)	Hydraulic Lifter (W)	Bone Carrier (W)	Implant
Soft	►	►	►	Ø2.8x7	Ø3.8x7	Ø3.8x7	►	Ø3.8x7	►	►	►	►
Normal	►	►	►	Ø3.1x7	Ø4.1x7	Ø4.1x7	►	Ø4.1x7	►	►	►	►
Stopper (W)			5	5	5	4	7	3	6			




# Drilling Sequence OneCAS Kit

Residual Bone 8mm  
Ø4.0




Bone Quality	Tissue Punch	Flattening Drill	Twist Drill (Ø2.2x7)	OneCAS Drill	OneCAS Drill	Depth Gauge	OneCAS Drill	Depth Gauge	Hydraulic Lifter	Bone Carrier	Implant
Soft	▶	▶	▶	Ø2.8x7	Ø2.8x10	▶	Ø2.8x10	▶	▶	▶	▶
Normal	▶	▶	▶	Ø3.1x7	Ø3.1x10	▶	Ø3.1x10	▶	▶	▶	▶
Stopper					2	2	1	1			

Residual Bone 8mm  
Ø4.5



Bone Quality	Tissue Punch	Flattening Drill	Twist Drill (Ø2.2x7)	OneCAS Drill	OneCAS Drill	OneCAS Drill	Depth Gauge	OneCAS Drill	Depth Gauge	Hydraulic Lifter	Bone Carrier	Implant
Soft	▶	▶	▶	Ø2.8x7	Ø3.3x7	Ø3.3x10	▶	Ø3.3x10	▶	▶	▶	▶
Normal	▶	▶	▶	Ø3.1x7	Ø3.6x7	Ø3.6x10	▶	Ø3.6x10	▶	▶	▶	▶
Stopper						2	2	1	1			

Residual Bone 8mm  
Ø5.0



Bone Quality	Tissue Punch(W)	Flattening Drill(W)	Twist Drill (Ø2.2x7)(W)	OneCAS Drill(W)	OneCAS Drill(W)	OneCAS Drill(W)	Depth Gauge(W)	OneCAS Drill(W)	Depth Gauge(W)	Hydraulic Lifter(W)	Bone Carrier(W)	Implant
Soft	▶	▶	▶	Ø2.8x7	Ø3.8x7	Ø3.8x10	▶	Ø3.8x10	▶	▶	▶	▶
Normal	▶	▶	▶	Ø3.1x7	Ø4.1x7	Ø4.1x10	▶	Ø4.1x10	▶	▶	▶	▶
Stopper (W)						2	2	1	1			

# OneEM Kit (HOMSK)

### Top panel components

**Torque Wrench**  
TQWCB



## Depth Gauge




**For**


EKIII


ETIII/IV


EM(MS)




# OneEM Kit Surgical Kit Instruments

OneGuide Template	
Description/Item code	Image/Guide
<ul style="list-style-type: none"> <li>Two sizes of guide holes depending on the diameter of the implant               <ul style="list-style-type: none"> <li>- D5.0 for implant diameters: F3.5/4.0/4.5</li> <li>- D5.8 for implant diameter: F5.0</li> </ul> </li> <li>Dual contact feature ensures excellent accuracy in positioning and stability</li> <li>Simple drilling sequence by adopting 122 Taper Kit concepts</li> <li>Packing unit: surgical guide               <ul style="list-style-type: none"> <li>(optional: SmartFit abutment, Temporary Crown)</li> </ul> </li> </ul>	

OneMS Drill							
Description/Item code							
<ul style="list-style-type: none"><li>Straight type drills for EM implant / ETIII Ø3.2 implants (EM Ø2.0 ~ Ø3.0 / ETIII Ø3.2 implants)</li><li>OneMS Cortical Drill is used for placing ETIII Ø3.2 implants in hard bones</li><li>Start with 8.5mm drill for stable drilling</li><li>Recommended speed: 800 ~ 1,200 rpm</li></ul>							
L	TL	D/Ø	Ø1.5	Ø1.8	Ø2.3	Ø2.7	
		GD	0.6	0.6	0.6	0.6	
8.5	37.5		HOMD1508	HOMD1808	HOMD2308	HOMD2708	
10	39		HOMD1510	HOMD1810	HOMD2310	HOMD2710	
11.5	40.5		HOMD1511	HOMD1811	HOMD2311	HOMD2711	
13	42		HOMD1513	HOMD1813	HOMD2313	HOMD2713	
Color			Yellow	Green	Blue	Pink	

OneMS Cortical Drill			
Description/Item code		Image	
<ul style="list-style-type: none"> <li>Used to remove cortical bone in hard bone</li> <li>Used for expanding the cortical bone after using the Straight Drill (ETIII Ø3.0 only)</li> <li>Recommended speed: 800 ~ 1,200 rpm</li> </ul>			
ETIII Ø3.0	EKIII Ø3.3		
HOGCD30	HOGCD30K		

OneMS Lance Drill			
Description/Item code		Image/Guide	
<ul style="list-style-type: none"> <li>Creates holes in the bone to facilitate initial drilling</li> <li>Bone density can be determined through drilling</li> <li>Single item (excluded from OneMS Kit)</li> </ul>			
L	TL	D/Ø	Ø1.5
8.5	37.5		HOMLD1508
10	39		HOMLD1510
11.5	40.5		HOMLD1511
13	42		HOMLD1513

# OneEM Kit Surgical Kit Instruments

Flattening Drill	
Description/Item code	Image
<ul style="list-style-type: none"> <li>Used for narrow or uneven ridges</li> <li>Cutting edge design allows for stable bone removal</li> </ul>	
HOGFD30	


Driver Separator	
Description/Item code	Image
<ul style="list-style-type: none"> <li>Used in case a driver is wedged with the implant during insertion</li> <li>Insert the Driver Separator into the driver groove to loosen it from the implant</li> </ul>	
HMSDS	


Tissue Punch	
Description/Item code	Image
<ul style="list-style-type: none"> <li>Used to remove gingiva to continue with flapless surgery</li> </ul>	
HOGTP30N	

OneMS Driver	
Description	Image
<ul style="list-style-type: none"> <li>Used with the torque wrench for adjustment of the final implant placement for MS implant narrow ridge</li> <li>Match the mark with the implant</li> </ul>	
<b>F2.0/2.5/3.0</b>	
NoMount Driver	HOGMDMSN
Implant Driver	HOGFDMSN

Implant Driver	
Description	Image
<ul style="list-style-type: none"> <li>Used with the torque wrench for adjustment of the final implant placement</li> <li>Yellow grooves aligns with abutment hex direction</li> <li>Match the groove of OneGuide template with the groove of driver</li> <li>*C = Connection</li> </ul>	
ET Ø3.2	EK Ø3.3
HOGFDM35	OGFDM35K

NoMount Driver	
Description	Image
<ul style="list-style-type: none"> <li>Used with the torque wrench for adjustment of the final implant placement</li> <li>Yellow grooves aligns with abutment hex direction</li> <li>Match the groove of OneGuide template with the groove of driver</li> <li>*C = Connection</li> </ul>	
ET Ø3.2	EK Ø3.3
HOGNMDM35	OGNMDM35K

Adapter	
Description/Item code	Image
<ul style="list-style-type: none"> <li>Use as a torque driver to connecting to the engine</li> </ul>	
TEADT	

OneMS Path Drill				
Description				
<ul style="list-style-type: none"><li>Used for correcting path deviation</li><li>Can establish implant placement path during surgery</li><li>For inclined bone cutting with a flat blade design</li><li>For the 13mm specification, depth adjustment follows marking lines: upper line (11.5mm), middle line (10mm), lower line (8.5mm)</li><li>Recommended speed: 1,200~1,500 rpm</li></ul>				
L	D/Ø	Ø1.5	Ø2.0	
7.0		HOMSD1507	HOMSD2007	
13.0		HOMSD1513	HOMSD2013	

# Drilling Sequence OneEM Drill

**EKIII | ETIII | EM**

(Length: 10mm)


For Implant 10 / 11.5 / 13 mm sequences, an Ø8.5 mm drill must be used at each step

Ex. Ø2.5 × 11.5mm MS implant

:Tissue punch ► Flattening drill ► Ø1.8 × 8.5mm (mandatory pilot drilling) ► Ø1.8 × 11.5mm ► NoMount driver ► Implant driver


※ Ø1.5 drill (optional): may be selectively used for initial drilling when higher precision is required

## EKIII Ø3.0




Bone Quality	Tissue Punch	Flattening Drill	OneMS Drill (Ø1.5)	OneMS Drill (Ø1.8)	OneMS Drill (Ø2.3)	OneMS Drill (Ø2.7)	F3.0 Cortical Drill	NoMount Driver	Implant Driver
Soft	►	(►)	(►)	►	►			►	►
Normal	►	(►)	(►)	►	►	►		►	►
Hard	►	(►)	(►)	►	►	►	►	►	►

## ETIII Ø3.0




Bone Quality	Tissue Punch	Flattening Drill	OneMS Drill (Ø1.5)	OneMS Drill (Ø1.8)	OneMS Drill (Ø2.3)	OneMS Drill (Ø2.7)	F3.0 Cortical Drill	NoMount Driver	Implant Driver
Soft	►	(►)	(►)	►	►			►	►
Normal	►	(►)	(►)	►	►	►		►	►
Hard	►	(►)	(►)	►	►	►	►	►	►

## EM Ø2.0




Bone Quality	Tissue Punch	Flattening Drill	OneMS Drill (Ø1.5)	NoMount Driver	Implant Driver		
					G/H 2.5	G/H 4.0	Denture
Soft	►	(►)	►	►		►	
Normal	►	(►)	►	►		►	
Hard	►	(►)	►	►		►	

## EM Ø2.5



Bone Quality	Tissue Punch	Flattening Drill	OneMS Drill (Ø1.5)	OneMS Drill (Ø1.8)	OneMS Drill (Ø2.3)	NoMount Driver	Implant Driver		
							G/H 2.5	G/H 4.0	Denture
Soft	►	(►)	(►)	►	-	►		►	
Normal	►	(►)	(►)	►	-	►		►	
Hard	►	(►)	(►)	►	(►)	►		►	

## EM Ø3.0



Bone Quality	Tissue Punch	Flattening Drill	OneMS Drill (Ø1.5)	OneMS Drill (Ø1.8)	OneMS Drill (Ø2.3)	OneMS Drill (Ø2.7)	NoMount Driver	Implant Driver		
								G/H 2.5	G/H 4.0	Denture
Soft	►	(►)	(►)	►	-	-	►		►	
Normal	►	(►)	(►)	►	-	-	►		►	
Hard	►	(►)	(►)	►	►	(►)	►		►	



# One485 Kit (HO485K)

For

EKIII


ETIII/IV

SSIII

**Torque Wrench**  
TW30B



**Depth Gauge**  
OSDG





# One485 Kit Surgical Kit Instruments

One485 Twist Drill				Image/Guide	
Description/Item code				Image/Guide	
<ul style="list-style-type: none"> <li>Ensures initial drilling stability and secures the exact depth of the guide hole</li> <li>Straight design allows drilling to 1 mm shorter than the planned implant length</li> <li>Available in 8 configurations according to OneGuide hole diameter <ul style="list-style-type: none"> <li>Regular hole (Ø5.1) / Wide hole (Ø5.8)</li> </ul> </li> <li>Recommended drilling speed: 800–1,200 rpm</li> </ul>				<p>Regular Hole (Ø5.1)</p> <p>Wide Hole (Ø5.8)</p>	
L	TL	Regular Hole(Ø5.1) F4.0 / F4.5	Wide Hole(Ø5.8) F5.0W		
4.0	30.2	-	O4D2204WC		
5.0	31.2	-	O4D2205WC		
6.0	32.2	O4D2206	O4D2206WC		
7.0	33.2	O4D2207	O4D2207WC		
8.5	34.7	O4D2208	O4D2208WC		


One485 Pilot Drill				Image/Guide	
Description/Item code				Image/Guide	
<ul style="list-style-type: none"> <li>Expands the hole diameter during intermediate drilling</li> <li>Available in 3 configurations according to OneGuide hole diameter <ul style="list-style-type: none"> <li>Regular hole (Ø5.1) / Wide hole (Ø5.8)</li> </ul> </li> <li>Implant length 4–5 mm: use 4 mm,</li> <li>Implant length 6–8.5 mm: use 6 mm</li> <li>Recommended drilling speed: 800–1,200 rpm</li> </ul>				<p>Regular Hole (Ø5.1)</p> <p>Wide Hole (Ø5.8)</p>	
L	TL	Regular Hole(Ø5.1) F4.5	Wide Hole(Ø5.8) F5.0W		
4.0	30.9	-	O4D3504WC		
6.0	32.9	O4D3506	O4D3506WC		

One485 Drill							Image/Guide	
Description/Item code							Image/Guide	
<ul style="list-style-type: none"> <li>Used as the final drill to achieve final osteotomy diameter and optimize insertion torque</li> <li>Lateral cutting edges: tapered drill geometry</li> <li>Available in 19 configurations according to OneGuide hole diameter <ul style="list-style-type: none"> <li>Regular hole (Ø5.1) / Wide hole (Ø5.8)</li> </ul> </li> <li>For F4.5 and F5.0 implants in hard bone, use the F4.5 and F5.0 hard drills</li> </ul>								
L	TL	Regular Hole(Ø5.1)			Wide Hole(Ø5.8)			
		F4.0	F4.5	F4.5 Hard	F5.0W	F5.0W Hard		
4.0	30.9	-	-		O4D5004WC	O4HD5004WC		
5.0	31.9	-	-		O4D5005WC	O4HD5005WC		
6.0	32.9	O4D4006	O4D4506	O4HD4506	O4D5006WC	O4HD5006WC		
7.0	33.9	O4D4007	O4D4507	O4HD4507	O4D5007WC	O4HD5007WC		
8.5	35.4	O4D4008	O4D4508	O4HD4508	O4D5008WC	O4HD5008WC		

# One485 Kit


Implant insertion depth adjustment method (Same concept applied as ETIII and EKIII)

## Ex. Short Implant / NoMount / Implant Driver Features



## Drilling Sequence One485 Drill EKIII | ETIII | SSIII

Extra Short  
Ø5.0 X 6mm




Ø5.0X6mm (Extra Short)	Bone Quality	Tissue Punch (W)	Flattening Drill (W)	Twist Drill (W)	Pilot Drill (W)	One485 Drill (W)	One485 Hard Drill (W)	NoMount Driver (W)	Implant Driver (W)
Ø5.0x 6/7/8.5mm	Normal	▶	(▶)	Ø2.2x6/7/8.5	Ø3.5x6	F5.0x6/7/8.5		Implant Placement (Up to 80%)	Implant Placement
	Hard	▶	(▶)	Ø2.2x6/7/8.5	Ø3.5x6		F5.0x6/7/8.5		

# One485 Kit


Implant insertion depth adjustment method (SSIII)

## NoMount / Implant Driver Features




## Drilling Sequence **One485 Drill** **EKIII | ETIII | SSIII**

### G/H 1.8 & 2.8



Implant	Bone Quality	Tissue Punch	Flattening Drill	Twist Drill	Pilot Drill	One485 Drill	One485 Drill	NoMount Driver	Implant Driver	
									G/H 1.8	G/H 2.8
Diameter 6/7/8.5mm	Normal	▶	(▶)	▶	▶	▶		Implant Placement (Up to 80%)	Implant Placement	
	Hard	▶	(▶)	▶	▶		▶			

### G/H 1.8&2.8



Implant	Bone Quality	Tissue Punch (W)	Flattening Drill (W)	Twist Drill (W)	Pilot Drill (W)	One485 Drill (W)	One485 Hard Drill (W)	NoMount Driver (W)	Implant Driver (W)	
									G/H 1.8	G/H 2.8
Diameter 6/7/8.5mm	Normal	▶	(▶)	▶	▶	▶		Implant Placement (Up to 80%)	Implant Placement	
	Hard	▶	(▶)	▶	▶		▶			

# 122 Taper Kit (HK122TPK)

For

EKIII

ETIII/IV


SSIII

Top panel components


**Torque Wrench**  
TQWCB





**Depth Gauge**  
ODG





# 122 Taper Kit Surgical Kit Instruments







122 Taper Drill									
Description				Image/Guide					
<ul style="list-style-type: none"> <li>Taper drill for Type III implants</li> <li>Specification by diameter and length</li> <li>Color coding displays implant diameter</li> <li>One step large-diameter drill is used to remove cortical bone from the hard bone</li> <li>122 Taper Kit single item (excluded from Taper Kit)</li> </ul>									
L	TL	D/Ø Y-Dim.	F3.5	F4.0	F4.5	F5.0	F5.5	F6.0	F7.0
			0.7	0.9	1.0	1.0	1.0	1.0	1.0
4.0	29.5		H122TPD3504	H122TPD4004	H122TPD4504	H122TPD5004	H122TPD5504	-	-
5.0	29.5		H122TPD3505	H122TPD4005	H122TPD4505	H122TPD5005	H122TPD5505	-	-
6.0	30.5		H122TPD3506	H122TPD4006	H122TPD4506	H122TPD5006	H122TPD5506	H122TPD6006	H122TPD7006
7.0	31.5		H122TPD3507	H122TPD4007	H122TPD4507	H122TPD5007	H122TPD5507	H122TPD6007	H122TPD7007
8.5	33		H122TPD3508	H122TPD4008	H122TPD4508	H122TPD5008	H122TPD5508	H122TPD6008	H122TPD7008
10	34.5		H122TPD3510	H122TPD4010	H122TPD4510	H122TPD5010	H122TPD5510	H122TPD6010	H122TPD7010
11.5	34.5		H122TPD3511	H122TPD4011	H122TPD4511	H122TPD5011	H122TPD5511	H122TPD6011	H122TPD7011
13	36		H122TPD3513	H122TPD4013	H122TPD4513	H122TPD5013	H122TPD5513	H122TPD6013	H122TPD7013
15	38		H122TPD3515	H122TPD4015	H122TPD4515	H122TPD5015	H122TPD5515	-	-
Color			Yellow	Green	Blue	Red	Yellow	Green	Blue

122 Taper Parallel Pin		
Description	Item Code	Image
<ul style="list-style-type: none"> <li>Parallel pin for 122 Taper Drill</li> <li>Used for checking position and direction of bone preparation</li> <li>Lower part for 2.2 drill, upper part for guide drill</li> <li>122 Taper Kit single item (excluded from Taper Kit)</li> <li>Other components same as Taper Kit</li> </ul>	HAPP2227	







Cortical Drill for Ultra-Wide			
Description	D/Ø	Item Code	Image
<ul style="list-style-type: none"> <li>Drill is used to remove cortical bone from hard bone (for ultra-wide)</li> <li>Dedicated drill by implant diameter</li> <li>It is recommended to drill to the bottom line of the marking line</li> </ul>	F6.0	HCD4C60	
	F7.0	HCD4C70	

# 122 Taper Kit Surgical Kit Instruments


Implant Driver for EK				
Description	*C	Length	Item Code	Image
<ul style="list-style-type: none"> <li>Connects directly to an EK implant for final adjustments to the implant's depth</li> <li>*C = Connection</li> </ul>	Regular	Short	HKSFDS	
		Long	HKSFDL	


Implant Driver for ET				
Description	*C	Length	Item Code	Image
<ul style="list-style-type: none"> <li>Connects directly to an ET implant for final adjustments to the implant's depth</li> <li>*C = Connection</li> </ul>	Mini	Short (17)	HGSMFDS	
		Long (24)	GSMFDL	
		Ex. Long (34)	HGSMFDE	
	Regular	Short (19)	HGSRFDS	
		Long (26.6)	GSRFDL	
		Ex. Long (33.6)	HGSRFDE	

NoMount Driver for EK				
Description	*C	Length	Item Code	Image
<ul style="list-style-type: none"> <li>Ø3.5 implant is combined with the bottom of the lower marking</li> <li>Ø4.0, Ø4.5, Ø5.0, Ø6.0 and Ø7.0 implants are combined with the upper part of the lower marking</li> <li>The distance between the two laser marking is 0.5mm</li> <li>*C = Connection</li> </ul>	Regular	Short (27.6)	HKSNMDCRS	
		Long (32.6)	HKSNMDCRL	
		Extra Long (37.6)	HKSNMDCRE	 <p>Ø4.0, Ø4.5, Ø5.0, Ø6.0, Ø7.0    Ø3.0, 3.5 Implant    Implant</p>


NoMount Torque Driver for ET				
Description	*C	Length	Item Code	Image
<ul style="list-style-type: none"> <li>Directly connects to an ET Implant for placement with a Torque Wrench</li> <li>Ensure correct and complete seating before applying torque; loose connection may cause implant fracture</li> <li>*C = Connection</li> </ul>	Mini	Short (19)	HGSNMT32S	
		Long (26.6)	HGSNMT32L	
		Ex. Long (33.6)	HGSNMT32E	
	Regular	Short (19)	HGSNMT35S	
		Long (26.6)	HGSNMT35L	
		Ex. Long (33.6)	HGSNMT35E	

# 122 Taper Kit Surgical Kit Instruments

LaceDrill (Guide Drill)			
Description	L	Item Code	Image
<ul style="list-style-type: none"> <li>Forming a hole to facilitate initial drilling</li> <li>Bone density determined through drilling</li> </ul>	<b>Long</b>	HAGDLC	

Twist Drill (Stopper Drill)										
Description				Image/Guide						
<ul style="list-style-type: none"> <li>Included in New Hanaro KIT</li> <li>Long stopper (6mm)</li> <li>Enabling a procedure without drill extension for posterior region</li> <li>The color coded stopper indicates the drill length</li> </ul>										

L	TL	D/Ø Y-Dim.	F2.2 0.6	F3.0 0.9	F3.3 1.0	F3.6 1.0	F3.8 1.0	F4.1 1.0	F4.3 1.0	F4.3 1.0
6	30.5		2D2206LC	3D3006LC	-	-	3D3806LC	-	-	-
7	31.5		2D2207LC01	3D3007LC01	-	-	3D3807LC01	-	-	-
8.5	33		2D2208LC01	3D3008LC01	-	-	3D3808LC01	-	-	-
10	34.5		2D2210LC01	3D3010LC01	-	-	3D3810LC01	-	-	-
11.5	34.5		2D2211LC01	3D3011LC01	3D3311LC01	3D3611LC01	3D3811LC01	3D4111LC01	3D4311LC01	3D4611LC01
13	36		2D2213LC01	3D3013LC01	-	-	3D3813LC01	-	-	-

Twist Drill (Non Stopper Drill)										
Description				Image/Guide						
<ul style="list-style-type: none"> <li>Included in New Hanaro KIT</li> <li>Used for limited access for the Stopper Drill into the oral cavity</li> <li>See the image provided in the Non-stopper Drill section for the sizes of the drill marking lines for short/long types</li> </ul>										

TL	D/Ø	F1.5	F2.0	F2.2	F2.7	F3.0	F3.3
33	Short	2D1518FNLC	2D2018FNLC	2D2218FNLC	3D2718FNLC	3D3018FNLC	3D3318FNLC
41		-	-	2D2215FNLC01	3D2715FNLC01	3D3015FNLC01	3D3315FNLC01
	D/Ø	F3.6	F3.8	F4.1	F4.3	F4.6	
33	Long	3D3618FNLC	3D3818FNLC	3D4118FNLC	3D4318FNLC	3D4618FNLC	
41		3D3615FNLC01	3D3815FNLC01	3D4115FNLC01	3D4315FNLC01	3D4615FNLC01	




# Drilling Sequence 122 Taper Drill

**EKIII | ETIII | SSIII**

(Length: 10mm)

Ø3.2mm for ET  
Ø3.3mm for EK




Bone Quality	Ø2.2 Drill	Ø2.7 Drill	F3.5 Taper Drill	Ø3.0 Implant
Soft	►	►		Implant Placement
Normal	►	► or ►		
Hard	►		►	

※ Lance drill can be used before Ø2.2 to mark the Osteotomy

※ The Ø2.7 Twist Drill is labeled as a 3.0mm drill in the kit case

Ø3.5mm




Bone Quality	Side Cut Drill	F3.5 Taper Drill	F4.0 Taper Drill	Ø3.5 Implant
Soft	►			Implant Placement
Normal	►	►		
Hard	►	►	►	

※ The recommended drill speed is 800–1,200 RPM, with an implant insertion torque of ≤40 Ncm.

※ The ET/EK implant is placed 1 mm subcrestally for normal and hard bone. However, soft bone may be placed at the bone level to ensure stability.

Ø4.0mm




Bone Quality	Lance Drill	F3.5 Taper Drill	F4.0 Taper Drill	F4.5 Taper Drill	Ø4.0 Implant
Soft	►	►			Implant Placement
Normal	►	►	►		
Hard	►	►		►	

※ The recommended drill speed is 800–1,200 RPM, with an implant insertion torque of ≤40 Ncm.

※ The ET/EK implant is placed 1 mm subcrestally for normal and hard bone. However, soft bone may be placed at the bone level to ensure stability.



## Ø4.5mm




Bone Quality	Lance Drill	F3.5 Taper Drill	F4.0 Taper Drill	F4.5 Taper Drill	F5.0 Taper Drill	Ø4.5 Implant
Soft	▶		▶			Implant Placement
Normal	▶	▶		▶		
Hard	▶	▶			▶	

※ The recommended drill speed is 800–1,200 RPM, with an implant insertion torque of  $\leq 40$  Ncm.

※ The ET/EK implant is placed 1 mm subcrestally for normal and hard bone. However, soft bone may be placed at the bone level to ensure stability.

## Ø5.0mm



Bone Quality	Lance Drill	F3.5 Taper Drill	F4.5 Taper Drill	F5.0 Taper Drill	F5.5 Taper Drill	Ø5.0 Implant
Soft	▶		▶			Implant Placement
Normal	▶	▶		▶		
Hard	▶	▶			▶	

※ The recommended drill speed is 800–1,200 RPM, with an implant insertion torque of  $\leq 40$  Ncm.


※ The ET/EK implant is placed 1 mm subcrestally for normal and hard bone. However, soft bone may be placed at the bone level to ensure stability.

## Drilling Sequence 122 Taper Drill

**ETIV**

(Length: 10mm)

### Ø4.0mm




Bone Quality	Ø2.2 Drill	F3.5 Taper Drill	Ø4.0 Implant
D4	▶		Implant Placement
Soft	▶	▶	

※ The recommended drill speed is 800–1,200 RPM, with an implant insertion torque of ≤40 Ncm.

※ The ET/EK implant is placed 1 mm subcrestally for normal and hard bone. However, soft bone may be placed at the bone level to ensure stability.

### Ø4.5mm




Bone Quality	Ø2.2 Drill	Ø3.0 Drill	F4.0 Taper Drill	Ø4.5 Implant
D4		▶		Implant Placement
Soft	▶		▶	

※ The recommended drill speed is 800–1,200 RPM, with an implant insertion torque of ≤40 Ncm.

※ The ET/EK implant is placed 1 mm subcrestally for normal and hard bone. However, soft bone may be placed at the bone level to ensure stability.

### Ø5.0mm



Bone Quality	Ø2.2 Drill	Ø3.0 Drill	F4.5 Taper Drill	Ø5.0 Implant
D4		▶		Implant Placement
Soft	▶		▶	

※ The recommended drill speed is 800–1,200 RPM, with an implant insertion torque of ≤40 Ncm.

※ The ET/EK implant is placed 1 mm subcrestally for normal and hard bone. However, soft bone may be placed at the bone level to ensure stability.



# Taper Kit (HKTAPEK)

For

EKIII

ETIII/IV


SSIII

Top panel components

**Torque Wrench**  
TQWCB



**Depth Gauge**  
ODG



# Taper Ultra Kit (HULTPK)


For **Ultra-wide**

Top panel components

**Open Wrench**  
SPOW



**Ratchet Wrench**  
RCWC




# Taper Kit Surgical Kit Instruments


Parallel Pin for Taper Drill			
Description	D/Ø	Item code	Image
<ul style="list-style-type: none"> <li>Parallel pin for taper drill</li> <li>Used for checking position and direction of bone preparation</li> <li>The lower part is for implant diameter drill and the upper part is for initial drill</li> <li>Color coding by implant diameter (F3.5: yellow, F4.0: green, F4.5: blue, F5.0: silver)</li> <li>122 Taper &amp; Taper Kit common components</li> </ul>	<b>F3.5</b>	HTPP3522	
	<b>F4.0</b>	HTPP4022	
	<b>F4.5</b>	HTPP4522	
	<b>F5.0</b>	HTPP5022	

Taper Drill							
Description			Image/Guide				
<ul style="list-style-type: none"> <li>Taper drill for taper(III type) implant by diameter and length</li> <li>Stopper drill with 1mm space</li> <li>Color coding displays implant diameter</li> <li>F3.5: yellow, F4.0: green, F4.5: blue, F5.0: red, F5.5: yellow</li> <li>Taper Kit single item</li> <li>(*excluded from 122 Taper Kit)</li> </ul>							
L	TL	D/Ø Y-Dim.	F3.5	F4.0	F4.5	F5.0	F5.5
5.0	29.5		HTPD3C3505	HTPD3C4005	HTPD3C4505	HTPD3C5005	-
6.0	30.5		HTPD3C3506	HTPD3C4006	HTPD3C4506	HTPD3C5006	HTPD3C5506
7.0	31.5		HTPD3C3507	HTPD3C4007	HTPD3C4507	HTPD3C5007	HTPD3C5507
8.5	33		HTPD3C3508	HTPD3C4008	HTPD3C4508	HTPD3C5008	HTPD3C5508
10	34.5		HTPD3C3510	HTPD3C4010	HTPD3C4510	HTPD3C5010	HTPD3C5510
11.5	34.5		HTPD3C3511	HTPD3C4011	HTPD3C4511	HTPD3C5011	HTPD3C5511
13	36		HTPD3C3513	HTPD3C4013	HTPD3C4513	HTPD3C5013	HTPD3C5513
15	38		HTPD3C3515	HTPD3C4015	HTPD3C4515	HTPD3C5015	HTPD3C5515
Color			Yellow	Green	Blue	Red	Yellow


LaceDrill (Guide Drill)			
Description	L	Item Code	Image
<ul style="list-style-type: none"> <li>Forming a hole to facilitate initial drilling</li> <li>Bone density determined through drilling</li> </ul>	<b>Long</b>	HAGDLC	

# Taper Kit Surgical Kit Instruments

Slidecut Drill		
Description	Item Code	Image
<ul style="list-style-type: none"> <li>Forming a hole to facilitate initial drilling</li> <li>Bone density determined through drilling</li> </ul>	HOSLMD20L	

Twist Drill (Stopper Drill)										
Description			Image/Guide							
<ul style="list-style-type: none"> <li>Included in New Hanaro KIT</li> <li>Long stopper (6mm)</li> <li>Enabling a procedure without drill extension for posterior region</li> <li>The color coded stopper indicates the drill length</li> </ul>										






L	TL	D/Ø Y-Dim.	F2.2 0.6	F3.0 0.9	F3.3 1.0	F3.6 1.0	F3.8 1.0	F4.1 1.0	F4.3 1.0	F4.3 1.0
6	30.5		2D2206LC	3D3006LC	-	-	3D3806LC	-	-	-
7	31.5		2D2207LC01	3D3007LC01	-	-	3D3807LC01	-	-	-
8.5	33		2D2208LC01	3D3008LC01	-	-	3D3808LC01	-	-	-
10	34.5		2D2210LC01	3D3010LC01	-	-	3D3810LC01	-	-	-
11.5	34.5		2D2211LC01	3D3011LC01	3D3311LC01	3D3611LC01	3D3811LC01	3D4111LC01	3D4311LC01	3D4611LC01
13	36		2D2213LC01	3D3013LC01	-	-	3D3813LC01	-	-	-

Twist Drill (Non Stopper Drill)										
Description			Image/Guide							
<ul style="list-style-type: none"> <li>Included in New Hanaro KIT</li> <li>Used for limited access for the Stopper Drill into the oral cavity</li> <li>See the image provided in the Non-stopper Drill section for the sizes of the drill marking lines for short/long types</li> </ul>										


TL	D/Ø	F1.5	F2.0	F2.2	F2.7	F3.0	F3.3
33	Short	2D1518FNLC	2D2018FNLC	2D2218FNLC	3D2718FNLC	3D3018FNLC	3D3318FNLC
41		-	-	2D2215FNLC01	3D2715FNLC01	3D3015FNLC01	3D3315FNLC01
	D/Ø	F3.6	F3.8	F4.1	F4.3	F4.6	
33	Long	3D3618FNLC	3D3818FNLC	3D4118FNLC	3D4318FNLC	3D4618FNLC	
41		3D3615FNLC01	3D3815FNLC01	3D4115FNLC01	3D4315FNLC01	3D4615FNLC01	





### Taper Cortical Drill for Taper Implant (ETIII, SSIII)

Description	D/Ø	Item code	Image
<ul style="list-style-type: none"> <li>The drill is used to remove cortical bone of the hard bone (used right after the use of Taper Drill)</li> <li>Dedicated drill for each implant diameter</li> <li>F3.5~5.0 drill marking line: bottom line 8.5mm or less, top line 10mm or more implant placement standard</li> <li>F5.5 drill marking line: bottom line 6mm or less, middle line 7mm, top line 8.5mm or more implant placement standard</li> <li>It is recommended to drill to the bottom of the marking line</li> <li>Taper Kit single item (excluded from 122 Taper Kit)</li> </ul>	<b>F3.5</b>	HTCD4C35	
	<b>F4.0</b>	HTCD4C40	
	<b>F4.5</b>	HTCD4C45	
	<b>F5.0</b>	HTCD4C50	
	<b>F5.5</b>	HTCD4C55	







### Taper Ultra Drill

Description	L	F6.0	F7.0	Image
<ul style="list-style-type: none"> <li>Taper drill for taper ultra-wide implant by diameter and length</li> <li>Stopper drill with 1mm space</li> <li>Color coding displays implant diameter</li> </ul>	<b>6.0</b>	HTPD3C6006	HTPD3C7006	 <p>Normal/ hard bone</p> <p>Soft bone</p>
	<b>7.0</b>	HTPD3C6007	HTPD3C7007	
	<b>8.5</b>	HTPD3C6008	HTPD3C7008	
	<b>10</b>	HTPD3C6010	HTPD3C7010	
	<b>11.5</b>	HTPD3C6011	HTPD3C7011	
	<b>13</b>	HTPD3C6013	HTPD3C7013	
	<b>Color</b>	<b>Green</b>	<b>Blue</b>	




### Implant Driver for EK




Description	C	Length	Item Code	Image
<ul style="list-style-type: none"> <li>Connects directly to an EK implant for final adjustments to the implant's depth</li> <li>*C = Connection</li> </ul>	<b>Regular</b>	<b>Short</b>	HKSFDS	 
		<b>Long</b>	HKSFDL	

### Implant Driver for ET

Description	C	Length	Item Code	Image
<ul style="list-style-type: none"> <li>Connects directly to an ET implant for final adjustments to the implant's depth</li> <li>*C = Connection</li> </ul>	<b>Mini</b>	<b>Short (17)</b>	HGSMFDS	  
		<b>Long (24)</b>	GSMFDL	
		<b>Ex. Long (34)</b>	HGSMFDE	
	<b>Regular</b>	<b>Short (19)</b>	HGSRFDS	  
		<b>Long (26.6)</b>	GSRFDL	
		<b>Ex. Long (33.6)</b>	HGSRFDE	

# Taper Kit Surgical Kit Instruments

NoMount Driver for EK				
Description	C	Length	Item Code	Image
<ul style="list-style-type: none"> <li>• Ø3.5 implant is combined with the bottom of the lower marking</li> <li>• Ø4.0, Ø4.5, Ø5.0, Ø6.0 and Ø7.0 implants are combined with the upper part of the lower marking</li> <li>• The distance between the two laser marking is 0.5mm</li> <li>• C= Connection</li> </ul>	Regular	Short (27.6)	HKSNMDCRS	
		Long (32.6)	HKSNMDCRL	
		Extra Long (37.6)	HKSNMDCRE	 <p>Ø4.0, Ø4.5, Ø5.0, Ø6.0, Ø7.0    Ø3.0, 3.5 Implant    Implant</p>


NoMount Torque Driver for ET				
Description	C	Length	Item Code	Image
<ul style="list-style-type: none"> <li>• Directly connects to an ET Implant for placement with a Torque Wrench</li> <li>• Ensure correct and complete seating before applying torque; loose connection may cause implant fracture</li> <li>• C= Connection</li> </ul>	Mini	Short (19)	HGSNMT32S	
		Long (26.6)	HGSNMT32L	
		Ex. Long (33.6)	HGSNMT32E	
	Regular	Short (19)	HGSNMT35S	
		Long (26.6)	HGSNMT35L	
		Ex. Long (33.6)	HGSNMT35E	

# Drilling Sequence Taper Drill

**EKIII | ETIII | SSIII**


(Length: 10mm)

## Ø3.2mm




Bone Quality	Ø2.2 Drill	Ø2.7 Drill	F3.0 Cortical Drill	Ø3.0 Implant
Soft	►			Implant Placement
Normal	►	►		
Hard	►	►	►	

## Ø3.5mm




Bone Quality	Ø2.2 Drill	Ø3.0 Drill	F3.5 Taper Drill	F3.5 Taper Cortical Drill	F3.5 Taper Implant Tap	Ø3.5 Implant
Soft	►	►				Implant Placement
Normal	►		►			
Hard	►		►	►		
Hard (Option)	►		►		►	

## Ø4.0mm




Bone Quality	Ø2.2 Drill	F3.5 Taper Drill	F4.0 Taper Drill	F4.0 Taper Cortical Drill	F4.0 Taper Implant Tap	Ø4.0 Implant
Soft	►	►				Implant Placement
Normal	►	►	►			
Hard	►	►	►	►		
Hard (Option)	►	►	►		►	

## Ø4.5mm




Bone Quality	Ø2.2 Drill	F3.5 Taper Drill	F4.0 Taper Drill	F4.5 Taper Drill	F4.5 Taper Cortical Drill	F4.5 Taper Implant Tap	Ø4.5 Implant
Soft	▶	▶	▶				Implant Placement
Normal	▶	▶		▶			
Hard	▶	▶		▶	▶		
Hard (Option)	▶	▶		▶		▶	

## Ø5.0mm



Bone Quality	Ø2.2 Drill	F3.5 Taper Drill	F4.5 Taper Drill	F5.0 Taper Drill	F5.0 Taper Cortical Drill	F5.0 Taper Implant Tap	Ø5.0 Implant
Soft	▶	▶	▶				Implant Placement
Normal	▶	▶	▶	▶			
Hard	▶	▶	▶	▶	▶		
Hard (Option)	▶	▶	▶	▶		▶	

## Ø5.5mm



Bone Quality	Ø2.2 Drill	F3.5 Taper Drill	F4.5 Taper Drill	F5.5 Taper Drill	F5.5 Taper Implant Tap	Ø5.5 Implant
Soft	▶	▶	▶			Implant Placement
Normal	▶	▶	▶	▶		
Hard	▶	▶	▶	▶	▶	

F5.5 taper cortical drill marking line: bottom line 6mm or less, middle line 7mm, top line 8.5mm or more

Recommended insertion torque ≤40Ncm. ET implant insertion depth in normal/hard bone is placed 1mm deeper than the bone level, and the soft bone is placed at the bone level to maintain initial stability.


Implant tap used in hard bone: 25rpm recommended with engine or use with torque wrench after fastening to mount extension (F5.0 implant tap: bottom line 7mm or less, top line 8.5mm or more).

## Drilling Sequence Taper Drill


**ETIV**

(Length: 10mm)


### Ø4.0mm



### Ø4.5mm



### Ø5.0mm




## Drilling Sequence Taper Drill

### ETIII Ultra-wide


(Length: 10mm)

#### Ø6.0mm



Bone Quality	Ø2.2 Drill	Ø2.0/3.0 Pilot Drill	Ø3.0 Drill	Ø3.0/3.8 Pilot Drill	Ø3.8 Drill	F6.0 Taper Drill	F6.0 Taper Drill	F6.0 Cortical Drill	Ø6.0 Implant
Soft	►	►	►	►		►			Implant Placement
Normal	►	►	►	►	►		►		
Hard	►	►	►	►	►		►	►	

#### Ø7.0mm



Bone Quality	Ø2.2 Drill	Ø2.0/3.0 Pilot Drill	Ø3.0 Drill	Ø3.0/3.8 Pilot Drill	Ø3.8 Drill	F6.0 Taper Drill	F7.0 Taper Drill	F7.0 Taper Drill	F7.0 Cortical Drill	Ø7.0 Implant
Soft	►	►	►	►		►	►			Implant Placement
Normal	►	►	►	►	►	►		►		
Hard	►	►	►	►	►	►		►	►	

Recommended placement torque less than 40Ncm.

ET implant placement depth. The normal/hard bone is placed 1mm deeper than bone level, and the soft bone is placed at the bone level to maintain initial stability.




## Drilling Sequence Taper Drill


### ETIV Ultra-wide

(Length: 10mm)

#### Ø6.0mm



#### Ø7.0mm



Recommended placement torque less than 40Ncm.

ET implant placement depth. The normal/hard bone is placed 1mm deeper than bone level, and the soft bone is placed at the bone level to maintain initial stability.





# 485 Kit (H485K)

For

EKIII

ETIII/IV

SSII/III

Ultra-Wide

Top panel components


**Torque Wrench**

TQWCB





**Depth Gauge**

ODG



# 485 Kit Surgical Kit Instruments

Guide Drill		
Description	Ø2.2	Image
<ul style="list-style-type: none"> <li>Drill for marking location of osteotomy to facilitate initial drilling</li> </ul>	OGD2027L	

485 Drill						
Description				Image/Guide		
<ul style="list-style-type: none"> <li>Included in 485 Kit</li> <li>A drill for placing short implants in alveolar bone lacking in vertical dimension</li> <li>Ø 2.2 drill: straight drill</li> <li>Except for Ø 2.2 drill, the top blade of the drill is in the shape of CAS Drill, and the side blade is in the shape of taper drill</li> <li>A stopper drill with 1mm margin</li> <li>Recommended drilling speed: 800~1,200rpm</li> </ul>						
L	Ø2.2	Pilot	F4.0	F4.5	F5.0	F5.5
4.0	H485D2204	H485D3504	H485D4004	H485D4504	H485D5004	H485D5504
5.0	H485D2205	-	H485D4005	H485D4505	H485D5005	H485D5505
6.0	H485D2206	H485D3506	H485D4006	H485D4506	H485D5006	H485D5506
7.0	H485D2207	-	H485D4007	H485D4507	H485D5007	H485D5507
8.5	H485D2208	-	H485D4008	H485D4508	H485D5008	H485D5508

# 485 Kit Surgical Kit Instruments

Implant Driver for EK				
Description	C	Length	Item Code	Image
<ul style="list-style-type: none"> <li>Connects directly to an EK implant for final adjustments to the implant's depth</li> <li>*C = Connection</li> </ul>	Regular	Short	HKSFDS	
		Long	HKSFDL	

Implant Driver for ET				
Description	C	Length	Item Code	Image
<ul style="list-style-type: none"> <li>Connects directly to an ET implant for final adjustments to the implant's depth</li> <li>*C = Connection</li> </ul>	Mini	Short (17)	HGSMFDS	
		Long (24)	GSMFDL	
		Ex. Long (34)	HGSMFDE	
	Regular	Short (19)	HGSRFDS	
		Long (26.6)	GSRFDL	
		Ex. Long (33.6)	HGSRFDE	

NoMount Driver for EK				
Description	C	Length	Item Code	Image
<ul style="list-style-type: none"> <li>Ø3.5 implant is combined with the bottom of the lower marking</li> <li>Ø4.0, Ø4.5, Ø5.0, Ø6.0 and Ø7.0 implants are combined with the upper part of the lower marking</li> <li>The distance between the two laser marking is 0.5mm</li> <li>C= Connection</li> </ul>	Regular	Short (27.6)	HKSNMDCRS	
		Long (32.6)	HKSNMDCRL	
		Extra Long (37.6)	HKSNMDCRE	


NoMount Torque Driver for ET				
Description	C	Length	Item Code	Image
<ul style="list-style-type: none"> <li>Directly connects to an ET Implant for placement with a Torque Wrench</li> <li>Ensure correct and complete seating before applying torque; loose connection may cause implant fracture</li> <li>C= Connection</li> </ul>	Mini	Short (19)	HGSNMT32S	
		Long (26.6)	HGSNMT32L	
		Ex. Long (33.6)	HGSNMT32E	
	Regular	Short (19)	HGSNMT35S	
		Long (26.6)	HGSNMT35L	
		Ex. Long (33.6)	HGSNMT35E	

## Drilling Sequence 485 Drill

**EKIII | ETIII/IV | SSII/III | ULTRA-WIDE**


(Length: 7mm)

### Ø4.0mm




Bone Quality	Twist Drill (Ø2.2)	Pilot Drill	485 Drill (F4.0)	485 Drill (F4.5)	Ø4.0 Implant
Normal	▶	▶	▶		Implant Placement
Hard	▶	▶		▶	

### Ø4.5mm



Bone Quality	Twist Drill (Ø2.2)	Pilot Drill	485 Drill (F4.5)	485 Drill (F5.0)	Ø4.5 Implant
Normal	▶	▶	▶		Implant Placement
Hard	▶	▶		▶	


### Ø5.0mm




Bone Quality	Twist Drill (Ø2.2)	Pilot Drill	485 Drill (F5.0)	485 Drill (F5.5)	Ø5.0 Implant
Normal	▶	▶	▶		Implant Placement
Hard	▶	▶		▶	


# Assist Kit (HOAK)

- Bone profilers are only sold in the packing unit of “Guide Screw + Bone Profiler”



# Assist Kit Surgical Instrument

Tissue Punch				Image/Guide
Description/Item code				
<ul style="list-style-type: none"><li>For flapless surgery</li><li>Laser marking to measure the height of gingiva, marked in 2mm increments</li><li>Packing unit: tissue punch + guide pin</li><li>Recommend to use a tissue punch smaller than the healing abutment by 0.7 to 1.5mm</li></ul>				
ET	SS	D/Ø	Item Code	
Ø 4.0/4.5	-	Ø3.3	HSTP33	
Ø 4.5/5.0	Ø 4.8	Ø3.8	HSTP38	
Ø 5.0	-	Ø4.3	HSTP43	
Ø 6.0	Ø 6.0	Ø4.8	HSTP48	
Ø 6.0	Ø 6.0	Ø5.3	HSTP53	
Application healing abutment standard				

Bone Profiler				Image/Guide
Description/Item code				
<ul style="list-style-type: none"><li>Used to remove bone around the implant after first or second stage surgery</li><li>Connect the Guide Screw to the implant in order to center the profiler</li><li>Guide screw protects the implant's platform from damage</li><li>Packing unit: bone profiler + guide screw</li><li>*C = Connection</li></ul>				
Guide Screw	D (Healing Abutment)	ET Mini / Regular	EK Regular	
ET (Mini + Regular)		HGSBGM	HGSBGS	
EK (3.0/3.5 + Regular)		HKS BGM	HKS BGS	
ET (Mini + Regular)	Ø4.0	HGSBP40	HKS BP40	
EK (3.0/3.5 + Regular)	Ø4.5	HGSBP45	HKS BP45	
Regular	Ø5.0	HGSBP50	HKS BP50	
	Ø6.0	HGSBP60	HKS BP60	
	Ø7.0	HGSBP70	HKS BP70	



### Trephine Drill

Description/Item code			Image/Guide
<ul style="list-style-type: none"> <li>Harvests bone or removes a failed implant</li> <li>Used to remove septal bone</li> <li>Can also be used as the initial drill for ultra-wide implants</li> </ul>			
D/Ø (Inner/Outer)	Short	Long	
3.7/4.5	HTD37S	HTD37	
4.2/5.0	HTD42S	HTD42	
4.7/5.5	HTD47S	HTD47	
5.2/6.0	HTD52S	HTD52	
5.7/6.5	HTD57S	HTD57	
6.2/7.0	HTD62S	HTD62	

### Crest Remover

Description/Item code			Image/Guide
<ul style="list-style-type: none"> <li>Marking the implant placement position after removing the narrow alveolar ridge horizontally</li> <li>Recommended drilling speed <ul style="list-style-type: none"> <li>- Angled type: 1,200~1,500rpm</li> <li>- Straight type: 15,000~30,000rpm</li> </ul> </li> </ul>			
L	29	45	
D Ø5.0	ERM50A	CERM50S	
D Ø7.0	-	CERM70A	

### AutoBone Collector®

Description/Item code			Image/Guide
<ul style="list-style-type: none"> <li>Used for autogenous bone collecting</li> <li>Comes in a Drill + Stopper set</li> <li>Recommended drilling speed: 300~600rpm</li> <li>Number of uses for the drill and stopper: 50 times</li> </ul> <p>※ Before initial drilling, connect the stopper to the first stage locking and harvest autogenous bone while drilling 4mm into the second stage locking (after harvesting, stop the drill and remove as it is with autogenous bone kept in the stopper)</p>			
D	Short (18.94)	Long (21.94)	
Ø3.0	ABC304S	ABC304L	
Ø4.0	ABC404S	ABC404L	
Ø5.0	ABC504S	ABC504L	
Ø6.0	ABC604S	ABC604L	



# Ultra Kit (HULTRK)

For


Ultra-Wide

Lower panel components


**Open Wrench**  
SPOW





**Ratchet Wrench**  
RCWC




# Ultra Kit Surgical Kit Instruments


Sidecut Drill		
Description	Item Code	Image
<ul style="list-style-type: none"> <li>Forming a hole to facilitate initial drilling</li> <li>Bone density determined through drilling</li> </ul>	HOSLMD20L	

Direct Drill			
Description	D1/D2	Item code	Image/Guide
<ul style="list-style-type: none"> <li>Direct drill: two-step drill that functions like a pilot and twist drill</li> <li>Final drilling is possible without using pilot drill</li> <li>Increases initial stability in an extraction socket due to the reduced dead space at the apex</li> </ul>	Ø4.6/5.2	H3D5213FNLC	
	Ø4.6/5.5	H3D5513FNLC	
	Ø5.5/6.2	H3D6213FNLC	
	Ø5.5/6.5	H3D6513FNLC	

Trial Pin for Ultra-wide			
Description	D	Item code	Image/Guide
<ul style="list-style-type: none"> <li>Measures the width and depth of a failed implant site</li> <li>Measures the drilling depth after using the direct drill as the final drill</li> <li>Also serves as a parallel pin</li> </ul>	Ø5.2	HUWFTP52	
	Ø5.5	HUWFTP55	
	Ø6.2	HUWFTP62	
	Ø6.5	HUWFTP65	

# Ultra Kit Surgical Kit Instruments

Cortical Drill for Ultra-wide			
Description	D/Ø	Item code	Image
<ul style="list-style-type: none"> <li>Trims cortical bone in hard bone cases (for ultra-wide type implants)</li> <li>Drill specifically designed for ultra-wide implant's unique diameter</li> <li>Drilling recommended until the bottom of the marker has been reached</li> </ul>	<b>F6.0</b>	HCD4C60	
	<b>F7.0</b>	HCD4C70	

Trepine Drill			
Description/Item code			Image/Guide
<ul style="list-style-type: none"> <li>Harvests bone or removes a failed implant</li> <li>Used to remove septal bone</li> <li>Can also be used as the initial drill for ultra-wide implants</li> </ul>			
D/Ø (Inner/Outer)	Short	Long	
3.7/4.5	HTD37S	HTD37S	
4.2/5.0	HTD42S	HTD42S	
4.7/5.5	HTD47S	HTD47S	
5.2/6.0	HTD52S	HTD52S	
5.7/6.5	HTD57S	HTD57S	
6.2/7.0	HTD62S	HTD62S	




## Drilling Sequence III Type Straight Drill

**EKIII | ETIII | SSIII**

(Length: 10mm)

### Ø3.2mm



Bone Quality	Ø2.2 Drill	Ø2.7 Drill	F3.0 Cortical Drill 2	Ø3.0 Implant
Soft	▶			Implant Placement
Normal	▶	▶		
Hard	▶	▶	▶	

### Ø3.5mm



Bone Quality	Ø2.2 Drill	Ø2.0/3.0 Pilot Drill	Ø3.0 Drill	F3.5 Cortical Drill 3	F3.5 Cortical Drill 3	Ø3.5 Implant
Soft	▶	▶	▶			Implant Placement
Normal	▶	▶	▶	▶		
Hard	▶	▶	▶		▶	

### Ø4.0mm



Bone Quality	Ø2.2 Drill	Ø2.0/3.0 Pilot Drill	Ø3.0 Drill	Ø3.3 Drill	F4.0 Cortical Drill 3	F4.0 Cortical Drill 3	Ø4.0 Implant
Soft	▶	▶	▶	▶			Implant Placement
Normal	▶	▶	▶	▶	▶		
Hard	▶	▶	▶	▶		▶	




## Drilling Sequence III Type Straight Drill

**EKIII | ETIII | SSIII**

(Length: 10mm)

### Ø4.5mm



Bone Quality	02.2 Drill	02.0/3.0 Pilot Drill	03.0 Drill	03.0/3.8 Pilot Drill	03.8 Drill	F4.5 Cortical Drill 3	F4.5 Cortical Drill 3	04.5 Implant
Soft	▶	▶	▶	▶	▶			Implant Placement
Normal	▶	▶	▶	▶	▶	▶		
Hard	▶	▶	▶	▶	▶		▶	

### Ø5.0mm



Bone Quality	02.2 Drill	02.0/3.0 Pilot Drill	03.0 Drill	03.0/3.8 Pilot Drill	03.8 Drill	04.3 Drill	F5.0 Cortical Drill 3	F5.0 Cortical Drill 3	05.0 Implant
Soft	▶	▶	▶	▶	▶				Implant Placement
Normal	▶	▶	▶	▶	▶	▶	▶		
Hard	▶	▶	▶	▶	▶	▶		▶	

### Ø5.5mm



Bone Quality	02.2 Drill	02.0/3.0 Pilot Drill	03.0 Drill	03.0/3.8 Pilot Drill	03.8 Drill	04.6 Drill	F5.5 Cortical Drill 3	F5.5 Cortical Drill 3	05.5 Implant
Soft	▶	▶	▶	▶	▶	▶			Implant Placement
Normal	▶	▶	▶	▶	▶	▶	▶		
Hard	▶	▶	▶	▶	▶	▶		▶	

Recommended insertion torque  $\leq 40\text{Ncm}$ .


ET implant insertion depth in normal/hard bone is placed 1mm deeper than the bone level, and the soft bone is placed at the bone level to maintain initial stability.

Drilling Sequence **IV Type** Straight Drill

**ETIV**


(Length: 10mm)

**Ø4.0mm**




Bone Quality	Ø2.2 Drill	Ø2.0/3.0 Pilot Drill	Ø3.0 Drill	Ø3.0 Drill (Half)	Ø4.0 Implant
D4	▶				Implant Placement
Soft	▶	▶	▶	▶	

**Ø4.5mm**



Bone Quality	Ø2.2 Drill	Ø2.0/3.0 Pilot Drill	Ø3.0 Drill	Ø3.0/3.8 Pilot Drill	Ø3.8 Drill	Ø4.1 Drill (Half)	Ø4.5 Implant
D4			▶				Implant Placement
Soft	▶	▶	▶	▶	▶	▶	

**Ø5.0mm**




Bone Quality	Ø2.2 Drill	Ø2.0/3.0 Pilot Drill	Ø3.0 Drill	Ø3.0/3.8 Pilot Drill	Ø3.8 Drill	Ø4.6 Drill (Half)	Ø5.0 Implant
D4			▶				Implant Placement
Soft	▶	▶	▶	▶	▶	▶	

## Drilling Sequence Ultra-wide Straight Drill

### ETIV Ultra-Wide


(Length: 10mm)

#### Ø6.0mm



Bone Quality	Ø2.2 Drill	Ø2.0/3.0 Pilot Drill	Ø3.0 Drill	Ø3.0/3.8 Pilot Drill	Ø3.8 Drill	Ø4.6 Drill	Ø5.2 Direct Drill	Ø6.0 Implant
D4	▶	▶			▶			Implant Placement
Soft	▶	▶	▶	▶	▶	▶	▶	


#### Ø7.0mm



Bone Quality	Ø2.2 Drill	Ø2.0/3.0 Pilot Drill	Ø3.0 Drill	Ø3.0/3.8 Pilot Drill	Ø3.8 Drill	Ø4.6 Drill	Ø5.5 Direct Drill	Ø6.2 Direct Drill	Ø7.0 Implant
D4	▶	▶			▶	▶			Implant Placement
Soft	▶	▶	▶	▶	▶	▶	▶	▶	

Recommended insertion torque ≤40Ncm


# Prosthetic Simple Kit (HPSK)



# Prosthetic Kit (HPRSTKA)


Top Panel Components


**Torque Wrench**  
TQWCB






# Prosthetic Kit



## Surgical Kit Instruments






Hand Driver					
Description	L	0.9 Hex	1.2 Hex	Image	
<ul style="list-style-type: none"> <li>Manual driver</li> <li>Tip holding function (except internal hex type)</li> <li>Internal hex type length: 11mm</li> </ul>	Ex.Short (8)	HAHD09MSH	HAHD12MSH		
	Short (13)	HAHD09SH	HAHD12SH		
	Middle (15)	-	HAHD12MH		
	Long (18)	HAHD09LH	HAHD12LH		
	Ex.Long (25)	-	HAHD12EH		






Machine Screw Driver					
Description	L	0.9 Hex	1.2 Hex	Image	
<ul style="list-style-type: none"> <li>1.2 hex driver for engine handpiece</li> <li>Tip holding function (except in internal hex type)</li> <li>Internal hex type length: 8mm</li> </ul>	Osstem Torque (5)	-	-		
	Short (5.6)	HAMSD09S	HAMSD12S		
	Long (11.6)	HAMSD09L	HAMSD12L		
	Ex.Long (17.6)	-	HAMSD12E		
<b>Application Driver Applied Products</b> (common for hand, machine screw and torque driver)		Cover screw (US mini)	Healing abutment, UCLA, Cemented abutment screw, Mount screw		

Torque Driver					
Description	L	0.5 Slot	0.9 Hex	1.2 Hex	Image
<ul style="list-style-type: none"> <li>Driver for torque wrench</li> <li>Tip holding function</li> <li>May bend or break if excessive torque is applied</li> <li>Damage is possible even at low torque if not fully engaged</li> <li>Apply vertical pressure when applying torque driver (Do not tilt)</li> <li>If the tip is bent or stripped, replace immediately</li> </ul>	Ex.Short (8)	-	-	HTRHD12MS	
	Short (13)	HTRSD05S	HTRHD09S	HTRHD12S	
	Middle (15)	-	-	HTRHD12M	
	Long (20)	HTRSD05L	HTRHD09L	HTRHD12L	
	Ex.Long (25)	HTRSD05E	-	HTRHD12SE	

Angled Torque Driver				
Description	L	1.2 Hex	1.2 Hex (Set)	Image
<ul style="list-style-type: none"> <li>Driver for torque wrench</li> <li>No holding function</li> <li>Recommended tightening torque: 30Ncm (excessive torque causes fracture)</li> <li>Do not remove tube to prevent fragmentation when broken</li> <li>Recommended number of use: 10 times</li> <li>Set: 3 per pack</li> </ul>	Short (13)	HATD12S	HATD12S3S	
	Long (20)	HATD12L	HATD12L3S	

Repair Torque Driver			
Description	L	1.2 Hex	Image
<ul style="list-style-type: none"> <li>Reduced diameter of shank compared to the Torque Driver (Ø2.1 → 1.6)</li> <li>The diameter of the screw hole can be minimized during prosthetic repair or SCRP procedures</li> </ul>	Short (13)	HTRHD12SR	
	Long (20)	HTRHD12LR	

Solid Abutment Driver		Regular		Wide	
Description	Type	Short (6)	Long (12)	Short (10)	Long (L)
<ul style="list-style-type: none"> <li>Driver specific for solid abutments</li> <li>Connect to the solid abutment by matching up the groove with the triangular indicator</li> <li>Recommended tightening torque: 30Ncm</li> </ul>	Square	 HSDSS	 HSDSL	 HSD60S	-
	Round	 HSDRS	 HSDRL	-	-

Excellent Solid Abutment Driver		Regular		Wide	
Description	Type	Short (6)	Long (12)	Short (10)	Long (L)
<ul style="list-style-type: none"> <li>Driver for Excellent Solid Abutments</li> <li>Connect to the solid abutment by matching up the groove with the triangular indicator</li> <li>Recommended tightening torque: 30Ncm</li> </ul>	Square	 HESDSS	 HESDSL	 HESD60S	-
	Round	 HESDRS	 HESDRL	-	-



# Prosthetic Kit

## Surgical Kit Instruments


O-ring Abutment Driver		
Description	Item Code	Image
<ul style="list-style-type: none"> <li>Driver for Stud Abutment</li> </ul>	HAORD	


Rigid Outer Driver				
Description	D/Ø (Abutment)	Short (16.5)	Long (21.5)	Image
<ul style="list-style-type: none"> <li>Driver for Rigid Abutments</li> <li>Recommended tightening torque: 30Ncm</li> </ul>	Ø4.0	HRDMS	HORDML	
	Ø4.5	HRD45S	HRD45L	
	Ø5.0	HRDRS	HRDRL	
	Ø6.0	HRDWS	HRDWL	


Octa Abutment Driver			
Description	Type	Short	Long
<ul style="list-style-type: none"> <li>Driver for Octa Abutments</li> <li>Recommended tightening torque: 30Ncm</li> </ul>	Square	12.5 HODSS	18.5 HDSL
	Round	13.4 HDRS	19.4 HDRL


Multi Abutment Machine Driver		
Description	Item Code	Image
<ul style="list-style-type: none"> <li>Machine driver for Multi-Abutments</li> </ul>	HMAMD	


Abutment Holder		
Description	Item Code	Image
<ul style="list-style-type: none"> <li>Used to hold abutments and help deliver them to hard-to-reach sites of the oral cavity</li> </ul>	HOABH	

Multi Abutment Outer Driver			
Description	Type	Item Code	Image
<ul style="list-style-type: none"> <li>Torque driver for Multi-Abutments</li> </ul>	Normal	HMAODP	
	Rescue	HMAOD	

Locator Abutment Driver			
Description	Type	Item Code	Image
<ul style="list-style-type: none"> <li>Torque driver for Locator Abutments</li> </ul>	Short	TWLDSK	
	Long	TWLDLK	



Torque Driver				
Description	Type	Short (10)	Long (15)	Image
<ul style="list-style-type: none"> <li>It may not be fastened or disconnected when connecting a normal handpiece</li> <li>Driver should be used after matching the groove or section of the outer triangle and abutment</li> <li>Solid, excellent solid driver is compatible only with Ø4.8</li> <li>1.2 hex type L is 5</li> </ul>	1.2 Hex	HTH12S	-	
	Rigid 4.0	HTR40S	HTR40L	
	Rigid 4.5	HTR45S	HTR45L	
	Rigid 5.0	HTR50S	HTR50L	
	Rigid 6.0	HTR60S	HTR60L	
	Solid	HTS48S	HTS48L	
	Excellent Solid	HTE48S	HTE48L	


Path Probe for ET			
Description	Connection	Item Code	Image
<ul style="list-style-type: none"> <li>Tool to check path and measure gingival height after ET implant placement</li> <li>*C = Connection</li> </ul>	Mini	GIPAP-3016A	
	Regular	GIPAP-3516A	


Path Probe for EK			
Description	Connection	Item Code	Image
<ul style="list-style-type: none"> <li>Tool to check path and measure gingival height after EK implant placement</li> <li>*C = Connection</li> </ul>	Regular	HKSPPR	


# Prosthetic Kit

## Surgical Kit Instruments

Finishing Reamer Set		
Description	Item Code	Image
<ul style="list-style-type: none"> <li>Tool used to remove the excess cast after plastic coping is set</li> </ul> <p><b>Reamer user guide</b></p> <ol style="list-style-type: none"> <li>1. Select a reamer tip that is the same size as Abutment size and connect it to the burn-out cylinder</li> <li>2. Firmly grasp the casting body and rotate the Reamer Bite with consistent force</li> <li>3. Ream the body until it is clean and free of the excess casting</li> </ol> 	HFRSC	

Reamer Bite		
Description	Item Code	Image
<ul style="list-style-type: none"> <li>Tool used to remove the excess cast after plastic coping is set</li> </ul>	HFRBC	

Reamer Tip for Rigid Abutment			
Description	D/Ø	Item Code	Image
<ul style="list-style-type: none"> <li>Tool used to remove the excess cast after plastic coping is set</li> </ul>	Ø 4.0	HGSRFRT400	
	Ø 4.5	HGSRFRT450	
	Ø 5.0	HGSRFRT500	
	Ø 6.0	HGSRFRT600	

Reamer Tip for Solid, Excellent Solid Abutment				
Description	Platform	Solid	Ex. Solid	Image
<ul style="list-style-type: none"> <li>Tool used to remove the excess cast after plastic coping is set</li> <li>For both solid Ø6.0 and excellent solid Ø4.8</li> <li>P = Platform</li> </ul>	Regular Ø4.8	HFRTS480	HFRTS480	
	Wide Ø6.0	HFRTS600	HFRTS600	



# ESR Kit Easy Screw Removal Kit (HESRK)

※ ESR (EK) Kit: HKESRK


For

EKIII

ETIII/IV

SSIII


Ultra-Wide













# ESR Full Kit Easy Screw Removal Full Kit (HESRFK\_US)

For





Nobel Biocare Active/Replace / Straumann Bone Level / Astra Osseo Speed TX  
3i Full OSSEOTITE Tapered Certain / 3i Tapered / BioHorizons® Internal





# ESR Full Kit Surgical Kit Instruments


Guide							
Description							
<ul style="list-style-type: none"><li>• Connected to the implant to center and prevent shaking of the reverse driver, SR drill and retap</li><li>• Choose spec of guide according to implant type and diameter (6 overseas companies' internal and submerged type products)</li><li>• Short or long used according to intermaxillary distance</li><li>• Used in combination with the ESR handle</li><li>• *C = Connection/the number of use: 10 times</li></ul>							
Osstem							
Type	Length	C	Mini	Regular	Wide	Image	
ET(TS)	Short		HGTMS	HOGTRS	-		
	Long		HOGTML	HOGTRL	-		
EK(KS)	Short		-	HKGRS	-		
	Long		-	HKGRL	-		
SS	Short		-	HGSRS	HGSRS		
	Long		-	HOGSRL	HOGSRL		
Nobel Biocare							
Type	Length	F	Ø3.5	Ø4.3	Ø5.0	Ø6.0	Image
Active	Short		HGNA01S	HGNA02S	HGNA02S	-	
	Long		HGNA01L	HGNA02L	HGNA02L	-	
Replace	Short		-	HGNR02S	HGNR03S	HGNR04S	
	Long		-	HGNR02L	HGNR03L	HGNR04L	
Type	Length	F	Ø3.3	Ø3.75	Ø4.0	Ø5.0	Image
MKIII	Short		HGUMS	HGURS	HGURS	HGUWS	
	Long		HHGUML	HOGURL	HOGURL	HOGUWL	
Straumann							
Type	Length	F	NC (3.3)	RC (4.1)	RC (4.8)		Image
Bone Level	Short		HGSB01S	HGSB02S	HGSB02S		
	Long		HGSB01L	HGSB02L	HGSB02L		
Type	Length	F	RN (3.3 / 4.1 / 4.8)		WN (4.8)		Image
Roxolid SLActive	Short		HGSTRS		HGSTRS		
	Long		HGSTRL		HGSTRL		
Astra							
Type	Length	F	Small (3.5 S / 4.0 S)		Large (4.5 / 5.0 / 5.0S)		Image
Osseo Speed TX	Short		HGAO01S		HGAO02S		
	Long		HGAO01L		HGAO02L		
Zimvie (Zimmer)							
Type	Length	F	Green (3.7 / 4.1 / 4.7)		Green (6.0)		Image
Tapered	Short		HGZB01S		HGZB02S		
	Long		HGZB01L		HGZB02L		





BioHorizons					
Type	Length	F	Yellow / Green	Blue	Image
Internal (Tapered Bone Level)	Short		HGZB01S	HGZB02S	
	Long		HGZB01L	HGZB02L	
Type	Length	F	Ø3.5	Ø4.0 / 5.0 / 6.0	
External	Short		HGUMS	HGURS	
	Long		HHGUML	HOGURL	
3i					
Type	Length	F	3.25	4.0 / 5.0 / 6.0	Image
Full Osseotite Tapered Certain	Short		HGIF01S	HGIF02S	
	Long		HGIF01L	HGIF02L	
Full Osseotite Tapered	Short		-	HGURS	
	Long		-	HOGURL	




Reverse Drill				
Description	Type	Short	Long	Image
<ul style="list-style-type: none"> <li>Used to remove a fractured screw</li> <li>To be used in conjunction with the guide</li> <li>If the red marking of the reverse driver is visible on the guide, remove the fractured screw using a screw holder</li> <li>Direction of rotation: Counterclockwise</li> <li>Recommended number of usage: 10 times</li> </ul>	M1.6	-	HARVDMLC	
	M1.8	HARVDRSC	HARVDRLC	
	M2.0	HARVDRSC	HARVDRLC	




Screw Removal Drill (SR Drill)				
Description	Type	Short	Long	Image
<ul style="list-style-type: none"> <li>Used to create a hole in the fractured screw</li> <li>Make sure to connect the guide, irrigate and suction to remove any debris</li> <li>Available in long and short lengths for different intermaxillary distances</li> <li>Drill until the colored marking on the drill is no longer visible through the guide</li> <li>Recommended speed: 1,200~1,500 rpm (counterclockwise)</li> <li>Recommended number of uses: 5 times</li> <li>Connect the guide before use/Do not apply excessive vertical force/Do not clean with hydrogen peroxide</li> </ul>	M1.6	HSRD08S	HOSRD08L	
	M1.8	HSRD09S	HSRD09L	
	M2.0	HSRD10S	HOSRD10L	




Torque Driver Handle		
Description	Item Code	Image
<ul style="list-style-type: none"> <li>Manual handle for SR Tip, AR Tip, screw holder</li> </ul>	HMSTH	


# ESR Full Kit Surgical Kit Instruments



Reverse Driver				
Description	F	Short	Long	Image
<ul style="list-style-type: none"> <li>Reverse Driver</li> <li>Used to remove a fractured screw</li> <li>To be used in conjunction with the guide</li> <li>Insert until the red band is in the Guide and turn counterclockwise to remove the screw</li> <li>Use manually/Rotate counterclockwise/</li> <li>Number of usages: 10 times</li> </ul>	Mini	-	HRVDML	
	Regular/ Wide	HRVDRS	HRVDRL	


Screw Removal Tip (SR Tip)				
Description	Type	Short	Long	Image
<ul style="list-style-type: none"> <li>Removes a fractured screw by engaging into the hole created by the Screw Removal Drill</li> <li>Rotation direction: counterclockwise</li> </ul>	M1.6	HRT16S	HORT16L	
	M1.8	HRT18S	HRT18L	
	M2.0	HRTS	HORTL	



Screw Holder			
Description	Type	Item Code	Image
<ul style="list-style-type: none"> <li>Grasps onto a protruding fractured screw to unscrew it</li> <li>Color-coded for easy recognition</li> <li>Rotation direction: counterclockwise</li> </ul>	M1.6	HOSHM	
	M1.8	HSHR18	
	M2.0	HOSHR	

Re-Tap			
Description	Type	Item Code	Image
<ul style="list-style-type: none"> <li>Re-threads the internal connection of a implant</li> <li>Connects to a torque wrench or ratchet wrench to re-thread by hand</li> </ul>	M1.6	HORTM	
	M1.8	HRTR18	
	M2.0	HORTR	

ESR Handle		
Description	Item Code	Image
<ul style="list-style-type: none"> <li>Stabilizes the Guide to the implant</li> </ul>	OARH	

Abutment Removal Tip (AR Tip)					
Description	F	Short	Long	Ex. Long	Image
<ul style="list-style-type: none"> <li>Removes fractured or jammed abutments and mounts from the implant</li> <li>Insert into fractured abutment hole, turn counterclockwise, and rock back and forth to loosen and remove with forceps</li> <li>Mini: it can be used to remove a screw with a stripped hex               <ul style="list-style-type: none"> <li>To remove the screw, engage the tip into the stripped hex and rotate counterclockwise</li> </ul> </li> </ul>	Mini	HARTMS	HOARTML	HARTMEL	
	Regular	HARTRS	HOARTRL	HARTREL	

Slot Driver		
Description	Item Code	Image
<ul style="list-style-type: none"> <li>Used to unscrew a screw, healing abutment, cover screw, or abutment screw with a stripped hex after creating a slot with a Ø0.8 bur</li> </ul>	HOTSD07	

Transfer Abutment Separate Tool			
Description		Item Code	Image
<ul style="list-style-type: none"> <li>Remove jammed abutment of non-hex type transfer abutment</li> <li>The tip is for mini platform abutments; the next step for regular platform</li> <li>Remove the abutment screw, insert Separate Tool Body into the abutment, tighten clockwise with Driver, and remove the abutment. If there is difficulty separating the abutment, attach a ratchet wrench for extra torque</li> </ul>	Driver	HTASD	
	Body	HTASB	
	Set	HTAST	

# EIR Kit Easy Implant Removal Kit (HKSFRK)

※ EIR (EK) Kit: HKSFRK

For

EKIII


ETIII/IV

SSIII

Ultra-Wide

Top panel components

**Implant Wrench**  
HFRDFE



# EIR Full Kit Easy Implant Removal Full Kit (HSFRFK\_US)

**For** Nobel Biocare Active/Replace / Straumann Bone Level / Astra Osseo Speed TX  
3i Full OSSEOTITE Tapered Certain / Zimmer Tapered / BioHorizons® Internal

Lower panel components

**Implant Wrench**  
HFRDFE



**Torque Wrench**  
HTW400B



HFRBM35S  
HFRBM35L  
HFRBR40S  
HFRBR40L  
HFRBW50S  
HFRBW50L

**Remover Body**

Normal  
HFRSM35  
HFRSR40  
HFRSW50


**Remover Screw**  
(ET/SS)

Fracture  
HFRSM35F  
HFRSR40F  
HFRSW50F

Normal  
HKSFRSM35  
HKSFRSR40  
HKSFRSW50

**Remover Screw**  
KS/EK


Fracture  
HKSFRSM35F  
HKSFRSR40F  
HKSFRSW50F





**Torque Extension**  
HOTE


**Screw Driver**  
HFRSDM23  
HFRSDR25  
HFRSDW30


# EIR Full Kit Surgical Kit Instruments


Remover Screw					
Description			Image		
<ul style="list-style-type: none"><li>When securely fastened to the implant, it serves as a supporting structure for the Remover Body</li><li>Use the proper type that matches the diameter of the implant to be removed (ET/SS/US, normal/fracture)</li><li>Fracture type is specifically for removing a fractured implant</li><li>Recommended tightening torque: regular/wide 100Ncm, mini 80Ncm</li></ul> ※ Disposable; do not reuse					
Osstem					
Type	Mode	Mini Ø3.5	Regular Ø4.0~4.5/P4.8	Wide Ø5.0/P6.0	
ET/SS	Normal	HFRSM35 (3N)	HFRSR40 (4N)	HFRSW50 (5N)	
	Fracture	HFRSM35F (3F)	HFRSR40F (4F)	HFRSW50F (5F)	
EK	Normal	HKSFRSM35 (KS3N)	HKSFRSR40 (KS4N)	HKSFRSW50 (KS5N)	
	Fracture	HKSFRSM35F (KS3F)	HKSFRSR40F (KS4N)	HKSFRSW50F (KS5N)	
Nobel Biocare					
Type	Mode	Mini Ø3.5	Regular Ø4.3	Wide Ø5.0/6.0	
Active	Normal	HFRSMNA35 (N1)	HFRSR40 (4N)	HFRSW50 (5N)	
	Fracture	HFRSMNA35F (N2)	HFRSR40F(4F)	HFRSW50F (5F)	
Replace	Normal	HFRSMNR35 (N3)	HFRSR40 (4N)	HFRSW50 (5N)	
	Fracture	HFRSMNR35F (N4)	HFRSR40F (4F)	HFRSW50F (5F)	
Straumann					
Type	Mode	Mini Ø3.3	Regular Ø4.1	Wide Ø4.8	
Bone Level	Normal	HFRSMS33 (S1)	HFRSRS41 (S2)	HFRSWS48 (S3)	
	Fracture	HFRSMS33F (S4)	HFRSRS41F (S5)	HFRSWS48F (S6)	
Astra					
Type	Mode	Mini Ø3.5	Regular Ø4.0	Regular Ø4.5	Wide Ø5.0
Osseo Speed TX	Normal	HFRSMNA35 (N1)	HFRSRA40 (A1)	HFRSR40 (4N)	HFRSW50 (5N)
	Fracture	HFRSMNA35F (N2)	HFRSRA40F (A2)	HFRSR40F (4F)	HFRSW50F (5F)
3i					
Type	Mode	Mini Ø3.25	Regular Ø4.0	Wide Ø5.0/6.0	
Full Osseotite Tapered Certain	Normal	HFRSMS33 (S1)	HFRSRI40 (S2)	HFRSWI50 (S3)	
	Fracture	HFRSMS33F (S4)	HFRSRI40F (S5)	HFRSWI50F (S6)	
Zimmer					
Type	Mode	Mini Ø3.7	Regular Ø4.1	Wide Ø4.7	Ultra-Wide Ø6.0
Tapered	Normal	HFRSMZ37 (Z1)	HFRSRZ41 (Z2)	HFRSWZ47 (Z3)	HFRSWZ60 (Z4)
	Fracture	HFRSMZ37F (Z4)	HFRSRZ41F (Z6)	HFRSWZ47F (Z7)	HFRSWZ47F (Z7)
Biohorizons					
Type	Mode	Mini Ø3.8	Regular Ø4.6	Wide Ø5.8	
Internal	Normal	HFRSRZ41 (i1)	HFRSWZ47 (i2)	HFRSWZ60 (i3)	
	Fracture	HFRSRZ41F (i4)	HFRSWB46F (i5)	HFRSWB46F (i6)	

Screw Driver			
Description	Type	Item Code	Image
<ul style="list-style-type: none"> <li>Connects and fastens the Remover Screw to the implant</li> <li>Recommended tightening torque: regular/wide 100Ncm, mini 80Ncm</li> </ul>	Mini	HFRSDM23	
	Regular	HFRSDR25	
	Wide	HFRSDW30	

Remover Body				
Description	Type	Short	Long	Image
<ul style="list-style-type: none"> <li>Connects to the Remover Screw to apply torque to remove the implant</li> <li>Select the correct type that matches the diameter of the implant to be removed</li> </ul> <p>※ Disposable; do not reuse</p>	Mini	HFRBM35S	HFRBM35L	
	Regular	HFRBR40S	HFRBR40L	
	Wide	HFRBW50S	HFRBW50L	
		HFRBW57S	HFRBW57L	
	Ultra-Wide	HFRBUW60S	HFRBUW60L	


Torque Extension		
Description	Item Code	Image
<ul style="list-style-type: none"> <li>Extends the length of the screw driver and remover body (by 10mm)</li> </ul>	HOTE	

Torque Wrench		
Description	Item Code	Image
<ul style="list-style-type: none"> <li>Tightens Screw Driver and removes the implant using the Remover Body</li> <li>Applies up to 400Ncm of torque (markers at 80/100/200/300/400Ncm)</li> <li>Torque by pulling the bar back until it reaches the desired torque value marking</li> <li>Clean and sterilize for storage</li> </ul>	HTW400B	


Implant Wrench		
Description	Item Code	Image
<ul style="list-style-type: none"> <li>Wrench used to separate the implant from the Remover Body</li> </ul>	HFRDFE	





# IM-Cure Kit (HICK)






# IM-Cure Kit Surgical Kit Instruments

Metal Probe		
Description	Item Code	Image
<ul style="list-style-type: none"> <li>Used to diagnose periodontal disease</li> <li>Measures pocket depth/size</li> <li>Marking lines of 1mm increments</li> </ul>	HICPM	


Plastic Probe		
Description	Item Code	Image
<ul style="list-style-type: none"> <li>Used to diagnose periodontal disease</li> <li>Measures pocket depth/size</li> <li>Marking lines of 1mm increments</li> <li>Plastic material prevents scratches on implant</li> <li>Flexible probe makes it ideal for curved shape of alveolar bone</li> <li>Autoclavable</li> </ul>	HICPP	



Curette			
Description	Type	Item Code	Image
<ul style="list-style-type: none"> <li>Tool to remove granulation tissue firmly attached to a specific area</li> <li>Gracey curette</li> <li>01-02: used for removal of anterior tissue</li> <li>11-12: used for removal of ganglion tissue</li> <li>13-14: used to remove the tissue from the distal part of posterior teeth</li> </ul>	01-02	HICC0102	
	11-12	HICC1112	
	13-14	HICC1314	



Protect Screw			
Description	Type	Item Code	Image
<ul style="list-style-type: none"> <li>Protect the implant's internal connection when using the SmartBrush 1 &amp; 2</li> <li>Torque Using a 1.2 hex driver, tighten to about 5Ncm</li> </ul>	Mini	OPISB2PM	
	Regular	OPISB2PR	



SmartBrush 1			
Description	Short	Long	Image
<ul style="list-style-type: none"> <li>Bristle designed for effective debridement of the implant surface</li> <li>After removing the patient's prosthesis and abutment, connect the protect screw before using</li> <li>Recommended speed: 1,200~1,500rpm</li> <li>Recommended use time: approximately 1 minute per screw thread (not recommended over 4 minutes)</li> <li>Must use irrigation and suction during polishing</li> </ul> <p>※ Disposable. Do not reuse</p>	OPISB1S	OPISB1L	 26  32

# IM-Cure Kit Surgical Kit Instruments

SmartBrush 2				
Description	D/Ø	Short	Long	Image
<ul style="list-style-type: none"> <li>Debride implant</li> <li>After removing the patient's prosthesis and abutment, connect the protect screw before using</li> <li>Must use irrigation and suction during polishing</li> <li>Recommended speed: 1,200~1,500rpm</li> <li>Recommended use time: 1~2 minutes</li> <li>Excessive use longer than 3 minutes may cause the product to break or bend</li> </ul>	F3.0/F3.5	OPISB23S	OPISB23L	
	F4.0/F4.5	OPISB24S	OPISB24L	
	F5.0/F5.5	OPISB25S	OPISB25L	
	F6.0	OPISB26S	OPISB26L	
	F7.0	OPISB27S	OPISB27L	

SmartScaler - Metal			
Description	Brand	Item Code	Image
<ul style="list-style-type: none"> <li>Used to remove substances such as tartar by connecting it to an ultrasonic scaler</li> <li>Secondary use after using SmartBrush 1 or SmartBrush 2</li> <li>Easy to bend tip for easy access</li> <li>EMS, KaVo, SATELEC specifications</li> </ul>	EMS	HICSME	  Bendable
	KaVo	HICSMK	
	SATELEC	HICSMS	

SmartScaler - Plastic				
Description	Brand	125°	100°	Image
<ul style="list-style-type: none"> <li>Used in combination with SmartScaler plastic tip</li> <li>Do not use for removal of debris on the implant surface</li> <li>EMS, KaVo, SATELEC specifications</li> <li>A = Angle</li> </ul>	EMS	HICSCE	HICSCE2	  125° 100°
	KaVo	HICSCK	HICSCK2	
	SATELEC	HICSCS	HICSCS2	

SmartScaler Plastic Tip		
Description	Item Code	Image
<ul style="list-style-type: none"> <li>Used to remove substances from abutment or crown by attaching the SmartScaler</li> <li>※ Do not use to implant surface</li> <li>Packing unit: 30 per set</li> </ul>	PSTP	  Tighten



# CAS Kit (HCRSNK)

For

EKIII

ETIII/IV

SSII/III

Ultra-Wide

Top panel components

Hydraulic Membrane  
Lifter Tube  
SNMT



Lower panel components


Bone Carrier Head  
SNBCH30



Bone Carrier  
SNBCS35



Bone Condenser  
SNBC1114



# CAS Kit Surgical Kit Instruments

CAS Drill				
Description	D/Ø	Short	Long	Image
<ul style="list-style-type: none"> <li>Specialized drill designed to penetrate the sinus floor without damaging the Schneiderian membrane by forming a conical bone lid</li> <li>Four blade drill design provides superior bone removing capability from low to high speeds, and can collect autogenous bone at low speeds</li> <li>Safely advance to the floor of the sinus using stoppers (1mm increments)</li> <li>Final drill diameter is based on bone density, regardless of the implant type (straight or tapered)</li> <li>Recommended speed: 400~800rpm (For beginner: 400rpm)</li> </ul>	Ø2.8	SNDR2813TS	SNDR2813TL	
	Ø3.1	SNDR3113TS	SNDR3113TL	
	Ø3.3	SNDR3313TS	SNDR3313TL	
	Ø3.6	SNDR3613TS	SNDR3613TL	
	Ø3.8	SNDR3813TS	SNDR3813TL	
	Ø4.1	SNDR4113TS	SNDR4113TL	


Guide Drill			
Description	D/Ø	Item Code	Image
<ul style="list-style-type: none"> <li>Drill to mark the implant placement site</li> <li>Side cutting blades can be used to modify the side walls of the extraction socket</li> <li>Line marking 2mm from the apex of drill</li> </ul>	Ø2.0/2.7	SNGD2027TL	



Ø2.2 Twist Drill			
Description	D/Ø	Item Code	Image
<ul style="list-style-type: none"> <li>Recommended to under-drill by 1mm below the floor of the sinus</li> <li>Use with stoppers for safe and controlled drilling</li> <li>Apex tip measures an additional 0.6mm</li> </ul>	Ø2.2	SNTD2213TL	


Hydraulic Membrane Lifter Set			
Description	D/Ø	Item Code	Image
<ul style="list-style-type: none"> <li>Hydraulic pressure is used to separate and lift the sinus membrane</li> <li>Securely fits Ø2.8~Ø4.1 CAS drilled osteotomies</li> </ul>	Ø2.6/6.0	HOCHML	


Stopper											
<ul style="list-style-type: none"> <li>Laser mark numbers indicates drilling depth</li> <li>Color-coded by length</li> <li>Recommended number of usage: 50 times</li> </ul>											
Length	2	3	4	5	6	7	8	9	10	11	12
Item Code	SNST2	SNST3	SNST4	SNST5	SNST6	SNST7	SNST8	SNST9	SNST10	SNST11	SNST12
Color	Purple	Blue	Yellow	Purple	Blue	Yellow	Purple	Blue	Yellow	Purple	Blue


# CAS Kit Surgical Kit Instruments

Bone Carrier		
Description	Item Code	Image
<ul style="list-style-type: none"> <li>Handle for the bone carrier head</li> <li>Connects the bone carrier head and tighten at the opposite end</li> <li>Connects both heads (SNBCH30 or SNBCH35)</li> </ul>	SNBCS35	


Bone Carrier Head			
Description	D/Ø	Item Code	Image
<ul style="list-style-type: none"> <li>Cone shaped with an extended tip that reaches the sinus cavity and prevents bone material from spilling out</li> <li>SNBCH30 for Ø3.1/3.3 CAS drilled osteotomy</li> <li>SNBCH35 for Ø3.6/3.8/4.1 CAS drilled osteotomy</li> <li>Fill the reservoir (up to the marker) with bone material in small quantities using the bone condenser. Repeat the process as necessary</li> </ul>	Ø3.1	SNBCH30	
	Ø3.6	SNBCH35	



Bone Condenser		
Description	Item Code	Image
<ul style="list-style-type: none"> <li>Safely pushes bone material through the bone carrier into the sinus cavity</li> <li>SNBCH30: use Ø1.1/SNBCH35: use Ø1.4</li> </ul>	SNBC1114	


Hydraulic Membrane Lifter Tube		
Description	Item Code	Image
<ul style="list-style-type: none"> <li>Tubing to connect the hydraulic membrane lifter to the saline filled syringe</li> </ul>	SNMT	

Membrane Lifter				
Description	D/Ø	Short	Long	Image
<ul style="list-style-type: none"> <li>Round shape with no cutting edges for safe membrane lifting</li> <li>Select Membrane Lifter diameter according to the CAS drill diameter used (head diameter: CAS drill diameter -0.2mm)</li> <li>Use CAS stoppers to control length</li> <li>Recommended speed: 30~50rpm (Recommended speed for beginner: 30rpm)</li> <li>Be sure to spray water when using</li> </ul>	Ø2.6	SNML2813TS	SNDR2813TL	
	Ø2.9	SNML3113TS	SNDR3113TL	
	Ø3.1	SNML3313TS	SNDR3313TL	
	Ø3.4	SNML3613TS	SNDR3613TL	
	Ø3.6	SNML3813TS	SNDR3813TL	
	Ø3.9	SNML4113TS	SNDR4113TL	




Depth Gauge		
Description	Item Code	Image
<ul style="list-style-type: none"> <li>Checks for access into the sinus cavity and measures the thickness of residual bone</li> </ul>	SNDG	

Bone Spreader			
Description	D/Ø	Item Code	Image
<ul style="list-style-type: none"> <li>Spreads bone graft material in the sinus cavity</li> <li>Used with CAS stoppers</li> <li>Recommended speed: 30rpm or less</li> </ul>	Ø2.0	SNBS2015T	
	Ø3.0	SNBS3015T	

Y-Connector		
Description	Item Code	Image
<ul style="list-style-type: none"> <li>Y-type connecting tool capable of simultaneous hydraulic pressure elevation in two osteotomies</li> </ul>	SNYCT	

# LAS Kit (HLRSNK)



# LAS Kit Surgical Kit Instruments

Dome Drill			
Description	D/Ø	Item Code	Image
<ul style="list-style-type: none"> <li>Dome Drill</li> <li>Forms a bone window and collects autogenous bone simultaneously</li> <li>Excellent penetration due to the combination of macro and micro cutting edges</li> <li>Stopper system safely controls the penetration depth</li> <li>Recommended speed: 1,200~1,500rpm</li> <li>Excessive drilling may cause damage to the membrane</li> </ul>	Ø5.5	LSDR554TD	
	Ø7.0	LSDR74TD	
	Wide Ø7.0	LSDR74WTD	

Core Drill			
Description	D/Ø	Item Code	Image
<ul style="list-style-type: none"> <li>Creates a bone lid while forming the lateral window</li> <li>Cutting edge design based on the CAS drill, enhancing safety</li> <li>Recommended speed: 1,200~1,500rpm</li> <li>Excessive drilling may cause damage to the membrane</li> </ul>	Ø5.5	LSDR554TC	
	Ø7.0	LSDR74TC	

Side Wall Drill						
Description	Item Code					Image
<ul style="list-style-type: none"> <li>Enlarges and trims the rough edges of the bone window</li> <li>Cutting blades start 1 mm above the bottom of the drill</li> <li>Recommended speed: 1,500rpm</li> </ul>	SWDR36T					
Height of side cutting blade (mm)	1.0	2.0	3.0	4.0	.05	
CAS Kit stopper (mm)	8.0	9.0	10	11	12	
<b>Side wall Drill + CAS Kit stopper</b> <ul style="list-style-type: none"> <li>Stoppers safely control the drilling depth</li> </ul>						

Bone Separator		
Description	Item Code	Image
<ul style="list-style-type: none"> <li>Removes the bone lid inside the core drill</li> </ul>	HST75	

Side Wall Drill							
Description	Length	0.5	1.0	1.5	2.0	2.5	3.0
<ul style="list-style-type: none"> <li>Laser marked numbers indicate the drilling depth</li> <li>Color-coded by length</li> <li>Can be used up to 50 times before replacement</li> </ul>							
	Item Code Color	LSNST0.5 Yellow	LSNST1.0 Purple	LSNST1.5 Blue	LSNST2.0 Yellow	LSNST2.5 Purple	LSNST2.0 Blue

# ESSET Kit (HESEK)

For

EKIII

ETIII/IV

SSII/III


Ultra-Wide

Lower panel components

**Torque Wrench**  
TQWCB



**Depth Gauge**  
ODG



# ESSET Kit Surgical Kit Instruments

Crest Remover				
Description	L	Ø5.0	Ø7.0	Image
<ul style="list-style-type: none"> <li>Grinds down narrow alveolar ridge, and creates an indentation for the implant's insertion site</li> <li>Angled type recommended speed: 1,200~1,500rpm</li> <li>Straight type recommended speed: 15,000~30,000rpm</li> </ul>	29	CERM50A	CERM70A	
	45	CERM50S	-	

Twist Drill				
Description	L	TL	Ø1.8	Image
<ul style="list-style-type: none"> <li>Marks the implant's insertion site</li> <li>Controls depth with the built-in stopper</li> <li>Recommended speed: 1,200 ~ 1,500 rpm</li> </ul>	8.5	33	2D1808LC01	
	10	34.5	2D1810LC01	
	11	36	2D1811LC01	

Saw				
Description	T	TL	Item Code	Image
<ul style="list-style-type: none"> <li>For ridge modification and splitting</li> <li>After vertical incision, move from mesial to distal</li> <li>Recommended speed: 1,200 ~ 1,500 rpm</li> <li>Recommended use: 10 times</li> <li>*T = Thickness</li> </ul>	0.3	Ø7.0	HSAW070	
		Ø10	HSAW100	
		Ø13	HSAW130	

Saw Protector			
Description	D	Item Code	Image
<ul style="list-style-type: none"> <li>Semi-circular saw cover protects user when using saws</li> <li>See through protector for maximum procedure visibility</li> <li>360° rotary saw cover for flexible operation</li> <li>Material: Plastic</li> </ul> <p>※ Single use only (Do not reuse)</p>	Ø7.0	HSP070	
	Ø10	HSP100	
	Ø13	HSP130	

# ESSET Kit Surgical Kit Instruments

Expansion Drill					
Description	Type	8.5	10	11.5	Image
<ul style="list-style-type: none"> <li>Gradually expands narrow alveolar ridge</li> <li>Use the Expansion Drills in numerical order based on the diameter of the implant</li> <li>F4.0: I → II → III/F4.5: I → II → III → IV</li> <li>Recommended speed: 25~35rpm</li> </ul>	I Ø1.6/2.8	EXP162808	EXP162810	EXP162811	
	II Ø2.2/3.6	EXP223608	EXP223610	EXP223611	
	III Ø2.8/4.4	EXP284408	EXP284410	EXP284411	
	IV Ø3.2/4.7	EXP324708	EXP324710	EXP324711	

Mount Extension		
Description	Item Code	Image
<ul style="list-style-type: none"> <li>Used to apply manual torque when inserting/removing the Expansion Drills</li> </ul>	HASMEL	

EXP Mount Driver			
Description	Length	Item Code	Image
<ul style="list-style-type: none"> <li>Used to inserting/removing the Expansion Drills with a handpiece and engine</li> </ul>	Short (22.1)	AESMDS	
	Long (28.6)	AESMDL	

Torque Wrench		
Description	Item Code	Image
<ul style="list-style-type: none"> <li>Use to apply torque to the Expansion Drill</li> </ul>	TQWCB	

Depth Gauge		
Description	Item Code	Image
<ul style="list-style-type: none"> <li>Releases excessive torque when hand piece does not move due to being stuck in bone during expansion drill removal process.</li> <li>Use with an open wrench to turn the hex of the Expansion drill</li> <li>Prevents from over torquing</li> </ul>	ODG	




# EM Kit (HMISLK)

Lower panel components

**Depth Gauge**  
HMDTGG



**Ratchet Wrench**  
CITQW-1185A





# EM Kit Surgical Kit Instruments

Drill for EM Implant							
Description	Length	Ø1.5	Ø1.8	Ø2.0	Ø2.3	Ø2.5	Image
<ul style="list-style-type: none"> <li>Laser markings on drill match implant length specifications (8/10/11.5/13/15 mm)</li> <li>In cortical bone, it is recommended to use the lance drill to drill to the final implant length</li> <li>Long type has a stopper at 13mm</li> </ul>	Short (35)	HOSLB15	-	-	-	-	
	Long (38)	HOMSLD15L	-	-	-	-	
	Short (33)	-	HOSMSD18S	HSMSD20S	HOSMSD23S	HOSMSD25S	
	Long (31)	-	HOSMSD18L	HSMSD20L	HOSMSD23L	HOSMSD25L	

Driver for Narrow Ridge & Provisional Type				
Description		Length	Ø3.4	Image
<ul style="list-style-type: none"> <li>Driver for EM (MS) implants: narrow ridge &amp; provisional</li> <li>The triangular marking is used in line with the implant</li> </ul>	Torque Driver	Short (16.5) Long (21.5)	HMSPTS HMSPTL	
	Machine Driver	Short (24.4) Long (29.4)	HMSPMS HMSPML	

Driver for Denture Type				
Description		Length	Ø3.8	Image
<ul style="list-style-type: none"> <li>Driver for EM (MS) Implant denture</li> <li>The triangular marking should be aligned with the implant</li> </ul>	Torque Driver	Short (13.5) Long (18.5)	HMSDTS HMSDTL	
	Machine Driver	Long (21.4)	HMSDMS	

Gauge for MS Implant			
Description	Length	Item code	Image
<ul style="list-style-type: none"> <li>Depth gauge <ul style="list-style-type: none"> <li>Left: to check the drilled depth</li> <li>Right: to bend the neck of the EM (MS) provisional type</li> <li>MS narrow narrow ridge type cannot be bent</li> </ul> </li> <li>Parallel pin is used to confirm the path of the implant after drilling <ul style="list-style-type: none"> <li>MSPP: lower diameter Ø1.5/upper diameter Ø1.8</li> <li>MSPPN: lower diameter Ø1.5/upper shape is same as the abutment portion of the MS narrow ridge</li> </ul> </li> </ul>	Depth Gauge	HMDTGG	
	Parallel Pin	HMSPPP HMSPPN	

Torque Driver Handle	
Description	Image/Item code
<ul style="list-style-type: none"> <li>Used for initial insertion by hand after fastening to torque driver</li> </ul>	 HMSTH

Driver Separator	
Description	Image/Item code
<ul style="list-style-type: none"> <li>If the MS implant is jammed in the Driver, the Driver Separator can be leveraged to separate the two</li> </ul>	 HMSDS


MS Removal Tool		
Description	Image	
<ul style="list-style-type: none"> <li>Easy removal of fractured EM Implant (Narrow Ridge)</li> <li>The tool is used by rotating in the reverse direction after assembling to the universal handle</li> <li>Options based on diameter of fractured implant</li> <li>For Ø 2.0, use orthodontic screw removal tool (code: OSRT20E)</li> </ul> <p>※ Disposable. Do not reuse</p>		
D/Ø	Ø2.5	Ø3.0
	HMRT25E	HMRT30E

Drilling Sequence **EM Drill**

**Narrow Ridge** | **Denture** | **Provisional**


(Length: 10mm)

**Ø1.8mm**




Bone Quality	Lance Drill	Ø1.8 Implant
Soft	▶	Implant Placement
Normal	▶	
Hard	▶	

**Ø2.0mm**




Bone Quality	Lance Drill	Drill (Ø1.8)	Ø2.0 Implant
Soft	▶		Implant Placement
Normal	▶		
Hard	▶	▶	

**Ø2.5mm**





Bone Quality	Drill (Ø1.8)	Drill (Ø2.0)	Ø2.5 Implant
Soft	▶		Implant Placement
Normal	▶		
Hard		▶	



Ø3.0mm





Bone Quality	Drill (Ø1.8)	Drill (Ø2.5)	Ø3.0 Implant
Soft	►		Implant Placement
Normal	►		
Hard		►	


# Surgical Instrument


Abutment Positioning Driver				
Description	H + G/H	Length	Item Code	Image
<ul style="list-style-type: none"> <li>Used for assembling the abutment in the prosthetic process after placing an implant</li> <li>※ For Transfer Abutment only</li> <li>Function to help convenient and stable mounting and tightening of the abutment being pushed away by gingiva</li> <li>Used according to the H and G/H lengths of the abutment to be removed as shown below</li> </ul>	5.0, 6.0, 7.0, 8.0, 9.0	Short (≤9)	HAPDS	
	10, 11, 12, 13, 14	Long (≥10)	HAPDL	


NoMount Torque Driver for ET				
Description	C	Length	Item Code	Image
<ul style="list-style-type: none"> <li>Directly connects to an ET Implant for placement with a Torque Wrench</li> <li>Ensure correct and complete seating before applying torque; loose connection may cause implant fracture</li> <li>C= Connection</li> </ul>	Mini	Short (19)	HGSNMT32S	
		Long (26.6)	HGSNMT32L	
		Ex. Long (33.6)	HGSNMT32E	
	Regular	Short (19)	HGSNMT35S	
		Long (26.6)	HGSNMT35L	
		Ex. Long (33.6)	HGSNMT35E	

NoMount Torque Driver for SS				
Description	C	Length	Item Code	Image
<ul style="list-style-type: none"> <li>Directly connects to a SS Implant for placement with a Torque Wrench</li> <li>Ensure correct and complete seating before applying torque; loose connection may cause implant fracture</li> <li>*C = Connection</li> </ul>	Regular / Wide	Short (16.8)	HSSNMT39S	
		Long (26.8)	HSSNMT39L	

NoMount Driver for EK				
Description	C	Length	Item Code	Image
<ul style="list-style-type: none"> <li>Ø3.5 implant is combined with the bottom of the lower marking</li> <li>Ø4.0, Ø4.5, Ø5.0, Ø6.0 and Ø7.0 implants are combined with the upper part of the lower marking</li> <li>The distance between the two laser marking is 0.5mm</li> <li>C= Connection</li> </ul>	Regular	Short (27.6)	HKSNMDCRS	 <p>Ø4.0, Ø4.5, Ø5.0, Ø6.0, Ø7.0    Ø3.0, 3.5 Implant    Implant</p>
		Long (32.6)	HKSNMDCRL	
		Extra Long (37.6)	HKSNMDCRE	

Implant Driver for ET				
Description	C	Length	Item Code	Image
<ul style="list-style-type: none"> <li>Connects directly to an ET implant for final adjustments to the implant's depth</li> <li>*C = Connection</li> </ul>	Mini	Short (17)	HGSMFDS	
		Long (24)	GSMFDL	
		Ex. Long (34)	HGSMFDE	
	Regular	Short (19)	HGSRFDS	
		Long (26.6)	GSRFDL	
		Ex. Long (33.6)	HGSRFDE	

Implant Driver for SS				
Description	C	Length	Item Code	Image
<ul style="list-style-type: none"> <li>Connects directly to a SS implant for final adjustments to the implant depth</li> <li>*C = Connection</li> </ul>	Regular/ Wide	Short (14)	HSSRFDS	
	Regular/ Wide	Long (24)	SSRFDL	
	Regular/ Wide	Ex. Long (31)	HSSRFDE	

Implant Driver for EK				
Description	C	Length	Item Code	Image
<ul style="list-style-type: none"> <li>Connects directly to an EK implant for final adjustments to the implant's depth</li> <li>*C = Connection</li> </ul>	Regular	Short	HKSFDS	
		Long	HKSFDL	

Torque Extension				
Description			Item Code	Image
<ul style="list-style-type: none"> <li>Extends the length of an instrument by 10mm</li> <li>Connects to the torque wrench</li> </ul>			HOTE	

# Surgical Instrument Kit

Simple Mount Driver			
Description	Length	Item Code	Image
<ul style="list-style-type: none"> <li>Connects to mounted implants for placement</li> <li>For use with a handpiece/implant motor</li> </ul>	Short (20.1)	HASMDS	
	Long (26.5)	HASMDL	


Simple Mount Extension			
Description	Length	Item Code	Image
<ul style="list-style-type: none"> <li>Extends the length of the simple mount driver and converts it for use with the Torque Wrench</li> </ul>	Short (14.5)	HASMES	
	Long (20.5)	HASMEL	

Simple Open Wrench		
Description	Item Code	Image
<ul style="list-style-type: none"> <li>Disengages the simple mount when bone quality is poor</li> <li>Easy insertion into the mouth with a neck angle of 30°</li> </ul>	SPOW	


Removal Tool for Implant Mount			
Description	Application	Item Code	Image
<ul style="list-style-type: none"> <li>Removes the mount if the mount becomes wedged in the implant</li> <li>Used with Driver Handle or Torque Wrench</li> <li>Removes the mount screw, insert the Removal Tool into the mount, and turn clockwise</li> <li>App = Application</li> </ul>	Mini (ET,US)	HERFM	
	Regular (ET,US) Wide (SS)	HHRFR	
	Wide (US)	HERFW	


Depth Gauge		
Description	Item Code	Image
<ul style="list-style-type: none"> <li>Measures drilling depth (7~15mm)</li> </ul>	ODG	


Positioning Guide			
Description	W/L	Item Code	Image
<ul style="list-style-type: none"> <li>Help set the drilling interval for implant insertion</li> <li>Insert after initial drilling</li> </ul>	2.5 / 21.5	HAPG201	
	6.0 / 17.5	HAPG202	
	11 / 17.5	HAPG203	

Tissue Height Gauge for ET		
Description	Item Code	Image
<ul style="list-style-type: none"> <li>Connects to the ET implant to measure the height of the gingiva in relation to the implant</li> </ul>	HGTSHG	


Ratchet Wrench		
Description	Item Code	Image
<ul style="list-style-type: none"> <li>It prevents wrench from backdriving</li> <li>Excessive torquing may cause damage to the bone or the inside of a implant</li> </ul>	CITQW-1185A	


L-Wrench		
Description	Item Code	Image
<ul style="list-style-type: none"> <li>1.2 hex driver for hard to reach areas like narrow intermaxillary areas</li> <li>Torque indication: when the wrench starts to bend (around 10°), it is possible to apply 5~8Ncm of torque</li> </ul>	HLWC	


Torque Wrench - Bar Type		
Description	Item Code	Image
<ul style="list-style-type: none"> <li>Adjusts the implant depth, and tightens abutments, screws, etc.</li> <li>Pull the bar back until the desired torque value is reached</li> </ul>	TQWCB	


Torque Wrench - Spring Type		
Description	Item Code	Image
<ul style="list-style-type: none"> <li>Applies a precise amount of torque (10/20/30Ncm) to the screw and abutment</li> <li>The neck of the torque wrench will bend when the exact amount of torque has been delivered</li> <li>Do not continue to torque after the neck has bent; excessive force may cause screw fracture etc.</li> </ul>	HTW30	


# Surgical Instrument Kit

Torque Connector		
Description	Item Code	Image
<ul style="list-style-type: none"> <li>Tool used convert a square driver connection to a bi-directional connection for the torque wrench</li> </ul>	HRC	


Machine Driver Connector		
Description	Item Code	Image
<ul style="list-style-type: none"> <li>Tool used convert a machine driver into a bi-directional connection for the torque wrench</li> </ul>	HMDC	


Driver Handle		
Description	Item Code	Image
<ul style="list-style-type: none"> <li>Connects to the Torque Driver</li> </ul>	HTIDHC	

Machine Driver Handle		
Description	Item Code	Image/Guide
<ul style="list-style-type: none"> <li>Tool used convert engine type surgical tools into a manual type</li> </ul>	HMDH	

Bone Mill		
Description	Item Code	Image
<ul style="list-style-type: none"> <li>Grinds harvested autogenous bone</li> </ul>	HABM	



Anterior Hand Driver for Implant		
Description	Item Code	Image
<ul style="list-style-type: none"> <li>Manually torque implants in the anterior area</li> <li>Connect to a NoMount torque driver or a implant driver</li> <li>Excessive torque may cause damage to the implant and/or driver</li> </ul>	HAHDI	

Torque Handle		
Description	Item Code	Image
<ul style="list-style-type: none"> <li>Connectable to a contra-angle handpiece (Hand-Piece gear ratio to 1:1)</li> <li>Used to connect healing abutments, cover screws, abutment screws, orthodontic screws, etc. (Note: after connecting, make sure that it is tightened with a torque wrench)</li> <li>Excessive torque may cause damage to the screw and/or hand piece</li> </ul>	HTQHD	





**HIOSSEN**  
IMPLANT

**Smiles that last a life time**