

Super Convenient Procedure in 2 Steps

# 122 Taper Kit

- Surgical KIT That is Optimized to The Taper Implant
- Reducing Chair Time with Ever Convenient Procedure
- Increasing User Convenience as Cortical Drilling is Not Necessary

Use for **ETIII /IV**


Normal Bone  
**Done in Two Steps!**



**HROSSEN**  
IMPLANT



## 1 Soft Bone Done in One Step




## Excellent Drilling Power

Maximizing drilling power (with distinguishing drill design and excellent manufacturing technology) which skips intermediate drilling and provides smooth and fast drilling result




## 2 Normal Bone Done in Two Steps



## Preventing Bone Heating

Excellent drilling power leads to skipping intermediate drilling and this prevents bone heating by maintaining low drilling temperature




122 Taper Drill

Taper Drill

Drilling temperature when using

## 2 Hard Bone Done in Two Steps



## Preventing Drill Bouncing

Multiple shaped taper drill also functions as pilot drill and this prevents bouncing during drilling

