



Crestal Approach Sinus Lift Augmentation
without Fear of Membrane Damage

CAS Kit


- Prevention of maxillary sinus membrane perforation
- Easy and secure with hydraulic pressure
- Applicable to various sinus floor morphologies



HROSSEN
IMPLANT




The formation of an inverted conical bone lid provides puncture prevention




The CAS drill is round in shape and provides safe drilling

Prevents Perforation of Sinus Membrane

- Applying a hydraulic lift to the cavity gently separates the membrane from the sinus floor
- Provides significant comfort to patient





Osteotome KIT



CAS KIT

Easy and Safe Membrane Elevation with Hydraulic Pressure

- Hydraulic pressure lifts membrane high and wide
- Elevated sinuses have excellent blood supply and osteogenesis capability



Applicable in Various Cases

