


Normal Bone  
**Just in  
Two Steps!**




Convenient procedure in only 2 Steps

# 122 Taper Kit

- Optimized surgical kit for Hiossen's ETIII & ETIV implant system
- Reduce chair time with just 2 drilling procedures
- Increase user convenience by eliminating the cortical drilling step

Use for **EKIII/ETIII/ETIV**


**HIOSSSEN**  
IMPLANT



**1 Soft Bone**  
Just in One Step

**Excellent Drilling Power**


Maximized cutting power with a distinguished drill design and innovative manufacturing technology that skips intermediate drilling and provides a smooth and fast drilling result.



**2 Normal Bone**  
Just in Two Steps

**Preventing Bone Heating**


Excellent cutting power leads to skipping cortical drilling which prevents bone heating by maintaining low drilling temperature



122 Taper Drill

Taper Drill

Drill temperature comparison between the 122 taper drills and other taper drills



**2 Hard Bone**  
Just in Two Steps

**Preventing Drill Bouncing**

Multi-level tapered drill also functions as a pilot drill and prevents bouncing during drilling

