

HIOSSEN DENTAL EQUIPMENT

PRODUCT CATALOG



Hiossen
IMPLANT



Hiossen Dental Equipment

Version: PC24HDELTR1.0

270 Sylvan Ave. Ste 1130, Englewood Cliffs, NJ 07632
www.hiossen.com Email: marketing@hiossen.com

@Hiossen @HiossenImplants @Hiossen



Hiossen Dental Equipment



Intro

About us	004
History	005

Hiossen Dental Equipment


SM3	006
SM5	007
IS4™	008
IS3™	009
E-Drive / Pluto / FireCR / SNAP	014
Medit Intraoral Scanner	015
Trios® 5 Intraoral Scanner	016
Trios® 4 Intraoral Scanner	017
Trios® 3 Intraoral Scanner	018
Trios® Move+	019
K5 Dental Unit Chair	020
K3 Dental Unit Chair	021
T2 Plus CBCT	022
T2 CBCT	023


Legacy Based on Technology

Every Hiossen implant undergoes domestic production, manufacturing, and packaging in Pennsylvania, USA. Through collaboration with our international distribution network, we have established a stringent quality management system, encompassing thorough inspection, quality tracking, in-house surface treatment, and sterilization processes. Hiossen maintains stringent quality certifications, including FDA and CE approvals. Exporting implants to over 40 countries worldwide, Hiossen is renowned in global markets for its cutting-edge technology and unwavering commitment to quality.

Continuously striving to offer cutting-edge technology, Hiossen is on a trajectory to rank among the foremost global companies in implant, restorative, and digital dentistry. Our commitment to reliability, convenience, and simplicity is underpinned by relentless inquiry, perseverance, and confidence. In partnership with our parent company, Osstem Implant, we have ascended to become one of the top three leaders in dental implants, offering a broad array of products and services.

Hiossen is driven by a shared passion for enhancing lives through the refinement and development of innovative products. We aim to streamline medical procedures, making them safer and more straight forward for both clinicians and patients alike. Our team is dedicated to achieving optimal surgical and restorative outcomes through meticulous testing, research, clinical studies, and services. This commitment is underscored by our unwavering dedication to enhancing the quality of life for edentulous patients.






SM3 Implant Engine

SM3 Unit	
Description	Image
<ul style="list-style-type: none"> • High Precision Calibration • Offers Advanced Technology • Compact Body and Large LCD Display • LED Optics for Safe, More Accurate Treatment • Well-Balanced, compact and light weight micro-motor 	
SM3US1	

Component: SM3 Contra-angle Hand Piece	
Description/Item code	Image
<ul style="list-style-type: none"> • Product Name: SM3 Contra-angle Hand Piece • Body Material: Titanium • Transmission: Gear Ratio 20:1 Reduction • Maximum Speed: 2,000 min-1 • Maximum Torque: 80 Ncm 	
S200EL	


Component : Irrigation Tube	
Description/Item code	Image
<ul style="list-style-type: none"> • SM3 Irrigation Tube (10 packs) - Disposable 	
Y900114	


SM5 Implant Engine

SM5 Unit	
Description	Image
<ul style="list-style-type: none"> Intelligent Power and Precise Torque One Touch Auto-Calibration Simple User Friendly Operation Enhanced workflow through easy programming 	
SM5US1	


Component: Hand Piece SURGmatic S201L (Optic)	
Description/Item code	Image
<ul style="list-style-type: none"> Product Name : S201L Reduction ratio 20:1 Removable head, easy to clean Small head with internal cooling system Speed: 15 – 2,000 rpm Maximum torque: 55 Ncm 	
1.012.1870	

Component: Hand Piece SURGmatic S11L	
Description/Item code	Image
<ul style="list-style-type: none"> Product Name : S11L Transmission : 1:1 Maximum Speed : 40,000rpm Maximum Torque : 5.5Ncm For 44.5 mm standard burs (can also be used with contra-angle burs) 	
1.012.1872	

Component: Irrigation Tube	
Description/Item code	Image
<ul style="list-style-type: none"> Tube Set Sterile S600 (10EA) - Disposable 	
1.009.8757	

Component: Irrigation Tube	
Description/Item code	Image
<ul style="list-style-type: none"> Tube Set Multiuse S600 (2EA) - Reusable 	
1.011.0633	

ISQ Stability Tester


IS4™	
Description	Image
<ul style="list-style-type: none"> • Precision Resonance Frequency Analysis (RFA) • Assess implant stability with the Implant Stability Quotient (ISQ) • Enhance implant success through accurate diagnostics • Achieve high-accuracy measurements for reliable results • Optimize loading protocols with precise data • User-Friendly design with single key operation 	
55250-HM	


Main Adapter
Image


55263


Main Adapter
Image

55267


IS4 Charging Station
Image

55245

IS4 Battery Replacement Kit
Image

55291

Multipeg™ Driver
Image

55245

Multipeg™	
Description/Item code	Image
<ul style="list-style-type: none"> • Up to 20 times without any deviation • Reusable, auto-clavable, calibrated, and compatible with all major implant systems • Optical Platform Fit 	

ISQ Stability Tester


IS3™	
Description	Image
<ul style="list-style-type: none"> • Resonance Frequency Analysis (RFA) • Implant Stability Quotient (ISQ) • Raise the implant success rate with accurate diagnosis • High accuracy measurement result • Determines optimal loading protocol • Hand held instrument with single key operation 	
55002-HM	


Main Adapter
Image

55093

Main Adapter
Image

55097

MulTipeg™ Driver
Image

55003

MulTipeg™	
Description/Item code	Image
<ul style="list-style-type: none"> • Up to 20 times without any deviation • Reusable, auto-clavable, calibrated, and compatible with all major implant systems • Optical Platform Fit 	

MulTipeg™ Assortment

		No	Type			No	Type			No	Type
AB Dental				Bio3				Dentalis			
Narrow Platform		55021	11	Advanced		55067	59	Regular		55051	43
Standard Platform		55049	41	Progressive		55067	59				
ADIN				BioHorizons				Dental Ratio			
Touareg S/OS		55034	26	Laser-Lok 3.0 Grey		55034	26	Tissue Level RP, WP		55014	4
Touareg Closefit UNP 2.75		55059	51	Internal Laser-Lok 3.5 Yellow		55049	41	Bone Level NC		55061	53
Touareg Closefit NP 3.0		55059	51	Internal Laser-Lok 4.5 Green		55051	43	Bone Level RC		55062	54
Touareg Closefit RP 3.5		55046	38	Internal Laser-Lok 5.7 Blue		55051	43				
Touareg CloseFit WP 4.3, 5.0		55067	59	External 4.0 Platform		55011	1	Dental Tech			
Swell 3.3		55034	26					Implassic FT3 3.25		55039	31
AIRES				Biomet/3i				Implassic FT3 3.75		55049	41
3.3		55049	41	External Hex		55011	1	Dentegris			
Alpha-Bio				Certain Internal Purple		55036	28	3.3		55034	26
Internal Hex		55034	26	Certain internal Purple Platform 3.4		55036	28	3.75		55057	49
Conical Hex 3.2		55065	57	Certain internal Blue, Yellow, Green		55036	28	4.5		55049	41
Conical Hex 3.75-5.0		55051	43	T3 Purple		55036	28	3.8/4.5		55049	41
Alfa Gate				T3 Blue, Yellow, Green		55066	58	5.5		55020	10
Bioactive, All diam		55034	26	Bionnovation				Dentis			
MAX, All diam		55034	26	Biomorse		55023	14	Cleanlant S-Clean 3.7 - 4.8		55020	10
M+ NP 3.5		55044	36	Biotech Dental				S-Clean OneQ Narrow 3.3		55064	56
M+ RP 4.3 - 5.0		55037	29	Kontakt, Kontakt-S 3.6 - 4.8		55045	37	Dentium			
Porous, All diam		55034	26	BlueSkyBio				Superline		55037	29
Alpha Dent				BIO Max NP Platform, All		55046	38	NR Line		55066	58
Active Conus		55067	59	BIO Max RP Platform 4.3-5.0		55037	29	Simple Line II GH 1,6		55014	4
Classic Conus		55067	59	BIO Int Hex 3.5 Platform 3.7-4.1		55049	41	DEVA			
Active Bio		55067	59	BIO Int Hex 3.5 Platform 4.7		55051	43	Bone Level 3.5		55064	56
Slim		55039	31	BIO Int Hex 4.5 Platform 4.7-7.0		55051	43	Bone Level 4.5		55051	43
Anker				BIO One-stage Narrow, Regular, Wide		55014	4	DIO			
SB III 3.5		55023	14	BIO One-stage Short Regular & Wide		55014	4	SM Submerged 3.8		55049	41
SB III 4.0-5.0		55020	10	BIO Quattro Narrow 3.3		55061	53	SM Submerged 5.0		55020	10
Ankylos				BIO Quattro Regular 4.1 - 5.6		55062	54	UF II Narrow		55064	56
C/X		55025	16	BIO Conus 12, 3.0 Platform 3.0		55064	56	UF II Regular		55020	10
Anthogyr				BIO Conus 12, 3.5/4.0 Platform 3.5-4.0		55046	38	Easy Implant			
Axiom Bone Level 3.4 - 5.2		55065	57	BIO Conus 12, 4.5 x 6mm		55046	38	Ocean		55046	38
Axiom Tri-lobe		55065	57	BIO Conus 12, 4.5/5.0 Platform 4.5 - 5.0		55017	7	Lilas		55037	29
Axiom Tissue Level Narrow, Regular		55050	42	BIO Conus 12, 5.0 x 6mm		55046	38	Edierre			
Argon				BIO Trilobe 3.5		55016	6	3-lobe		55021	11
K3pro Rapid 3.0 - 3.5		55026	17	BIO Trilobe 4.3-6.0		55017	7	4-lobe		55025	16
K3pro Rapid Platf 3mm		55019	9	Bredent				Euroteknika			
K3pro Rapid Platf 3mm Short		55044	36	SKY		55025	16	Internal Platform 3.0		55064	56
Astra Tech				CopaSKY		55027	19	Internal Platform 3.5		55081	73
Osseospeed EV 3.0		55058	50	Mini2SKY		55033	25	Naturactis/Naturall+ 3.0		55064	56
Osseospeed EV 3.6		55046	38	Fast & fixed Abutment		55083	75	Naturactis 3.5-5.0		55081	73
Osseospeed EV 4.2		55057	49	BTI				Uneva		55011	1
Osseospeed EV 4.8		55055	47	Interna Universal		55051	43	Galimplant			
Osseospeed EV 5.4		55048	40	Interna Wide		55051	43	iPx 4012		55046	38
Osseospeed TX 3.0 S		55064	56	BTK				GC Tech			
Osseospeed TX 3.5-4.0 S		55046	38	External Hex EN		55038	30	Aadva Narrow		55081	73
Osseospeed TX 4.5-5.0 S		55017	7	External Hex ER		55011	1	Aadva Regular		55066	58
Avinent				External Hex EW		55011	1	Geass			
External, Coral & Ocean 3.5		55038	30	Internal Hex IR		55049	41	Way Extra		55081	73
External, Coral & Ocean 4.1		55011	1	Internal Hex IM		55051	43	Global D			
External, Coral 5.1		55013	3	Internal Hex IW		55051	43	In-Kone 4.0-5.0		55024	15
Internal, Coral & Ocean 3.5, 4.1		55051	43	Hexagonal Morse Taper KR		55081	73	Geass			
Conical, Ocean 3.5		55066	58	Hexagonal Morse Taper KW		55056	48	Way Extra		55081	73
Conical, Ocean 4.1		55004	18	Octagonal Morse Taper SR		55014	4	Global D			
B&B Implant				Camlog				GMI			
Duravit Slim		55021	11	Screw-Line, Root-Line 3.3		55036	28	Avantgard RP 3.75-4.25		55066	58
Duravit 3P		55035	27	Screw-Line, Root-Line 3.8, 4.3		55031	23	Avantgard WP 4.75-5.75		55075	67
Duravit Evolution		55035	27	Cowellmedi				Frontier RP/WP		55051	43
BEGO				Inno Int Octa		55014	4	Phoenix NP 3.3		55038	30
Semados Mini		55012	2	Inno Sub Hex		55037	29	Phoenix Std 3.3 - Wide 5.0		55011	1
Semados S / RI 3.25-5.5		55034	26	C-Tech				Monolith 3.0		55043	35
Semados RS / RSX 3.0		55082	74	EL		55067	59	Dentack			
Semados RS / RSX 3.75-5.5		55034	26	Dentack				Internal Hex		55056	48
Semados SC / SCX 3.25-5.5		55034	26	Internal Hex				Dentalis			
				Regular		55051	43	Regular			

MulTipeg™ Assortment

	No	Type		No	Type		No	Type
GP Implants			Implantvision			Medentis		
SSI Spiral Shape implant	55051	43	Internal Hex 3.7-5.7	55034	26	ICX Premium 3.75 - 4.8	55062	54
GT Medical			Innobiosurg			ICX Plus 3.45	55039	31
HXE	55011	1	MagiCore Regular	55018	8	ICX Narrow	55058	50
HXE RP	55012	2	Intoss Anchor			Medical Instinct		
HXE WP	55013	3	Internal Hex 3.7-5.7	55088	80	Hex 3.0	55071	63
HXI 4.1	55037	29				Hex 3.4, 4.0, 5.0	55052	44
002 HXI 4.7	55049	41	Intra-lock			Cone 3.0, 3.4, 4.0	55071	63
OCT 4.3	55014	4	Intrahex 3.4 - 4.5	55051	43	Cone 5.0	55072	64
			Gold & Blue 3.4 - 6.0	55051	43			
Hahn			InDex 3.4 - 3.75	55015	5	Medigma Biomedical		
Tapered 3.0	55058	50	InDex 4.0 - 4.3	55020	10	Mars 3,5-6,0	55057	49
Tapered 3.5	55046	38	FlatTop 3.75	55012	2	FixTite 3,5-6	55057	49
Tapered 4.3	55046	38				Fix-a-Dent 3,75-6	55057	49
Tapered 5.0	55037	29	IRES					
Tapered 7.0	55048	40	External Hex Narrow 3.3	55012	2	Medimecca		
			External Hex Regular 3.7 - 5.2	55011	1	PT SLA Chaorum	55047	39
Hiossen			Internal Hex Narrow 3.3	55087	79			
HS System	55028	20	Internal Hex Regular 3.7 - 5.2	55049	41	Megagen		
ET System	55037	29	Internal Octagon Regular 3.7 - 5.2	55014	4	Anyridge	55051	43
ET 3.2	55047	39				Anyone Internal	55037	29
ET mini	55015	5	J Dental Care			Mini	55064	56
EK 3.2/3.5	55046	38	JDEvolution 3.7-6.0	55067	59	MUA	55080	72
EK 4.0	55044	36	JD Evolution S 3.25	55044	36			
			JD ICON ultra, S 2.75	55064	56	Microdent		
			Conical Abutment	55080	72	Universal	55011	1
Hi-Tec			Keystone			Microdent 2.8 - 3.25	55068	60
Tapered Self Thread NTI 3.3, SLT 3.75	55034	26	Prima / Connex Narrow	55049	41	Microdent 3.3 - 5.0	55011	1
Tapered Self Thread GTI 4.2, WTI 5.0	55051	43	Prima / Connex Regular	55067	59	Microdent 5.5 - 8.0	55013	3
Spark TPI	55034	26	Prima / Connex Wide	55041	33	Genius 3.0	55025	16
Logic Plus	55034	26	Prima / Genesis Narrow	55049	41	Genius 3.5 - 5.0	55051	43
Logic LGI 3.0-3.5	55046	38	Prima / Genesis Regular	55067	59	Ektos 3.2	55023	14
Logic LGI 4.3-5.0	55037	29	Prima / Genesis Wide	55041	33	Ektos 3.7 - 8.0	55051	43
			Restore Narrow	55079	71	Trylogic 3.5 - 4.2	55051	43
IDI Evolution			Restore Regular	55011	1	Trylogic 5.0	55075	67
Standard implant 3.75-4.0	55004	18	Restore Wide	55013	3			
B1One High	55080	72				MIS		
B1One Low	55033	25	Klockner			Conical connection Narrow (C1, V3)	55081	73
			VEGA Micro 3.0	55074	66	Conical connection Standard (C1, V3)	55087	79
IDI Implants			VEGA Narrow 3.5	55065	57	Conical connection Wide (C1)	55041	33
ID Cam ST 4.2	55010	13	VEGA Regular 4.0, 4.5	55034	26	Seven and M4 Narrow	55027	19
ID All	55010	13	KL Narrow 3.3	55012	2	Seven and M4 Standard	55049	41
IDS			KL Regular 3.7, 4.2	55011	1	Seven Wide 5.0	55051	43
Reflect Aspire	55081	73	KL Wide 4.7	55013	3	Connect abutment NP, SP, WP	55115	81
Reflect Certus	55019	9	SK2	55038	30			
Reflect Tapered	55034	26	Essential Cone	55067	59	Neo Biotech		
Reflect Recover	55020	10	S3M 3.1	55079	71	IS II 4.0	55020	10
Reflect Rapid	55020	10				IS III 4.0	55020	10
Implance			Konmet			Neodent		
Tissue level	55014	4	Classic 4.0, 4.8	55014	4	Drive CM	55025	16
Bone level Regular	55037	29	Cylindrical 3.5-4.0	55065	57	Titamax WS Cortical	55056	48
Bone level 3.3	55065	57	Cylindrical 4.5-5.5	55057	49	GM	55066	58
			Short Cylindrical 5.0-6.0	55078	70			
Implant One			Conical 4.4/3.4	55065	57	Neodent		
100 Series Std thread	55029	21	Conical 4.8/3.8, 5.5/4.5	55057	49	Drive CM	55025	16
200 Series Std thread	55029	21				Titamax WS Cortical	55056	48
300 Series Std thread (blue)	55029	21	Lasak			GM	55066	58
400 Series Std thread (purple)	55035	27	BioniQ QR	55066	58			
500 series Std thread (green)	55035	27	BioniQ QN	55021	11	Neoss		
300 Series Wide thread (blue)	55029	21				SP (Ø 3.5 implants and wider)	55004	18
400 Series Wide thread (purple)	55035	27	Leader			NP (Ø 3.25)	55027	19
500 Series Wide thread (green)	55035	27	Internal Hex, 3.3 - 5.5 mm	55034	26	Access Abutment	55060	52
			Cono-In	55067	59			
Implant Direct			Maxillent			Nobel Biocare		
Legacy 1/2/3/4 3.0/3.2	55012	2	iRaise	55049	41	Brånemark System Ext hex RP	55011	1
Legacy 1/2/3/4 3.5/3.7 - 5.7/7.0	55051	43				Brånemark System Ext hex NP	55012	2
RePlant 3.5	55022	12	Medentika			Brånemark System Ext hex WP	55013	3
RePlant 4.3 - 5.0	55010	13	Microcone NI	55021	11	Tri-lobe Connection NP	55022	12
RePlant 6.0	55032	24	Microcone RI	55081	73	Tri-lobe Connection RP, WP	55010	13
ReActive 3.5 Platform	55022	12	Quattrocone RI	55081	73	Tri-lobe Connection 6.0	55032	24
ReActive 4.3 - 5.0 Platform	55010	13	Quattrocone 30 Al	55081	73	Conical Connection 3.0	55058	50
Interactive 3.0 Platform	55046	38	Procone 3.3	55119	85	Conical Connection Narrow	55046	38
Interactive 3.4 Platform	55037	29	Procone 3.8	55031	23	Conical Connection Regular	55037	29
Screwindirect 5.0 PF 3.0 - 5.7	55070	62	Procone 4.3	55031	23	Conical Connection Wide 5.5	55048	40
			Procone 5.0	55032	24	MUA NP & RP	55080	72
Implant Swiss						On1 NP, RP, WP	55042	34
Tissue level	55014	4						
Bone level Regular	55037	29						
Bone level 3.3	55065	57						

MulTipeg™ Assortment

	No	Type		No	Type		No	Type
Noris Medical			Quest			Sweden & Martina		
Tuff 3.3 - 6.0	55034	26	Hylock NP	55030	22	Syra, Syra SP 3.9-5.0	55011	1
			Hylock RP	55066	58	Outlink2 3.75, 4.10, 5.0	55011	1
Novem			Resista			Shelta 3.8-5.0	55025	16
3.4	55039	31	WIS 3.0-5.0	55049	41	Premium Straight & SP 3.3-3.8	55025	16
3.8	55016	6				Kohno Straight & SP 3.3-3.8	55025	16
4.2, 5.0	55057	49	Rex Implants			Premium Straight & SP 4.25-6.0	55004	18
NucleOss			External hex TL 1.8	55060	52	Kohno Straight & SP 4.25-6.0	55004	18
Tpure Bone Level	55004	18	Internal hex TL 2.9	55034	26	Prama 3.8 - 5.0	55079	71
T6 Bone Level	55074	66				DAT Connection	55035	27
OCO Biomedical			Ritter			Tag Dental		
Engage 3.25	55027	19	Spiral SB/LA Narrow PF 2.9	55081	73	Axis	55051	43
Engage 4.0-6.0	55051	43	Spiral SB/LA Standard PF 3.75	55049	41	massif	55051	43
ERI 3.25	55069	61	RS Implant			Tatum		
ERI 4.0-5.0	55014	4	3.5 - 6.0	55020	10	Integrity 3.7 - 5.0	55015	5
TSI 3.25	55068	60				Integrity 6.0 - 8.0	55040	32
TSI 4.0-5.0	55011	1	Schutz Dental			T implant	55063	55
OsseoFuse			Impla cylindrical	55052	44	Simply Implants		
Hexaplus, Sinusplus	55081	73	IMPLA micro Retention Cone 4.2	55052	44	Internal Hex	55051	43
Osstem			SGS Implant			SIN		
US Regular	55011	1	P1, P7 Platform	55049	41	Strong SW CM	55023	14
US Mini	55012	2	P1D, P7D	55049	41	Strong SW Int Hex	55035	27
US Wide	55013	3	P7N	55059	51	Unitite CM 2.9	55084	76
US Wide PS	55011	1	SIC invent			Unitite CM 3.5	55077	69
SS Regular, Wide	55028	20	SICace 3.4 - 5.0	55052	44	Unitite CM Compact 4.0	55085	77
TS Regular, Reg Short, Diam 6 & 7	55037	29	SICmax, SICtapered 3.7 - 5.2	55052	44	Tryon 3.25 - 5.0	55011	1
TS Mini 3.5	55015	5	SICvantage® max, SICvantage tapered 3.0	55026	17	Southern Implants		
TS Mini 3.0	55047	39	SICvantage® max, SICvantage tapered 3.7	55116	82	Tri-Nex 4.3-5.0	55010	13
			SICvantage® max, SICvantage tapered 4.2-5.2	55054	46	Tri-Nex 3.5	55022	12
Osteo Plus			MUA	55120	86	External Hex	55011	1
Iridium PF 3.8 & 4.2	55079	71	Simply Implants			MAX Ext Hex	55011	1
She 3.45, 3.75	55011	1	Internal Hex	55051	43	Co-Axis Ext Hex	55011	1
Shi 3.75, 4.2	55049	41				Piccolo Ext Hex	55050	42
Shi 5.0	55051	43	SIN			IT-Connection	55014	4
OT-Medical			Strong SW CM	55023	14	Deep Conical 3.0	55064	56
OT-F2 3.4, 3.8	55053	45	Strong SW Int Hex	55035	27	Deep Conical 3.5-4.0	55046	38
OT-F2 4.1, 5.0	55073	65	Unitite CM 2.9	55084	76	Deep Conical Tapered 5.0	55017	7
Overmed			Unitite CM 3.5	55077	69	Spiral Tech		
Easyline Easydip	55020	10	Unitite CM Compact 4.0	55085	77	ESi, Hex connection	55051	43
			Tryon 3.25 - 5.0	55011	1	Ultimate, Hex connection	55051	43
Oxy Implant			Southern Implants			ESi Conical 3.5	55044	36
Internal	55020	10	Tri-Nex 4.3-5.0	55010	13	ESi Conical 4.3, 5.0	55037	29
Kone Conical Narrow	55021	11	Tri-Nex 3.5	55022	12	Straumann		
Kone Conical Regular	55051	43	External Hex	55011	1	Tissue level, Regular & Wide Neck	55014	4
PSK Line Regular	55051	43	MAX Ext Hex	55011	1	Tissue level, Narrow Neck Crossfit (NNC)	55076	68
External Mini	55079	71	Co-Axis Ext Hex	55011	1	Bone level / Bone Level Tapered Narrow	55061	53
External Regular	55011	1	Piccolo Ext Hex	55050	42	Bone level / Bone Level Tapered Regular	55062	54
External Maxi	55011	1	IT-Connection	55014	4	Bone level / Bone Level Tapered SC 2.9	55064	56
Paltop			Deep Conical 3.0	55064	56	Bone Level X (BLX) all platform	55027	19
Standard 3.75-5.0	55034	26	Deep Conical 3.5-4.0	55046	38	Screw retained Abutment NP, RP	55080	72
Narrow	55027	19	Deep Conical Tapered 5.0	55017	7	Surgikor		
PCA Conical	55066	58	Spiral Tech			Fixation 3.0	55058	50
Phibo			ESi, Hex connection	55051	43	Fixation 3.5	55066	58
Aurea EVO NP	55027	19	Ultimate, Hex connection	55051	43	Fixation 4.3	55037	29
Aurea EVO RP, WP	55051	43	ESi Conical 3.5	55044	36	Solution 3.5, 4.0	55066	58
TSA S3	55074	66	ESi Conical 4.3, 5.0	55037	29	Solution 4.5, 6.0	55037	29
TSA S4, S5	55034	26	Straumann			Immediate, all	55049	41
TSH S3, S4, S5	55011	1	Tissue level, Regular & Wide Neck	55014	4	Versatile, all	55049	41
P-I Brånemark			Tissue level, Narrow Neck Crossfit (NNC)	55076	68	Sweden & Martina		
HEX 3.3	55115	81	Bone level / Bone Level Tapered Narrow	55061	53	Syra, Syra SP 3.9-5.0	55011	1
HEX 3.75 -5.0	55117	83	Bone level / Bone Level Tapered Regular	55062	54	Outlink2 3.75, 4.10, 5.0	55011	1
AMP 3.3	55039	31	Bone level / Bone Level Tapered SC 2.9	55064	56	Shelta 3.8-5.0	55025	16
AMP 3.75 - 4.8	55118	84	Bone Level X (BLX) all platform	55027	19	Premium Straight & SP 3.3-3.8	55025	16
MT 3.3 - 4.8	55046	38	Screw retained Abutment NP, RP	55080	72	Kohno Straight & SP 3.3-3.8	55025	16
Prodent Italia			Surgikor			Premium Straight & SP 4.25-6.0	55004	18
Twinner & Prime, Conomet, Collar: Orange	55027	19	Fixation 3.0	55058	50	Kohno Straight & SP 4.25-6.0	55004	18
Twinner & Prime, Conomet, Collar: Fuchsia	55051	43	Fixation 3.5	55066	58	Prama 3.8 - 5.0	55079	71
Twinner & Prime, Conomet, Collar: Green	55067	59	Fixation 4.3	55037	29	DAT Connection	55035	27
Twinner & Prime, Conomet, Collar: Yellow	55020	10	Solution 3.5, 4.0	55066	58			
Prime, Conomet, Collar: Blue	55020	10	Solution 4.5, 6.0	55037	29			
			Immediate, all	55049	41			
			Versatile, all	55049	41			

MulTipeg™ Assortment

	No	Type		No	Type		No	Type
Sweden & Martina			Ticare			Vulkan		
Syra, Syra SP 3.9-5.0	55011	1	Osseus Mini, Standard	55011	1	Internal Hex	55034	26
Outlink2 3.75, 4.10, 5.0	55011	1	Osseus Maxi	55013	3	External Hex NP	55012	2
Shelta 3.8-5.0	55025	16	Osseus Maxi ps	55011	1	External Hex RP	55011	1
Premium Straight & SP 3.3-3.8	55025	16	Inhex Mini	55021	11	Conical con MIP	55058	50
Kohno Straight & SP 3.3-3.8	55025	16	Inhex Standard	55046	38	Conical con NP	55046	38
Premium Straight & SP 4.25-6.0	55004	18	Inhex Maxi	55017	7			
Kohno Straight & SP 4.25-6.0	55004	18	Osseus Quattro Standard	55011	1	Winsix		
Prama 3.8 - 5.0	55079	71	Inhex Quattro Standard	55046	38	Internal Hex	55020	10
DAT Connection	55035	27						
			TRATE			Xgate		
Tag Dental			ROOTFORM	55065	57	Conical 3.3 - 5.0	55034	26
Axis	55051	43						
Massif	55051	43	TRI Dental implants			Xive		
			Tissue level Octa	55014	4	S-Plus 3.0	55030	22
Tatum			Tri-Vent, all	55034	26	S-Plus 3.4	55061	53
Integrity 3.7 - 5.0	55015	5	Narrow	55025	16	S-Plus 3.8-5.5	55062	54
Integrity 6.0 - 8.0	55040	32				TG Plus 3.4-4.5	55018	8
T Implant	55063	55	Trinon					
			Q2-Implant	55011	1	Zeramex		
TBR			QK-Implant	55014	4	XT Implant 4.2, 5.5	55086	78
Infinity PF 3.5, PF 4.0	55004	18	QZA 3.5 Platform	55049	41			
M Implant	55044	36	QZA 4.5, 5.7 Platform	55051	43	Zimmer		
						Screw-Vent, Tapered Screw-Vent 3.7-4.1	55049	41
Thommen			URIS			Screw-Vent, Tapered Screw-Vent 4.7-5.7	55051	43
SPI 3.5	55021	11	Omni Narrow 3.0-3.5	55115	81	Tapered Swiss Platt 3.8	55070	62
SPI 4.0	55065	57	Omni Regular 4.0-6.5	55020	10	Tapered Swiss Platt 4.8	55014	4
SPI 4.5-6.0	55066	58				Spline Twist Platt 4 & 5	55038	30
						Eztetic NP 3.1	55065	57
						Zuga		
						All	55027	19

E-Driver Orthodontic/Prosthetic screws




E-Driver	
Description	Image
<ul style="list-style-type: none"> • Strong and accurate application of torque (Range: 5~35Ncm) • Adjustable drilling speed (Range: 15~60rpm) • Minimizes OrthAnchor fractures and ensures accurate placement path • Ease of abutment tightening and minimizing the chance of screw loosening 	
OSM-TORQ	

Image Sensor

Pluto	
Description	Image
<ul style="list-style-type: none"> • 16bit clear image • Thinnest intraoral sensor in the market - 4.4mm • Ergonomically designed to provide optimal patient comfort • Quick image acquisition - 4 sec 	
PLUTO0001X PLUTO0002X	

FireCR	
Description	Image
<ul style="list-style-type: none"> • Clear image of 21 lp/mm • Thin and flexible IP of various size allows for patient comfort • Non-contact scan provides up to 1.8 times longer IP lifetime • Quick image acquisition minimizes patient chair time 	
FIRECR	

SNAP	
Description	Image
<ul style="list-style-type: none"> • Autofocus provides clear and vivid images with less distortion • Slim head design enables accessibility to wisdom teeth with ease • 160 degree contra angle provides ease of use • Two-way shooting buttons allows for comfortable handling of the device 	
SNAP004	


Medit Intraoral Scanner


Medit i900	
Description	Image
<ul style="list-style-type: none"> • Ultra light: 165g • Touch band and touch panel • Free applications • No Subscriptions • 3 year warranty 	
MO1-i900	


Medit i700 Wired	
Description	Image
<ul style="list-style-type: none"> • 180 Degree reversible tip • Free applications • No subscription • Remote Control Mode • UV-C Disinfection • Light scanner for comfort: 245g 	
MD-IS0200	

Medit i600	
Description	Image
<ul style="list-style-type: none"> • 180 Degree reversible tip • Free applications • No subscriptions • Light scanner for comfort: 245g 	
MD-IS0100	

TRIOS® 5 Intraoral Scanner


Trios® 5	
Description	Image
<ul style="list-style-type: none"> • FDA Approved Caries Detection Aid • Next-level ergonomics • TriosShare • ScanAssist • Remote Control Mouse Mode • Improved battery life • LED and Haptic feedback ※ Engagement apps included 	
22003185	

TRIOS 5 Protection Tip
Image

22002382


TRIOS 5 Battery Charger
Image

22002384

TRIOS 5 Battery 3-pack
Image

22002385

Trios 5 Sleeves
Image

22002380

TRIOS® 4 Wireless Pod Intraoral Scanner

Trios® 4 Wireless Pod	
Description	Image
<ul style="list-style-type: none"> • Wired and wireless (2-in-1 technology) • Smart tips (Fast heating for fogless scanning) • Remote Control Mode • Reversible tips • Built-in counter (For when to buy new tips and to calibrate) <p>※ Engagement apps included</p>	
22003127	

TRIOS 4 Protection Tip
Image


82610101

Trios 4 Tips (Box of 3)
Image

82610100

TRIOS 4 Battery 3-pack
Image

22002321


Trios 4 Calibration Kit
Image

82610103

TRIOS® 3 Intraoral Scanner

Trios® 3 (2024)	
Description	Image
<ul style="list-style-type: none"> Award-winning TRIOS 3 technology Wired Only Reversible tips ※ Engagement apps included 	
22004000	

TRIOS 3 Protection Tip
Image

22002382

TRIOS 3 Calibration Tip
Image

22002384

TRIOS 3 Calibration Kit
Image

22002385

TRIOS® Move+

Trios® MOVE+ (Standalone Accessory)	
Description	Image
<ul style="list-style-type: none">• Navigate directly on the LCD touch screen: 15.6" (MOVE+): USB 2.0 port for import and export of your scans• Maneuver the arm and swivel screen• Comes with a powerful PC• Designed for fast and smooth scanning with TRIOS• Unrivalled open system• Secure access to the most up-to-date technology, materials, partners and prices	
22002375	


K5 Dental Unit Chair

K5 Dental Unit Chair

Description

- Patient Chair - High performance electric motor with noise filter
- Cuspidor Unit - Water and air delivery and intake & outtake unit connection to junction box with 90° rotatable spittoon
- Doctor Table/Delivery System - Holds various handpieces and master control unit with LCD display
- Assistant Delivery System - Controls the chair's movement/Holds instruments
- LED Dental Light - Touchless on/off sensor
- Improved Dr. Stool design - Allows for better customization of lumbar support, this helps to reduce strain on the lower back and improve overall comfort during long procedures
- Improved chair cushioning - The upgraded cushioning materials offer better durability to maintain their shape and support over time
- Wireless Foot Control - Controls chair positioning and handpieces
- Warm water supply - Supplies 77° F water to help alleviate discomfort and soothe symptoms during treatment
- Sterilized Water Supply Device - Supplies sterilized water to effectively suppress bacteria and biofilm in the water pipe
- Wide range of color options to match any clinic's interior settings

※ 5 Year Limited Warranty, please refer to the full warranty statement



Brown

Dark Blue

Dark Violet

Classic Blue

Gold

Coral

Royal Red

Light Green

Ivory

K3™ Next Dental Unit Chair

K3™ Next Dental Unit Chair

Description

- Patient Chair - High performance hydraulic motor with noise filter
- Cuspidor Unit - Water and air delivery and intake & outtake unit connection to junction box
- Doctor Table/Delivery System - Holds various handpieces and master control unit with LCD display
- Assistant Delivery System - Controls the chair's movement/Holds instruments
- LED Dental Light - Touchless on/off sensor
- Foot Control - Controls chair positioning and handpieces
- Wide range of color options to match any clinic's interior settings

※ 5 Year Limited Warranty, please refer to the full warranty statement



Brown

Dark
Blue

Dark
Violet

Blue

Gold

Red

Ivory

T2 Plus Dental CBCT

T2 Plus Dental CBCT

Description

- Faster scan time, higher resolution and AI-powered 3D diagnostics
- 3 in 1 system: CT, Pano & Cephalo
- Multiple FOV for a Wide Range of Applications:
* 5x5, 8x9, 10x9, 12x9, 15x9 (Standard) & 15x15 (Stitching)
- High Resolution Voxel Size to Provide Clear Images
- Simplified TMJ Imaging
- User-Friendly Operation
- Automated voice guidance to reduce re-imaging
- Improved radiolucent temple support for stable patient positioning
- ※ 2 Year Full Warranty, please refer to the full warranty statement
- ※ Up to 10 Year Extended Warranty is available for purchase



T2 Dental CBCT

T2 Dental CBCT

Description

- 3 in 1 system: CT, Pano & Cephalo
- Multiple FOV for a Wide Range of Applications:
 - * 5x5, 8x9, 10x9, 12x9, 15x9 (Standard) & 15x15 (Stitching)
- High Resolution Voxel Size to Provide Clear Images
- Simplified TMJ Imaging
- User-Friendly Operation
- Automated voice guidance to reduce re-imaging
- Improved radiolucent temple support for stable patient positioning
- ※ 2 Year Full Warranty, please refer to the full warranty statement
- ※ Up to 10 Year Extended Warranty is available for purchase



HIOSSEN IMPLANT

Smiles that last a life time

Hiossen Implant Inc.

270 Sylvan Ave. Ste 1130, Englewood Cliffs, NJ 07632
www.hiossen.com Email: marketing@hiossen.com

© @Hiossen  @HiossenImplants  @Hiossen

