

Osseodensification Kit for Soft Bone


Bone Compaction Kit

- Compacts bone chips into the surrounding bone
- Enables safe sinus lift without hydraulic elevation

HIOSSEN
IMPLANT

Enhancing bone density through compression

- Unlike a surgical taper drill, a compaction drill compresses the lateral bone during the drilling process




122 Taper Drill

Bone Compaction Drill

Increasing initial fixation in D3 and D4 bones through bone densification

- Torque required for placing a $\text{Ø}4.5 \times 10\text{mm}$ implant in soft bone




122 Taper Drill

Bone Compaction Drill

Engineered to retain bone chips during drilling

- Creates bone chips and compresses the bone laterally while drilling clockwise



Maxillary sinus lift without hydraulic elevation

- Unlike our CAS drill, the compaction drill does not require hydraulic pressure to elevate the sinus membrane, reducing the risk of membrane perforation




Bone Compaction Kit (HBCK)

Implant Systems

EKIII/IV

ETIII/IV

SSIII



Ø2.2 Initial Drill

3
HBCDI03
3mm Initial Drill (Optional)

8
HBCDI04
HBCDI05
HBCDI06
HBCDI07
HBCDI08
HBCDI10

F3.5 Compaction Drill

4
HBCD3503
3mm Initial Drill (Optional)

9
HBCD3504
HBCD3505
HBCD3506
HBCD3507
HBCD3508
HBCD3510

F4.0 Compaction Drill

5
HBCD4003
3mm Initial Drill (Optional)

10
HBCD4004
HBCD4005
HBCD4006
HBCD4007
HBCD4008
HBCD4010

F4.5 Compaction Drill

6
HBCD4503
3mm Initial Drill (Optional)

11
HBCD4504
HBCD4505
HBCD4506
HBCD4507
HBCD4508
HBCD4510

F5.0 Compaction Drill

7
HBCD5003
3mm Initial Drill (Optional)

12
HBCD5004
HBCD5005
HBCD5006
HBCD5007
HBCD5008
HBCD5010

HIOSSEN

IMPLANT

Smiles that last a lifetime






Please contact your local sales representative or visit our website today to learn more about Hiossen and its products.



All Hiossen Implants are processed and
Manufactured in the USA

270 Sylvan Ave. Ste 1130, Englewood Cliffs, NJ 07632
Marketing@hiossen.com | 888.678.0001 | www.hiossen.com

 @Hiossen  @HiossenImplants  @Hiossen

Copyright © 2024 Hiossen Inc. All rights reserved.

All information included herein is subject to change without notice. No part of this brochure may be reproduced in any manner without the prior written permission of Hiossen Inc. Hiossen and the Hiossen logo are trademarks of Hiossen Inc. All other brand or product names are trademarks of their respective companies or organizations.