

Medit Link & Apps

Medit Link is a **fully open system** for you.



Dashboard



Workflow Management



Cloud Service



Patients Communication

Workflow that is patient-centered.

Value created together with partners.

Data Management that has never been more convenient.


Experience of the real digital dentistry with Medit Link.

Scan with **free intelligent features** and **apps specialized** for your dental procedures.

Improve consultation experience with your patient to **provide better treatment options.**



Temporaries



Model Builder



Crown Fit



Compare



Smile Design



DCM Converter



Occlusion Analyzer



Ortho Simulation



Calibration Wizard

About Medit

Medit specializes in 3D measurement and CAD/CAM solutions for dental clinics and labs, including intraoral scanners, based on its own patented state-of-the-art technology. We also develop platform solutions for digital dentistry, supporting collaborative workflows.

Our goal is to provide innovative technology and the highest quality products to ensure mutual growth for all partners.



User Group



Medit Academy



Help Center



Contact Medit



www.medit.com


Unleash your clinic's potential

MEDIT *i700*



Medit i700 makes the scanning experience a comfortable one for both the dentist and patient.

With both **powerful hardware** and **intelligent software**, the full-feature Medit i700 is the key to **unleash your clinic's full potential**.



Specifications

Category	Medit i700	
Scanning technology	Scanning frame	Up to 70 FPS
	Imaging technology	3D-in-motion video technology 3D full color streaming capture
	Light source	LED
Accuracy	Anti-fogging technology	Adaptive anti-fogging
Handpiece	Full-arch	10.9µm ± 0.98
	Dimensions	248 x 44 x 47.4 mm
	Weight	245 g
Tip	Tip size	22.2 x 15.9 mm
	Mirror angle	45-degree angle <small>Easier to scan distal molar area</small>
	Scan area	15 x 13 mm
	Autoclavable	Up to 100 times Autoclave 121°C 30 mins Autoclave 134°C 4 mins
	Reversible tip	Yes
Special mode	Remote control mode	Yes
	UV-C disinfection	Yes
Cable	Length	2.0 m / Handpiece Detachable <small>2.5m (Sold Separately)</small>
	Connectivity	USB 3.1 Gen1 (C Power Delivery)

Please scan the QR code for up-to-date [system requirements](#). >>

