



Precise V.P.S. Impression Material

Hysil Plus


- Hydrophilic impression material for precise impression taking
- Easily removable impression material with high tensile strength
- Short Intra-oral Setting Time

HIOSSSEN
IMPLANT



Impression material with hydrophilic properties for accurate impression taking


- Angle testing with Hysil Plus Light Body compared to brands A & B
- It exhibits hydrophilicity that is more than double that of its competitors, allowing for accurate impressions




HySil Plus shows twice the hydrophilicity compared to Products A and B

Impression material with high tensile strength to facilitate effortless removal

- Facilitate easy removal and prevent tearing, even in the finer marginal areas



Capable of capturing impressions for various types of margins



Inlay



Single Crown




Partial Denture



Transfer Impression
(Fixture Level)

Intra-oral Setting Time

- Sufficient working time of 3min 50 second
- Short intraoral setting time of 2min 30sec for fast and precise impression taking



Hiossen Impression Materials Line ups

| HySil Plus Impression Materials

Regular Type

HySil Plus Heavy



ESS50HS (50ml X 4)
ESS50HSB (50ml X 15)

HySil Plus Mono



ESS50MS (50ml X 4)
ESS50MSB (50ml X 15)

HySil Plus Light



ESS50LS (50ml X 4)
ESS50LSB (50ml X 15)

HySil Plus Putty



ESS400P

HySil Plus Bite



ESS50B

Auto Type

HySil Plus Heavy Auto



ESS380HS

HySil Plus Mono Auto



ESS380MS