

# Get a **FREE** Night Guard

## Digital Center Special Offer!

**SEP 1 – NOV 30, 2025**

- For Hiossen Implant Package customers
- **1 FREE** custom night guard (*value \$100*)
- Choice: Soft (*comfortable*) or Hard (*durable*)
- Fabricated by Digital Center with precision

### Benefit

Protect patients' teeth  
& your treatment outcomes

*Ask your Hiossen rep today!*




#### Disclaimer:

Offer valid for Hiossen Implant Package customers only. One free night guard per customer. Valid September 1 – November 30, 2025. Non-transferable, no cash value.

**HÍOSSEN**  
IMPLANT


# Night Guard

## How to Claim



The dashboard shows the HiOSSEN account number (07632059) and digital site test information. It features four main buttons: 'Create Case' (red), 'Case Manager', 'Message Center', and 'Profile'.

**Create a Case:** Go to HiOSSEN OneDigital Website, create a new case for your patient.




The 'Create Case' form is shown at the 'Teeth Selection' step. It includes fields for Patient Name (First and Last name) and a diagram of teeth numbered 1-16 for the upper arch and 32-17 for the lower arch. Buttons for 'Upper Arch' and 'Lower Arch' are present, along with a 'Next' button.

**Choose the Arch:** Select the arch (upper or lower) for the night guard.



The 'Create Case' form is shown at the 'Arch' selection step. It displays the 'Night Guard' option checked under the 'Arch' section. The 'Next' button is visible at the bottom right.

**Check the Night Guard Option:** Locate the Night Guard option box and check it.



The 'Create Case' form is shown at the 'Instructions' step. It includes a text field for 'Free night guard' and checkboxes for 'Enclosed' and 'Additional' options. The 'Next' button is visible at the bottom right.

**Add Special Instructions:** In the instruction field, type: "Free night guard".

**Complete Shipping Information:** Enter and confirm the shipping details to finalize the case.

# Night Guard

## Use and Care Instruction

### Use instruction

#### Inserting the night guard

- Insert over your teeth. Gently press the nightguard into place.
- For thermoplastic night guards: Place in warm to hot water for 15 seconds to soften it before inserting.

#### Removing the night guard

- Remove by gently pulling evenly on both sides.
- Avoid using excessive force.

### Care instruction

#### Cleaning

- Rinse with cool water and gently brush with a soft toothbrush.
- Soak weekly in denture cleaner or non-alcoholic mouthwash for 20 mins.
- Rinse before use.

#### Storage

- Store in a clean, dry protective case to prevent bacterial buildup.
- Keep away from sunlight, hot water, or high temperatures.