

INTRODUCING ZCORETM EXPAND

expandable collagen-enriched porcine xenograft

- Prevents risks of membrane dislodging or particulate washout
- Resorption Profile: implant placement at 18 - 20 weeks post extraction grafting




Expandable Xenograft Collagen

- Composed of 65% porcine xenograft particulate and 35% porcine collagen by volume
- Made from porcine collagen and porcine graft material, sharing the same material origin
- Commonly used for extraction site grafting and sinus augmentation

Zcore Expand Function



Hydrates and expands almost immediately when introduced to the patient's blood or sterile saline



Once hydrated, it increases in diameter to fill the void space in a socket or sinus defect

65% Porcine
xenograft
particulate

35% Type 1
expandable
porcine collagen

100X



Zcore Expand Specification (3 Types)



Socket

5 mm x 17 mm | **10 mm x 17 mm EXPANDED**

| ZXSOCKET-1 (1 per box) | ZXSOCKET-5 (5 boxes bundle)



Small Sinus

13 mm diam. x 10 mm | **17 mm x 10 mm EXPANDED**

| ZX SINUSS (1 per box)



Large Sinus

17 mm diam. x 12 mm | **22 mm x 12 mm EXPANDED**

| ZX SINUSL (1 per box)



INTRODUCING ZCORETM FORM

moldable collagen-enriched porcine xenograft



- A composite of osteoconductive bone mineral and collagen

Xenograft Collagen Bone Block

- Composed of 80% porcine xenograft particulate and 20% porcine collagen by volume sharing the same material origin
- Its moldable consistency allows it to conform to the shape of the defect
- Offers easier and more convenient handling compared to particulate grafts


Zcore Form Function



Hydrates almost immediately when introduced to the patient's blood or sterile saline



Once hydrated, it becomes moldable and can take the shape of a variety of defect shapes and sizes



ZcoreTM Form Moldable Collagen-Enriched Porcine Xenograft



0.5 cc

9 mm diam. x 8 mm
| ZF050 (1 per box)




1.0 cc

11 mm diam. x 12 mm
| ZF100 (1 per box)



2.0 cc

11 mm diam. x 22 mm
| ZF200 (1 per box)



PF26ZCOLTR1.0

Copyright © 2026 Hiossen Inc. All rights reserved.
105 Challenger Rd. 7th Floor, Ridgefield Park, NJ 07660
Marketing@hiossen.com | 888.678.0001 | www.hiossen.com

@Hiossen @HiossenImplants @Hiossen

Hiossen
IMPLANT

Smiles that last a lifetime